

PICam™ 4.x Programmer's Manual

PICam™ 4.x

Revision History

Issue	Date	List of Changes
Issue 1	September 13, 2016	This is the initial release of this document.

©Copyright 2016

Princeton Instruments, a division of Roper Scientific, Inc.
3660 Quakerbridge Rd
Trenton, NJ 08619
TEL: 800-874-9789 / 609-587-9797
FAX: 609-587-1970

All rights reserved. No part of this publication may be reproduced by any means without the written permission of Princeton Instruments, a division of Roper Scientific, Inc. ("Princeton Instruments").

Printed in the United States of America.

PICam is a trademark of Roper Scientific, Inc.

Roper Scientific is a registered trademark of Roper Scientific, Inc.

Windows, Windows Vista, Windows 7, Windows 8, and Windows 10 are registered trademarks of Microsoft Corporation in the United States and/or other countries.

Red Hat and Enterprise Linux are registered trademarks of Red Hat, Inc.

The information in this publication is believed to be accurate as of the publication release date. However, Roper Scientific, Inc. does not assume any responsibility for any consequences including any damages resulting from the use thereof. The information contained herein is subject to change without notice. Revision of this publication may be issued to incorporate such change.

Table of Contents

Chapter 1:	About this Manual	7
1.1	Manual Organization	7
1.2	Contact Information	8
Chapter 2:	Introduction to PICam™	9
2.1	System Overview	9
2.2	Hardware Support	9
2.2.1	Camera Firmware [GigE Cameras Only]	9
2.3	Supported Operating Systems	9
2.3.1	WoW64 Support	10
2.4	Sample Code	10
2.5	Naming Conventions	11
2.6	Concepts	12
2.6.1	Camera Handles	13
2.7	Defined Data Types	14
2.8	Include Files	15
2.8.1	Optional and Advanced Files	15
Chapter 3:	General Library APIs	17
3.1	Data Type Definitions	17
3.2	Programmers' Reference for General Use Library APIs	18
Chapter 4:	Camera Identification APIs	25
4.1	Data Type Definitions	25
4.2	Structure Definitions	26
4.3	Programmers' Reference for Camera APIs	27
4.3.1	Camera Identification APIs	28
4.3.2	Camera Access APIs	33
4.3.3	Camera Information APIs	41
4.3.4	Demo Camera Identification APIs	43
Chapter 5:	Camera Configuration APIs	49
5.1	Data Type Definitions	50
5.1.1	Camera Parameter Enumerations	50
5.1.2	Camera Parameter Access Enumerations	81
5.1.3	Camera Parameter Constraint Enumerations	82
5.2	Structure Definitions	85
5.2.1	Camera Parameter Structures	85
5.2.2	Camera Parameter Constraints	89
5.3	Programmers' Reference for Camera Configuration APIs	95
5.3.1	Camera Parameter Value APIs	97
5.3.2	Camera Parameter Information APIs	137
5.3.3	Camera Parameter Constraints APIs	145
5.3.4	Camera Parameter Commitment APIs	155

Chapter 6:	Camera Data Acquisition APIs	157
6.1	Data Format	158
6.2	Data Type Definitions	159
6.2.1	Data Acquisition Enumerations	159
6.3	Data Acquisition Data Structures	160
6.4	Programmers' Reference for Acquisition Control APIs	162
Chapter 7:	Advanced Function APIs	169
7.1	Data Type Definitions	170
7.1.1	Camera Plug and Play Discovery Data Enumerations	170
7.1.2	Camera Access Enumerations	171
7.1.3	Camera Parameter Information Enumerations	171
7.1.4	Camera Data Acquisition Enumerations	172
7.2	Data Structures	173
7.2.1	Camera Information Data Structures	173
7.2.2	Camera Parameter Validation Data Structures	176
7.2.3	Camera Data Acquisition Data Structures	180
7.3	Callback Functions	181
7.3.1	Camera Discovery Callbacks	181
7.3.2	Camera Parameter Value Callbacks	182
7.3.3	Camera Parameter Constraints Callbacks	187
7.3.4	Camera Data Acquisition Callbacks	190
7.4	Programmers' Reference for Advanced APIs	191
7.4.1	Camera Discovery APIs	193
7.4.2	Camera Access APIs	198
7.4.3	Camera Information APIs	204
7.4.4	Camera Parameter Value APIs	208
7.4.5	Camera Parameter Information APIs	227
7.4.6	Camera Parameter Constraints APIs	233
7.4.7	Camera Commitment APIs	248
7.4.8	Acquisition Setup APIs	257
7.4.9	Acquisition Notification APIs	259
7.4.10	Acquisition State Notification APIs	261
7.4.11	Acquisition Control APIs	264
Chapter 8:	EM Calibration APIs	267
8.1	EM Calibration Applications	267
8.2	Structure Definitions	268
8.2.1	EM Calibration Structures	268
8.3	Callback Functions	269
8.3.1	EM Calibration	269
8.4	Programmers' Reference for EM Calibration APIs	270
8.4.1	EM Calibration Access APIs	271
8.4.2	EM Calibration Parameter Value APIs	275
8.4.3	EM Calibration Parameter Constraints APIs	280
8.4.4	EM Calibration APIs	281
Appendix A:	Available Parameters	283
Appendix B:	EM Gain Calibration Code Sample	289
B.1	EM Gain Calibration Procedure	289

Appendix C: Firmware Upgrade/Restore	291
C.1 Firmware Upgrade Procedure	291
C.2 Restore Firmware	293
C.2.1 Precautions	293
C.2.2 Procedure	294
Appendix D: Debugging GigE Cameras	297
D.1 Debugging	297
D.1.1 Timeout Period Considerations	298
D.1.2 Following Debugging	298
D.2 Timeout Configuration	298

List of Figures

Figure 2-1: Basic PICam Structure	12
Figure 2-2: Block Diagram of Handle Hierarchy	13
Figure 6-1: Data Format Diagram	158
Figure 7-1: PICam Structure - Advanced	169
Figure B-1: Typical EM Gain Calibration Dialog	290
Figure C-1: Firmware Upgrade: Typical IP Engine Selection Dialog	291
Figure C-2: Firmware Upgrade: Selecting Device to be Upgraded	292
Figure C-3: Firmware Upgrade: Typical Updating Dialog	292
Figure C-4: Firmware Upgrade: Upgrade Complete	293
Figure C-5: Firmware Restore: Typical IP Engine Selection Dialog	294
Figure C-6: Firmware Restore: Selecting Device to be Restored	294
Figure C-7: Firmware Restore: Typical Updating Dialog	295
Figure C-8: Firmware Restore: Complete	295

List of Tables

Table 2-1: List of Sample Code Files Provided	10
Table 2-2: Data Type Definitions	14
Table 2-3: Sized Data Type Definitions	14
Table 3-1: Data Enumeration Definitions for General Library APIs	17
Table 4-1: Data Definitions for Camera APIs	25
Table 5-1: PicamValueType Enumerator Definitions	51
Table 5-2: PicamConstraintType Enumerator Definitions	52
Table 5-3: PicamParameter Enumerator Definitions	53
Table 5-4: PicamActiveShutter Enumerator Definitions	62
Table 5-5: PicamAdcAnalogGain Enumerator Definitions	62
Table 5-6: PicamAdcQuality Enumerator Definitions	63
Table 5-7: PicamCcdCharacteristicsMask Enumerator Definitions	64
Table 5-8: PicamCoolingFanStatus Enumerator Definitions	65
Table 5-9: PicamEMICcdGainControlMode Enumerator Definitions	65
Table 5-10: PicamGateTrackingMask Enumerator Definitions	66
Table 5-11: PicamGatingMode Enumerator Definitions	66
Table 5-12: PicamGatingSpeed Enumerator Definitions	67
Table 5-13: PicamIntensifierOptionsMask Enumerator Definitions	67

Table 5-14: PicamIntensifierStatus Enumerator Definitions	68
Table 5-15: PicamModulationTrackingMask Enumerator Definitions.....	68
Table 5-16: PicamOrientationMask Enumerator Definitions	69
Table 5-17: PicamOutputSignal Enumerator Definitions	69
Table 5-18: PicamPhosphorType Enumerator Definitions	70
Table 5-19: PicamPhotocathodeSensitivity Enumerator Definitions	71
Table 5-20: PicamPhotonDetectionMode Enumerator Definitions.....	72
Table 5-21: PicamPixelFormat Enumerator Definitions.....	72
Table 5-22: PicamReadoutControlMode Enumerator Definitions	73
Table 5-23: PicamSensorTemperatureStatus Enumerator Definitions	73
Table 5-24: PicamSensorType Enumerator Definitions.....	74
Table 5-25: PicamShutterStatus Enumerator Definitions	74
Table 5-26: PicamShutterTimingMode Enumerator Definitions	75
Table 5-27: PicamShutterType Enumerator Definitions.....	75
Table 5-28: PicamTimeStampsMask Enumerator Definitions	76
Table 5-29: PicamTriggerCoupling Enumerator Definitions	77
Table 5-30: PicamTriggerDetermination Enumerator Definitions.....	78
Table 5-31: PicamTriggerResponse Enumerator Definitions	79
Table 5-32: PicamTriggerSource Enumerator Definitions	79
Table 5-33: PicamTriggerTermination Enumerator Definitions.....	80
Table 5-34: PicamValueAccess Enumerator Definitions.....	81
Table 5-35: PicamConstraintScope Enumerator Definitions	82
Table 5-36: PicamConstraintSeverity Enumerator Definitions	82
Table 5-37: PicamConstraintCategory Enumerator Definitions	83
Table 5-38: PicamRoisConstraintRulesMask Enumerator Definitions	84
Table 6-1: PicamAcquisitionErrorsMask Enumerator Definitions.....	159
Table 7-1: PicamDiscoveryAction Enumerator Definitions	170
Table 7-2: PicamHandleType Enumerator Definitions.....	171
Table 7-3: PicamDynamicsMask Enumerator Definitions	171
Table 7-4: PicamAcquisitionState Enumerator Definitions	172
Table 7-5: PicamAcquisitionStateErrorsMask Enumerator Definitions ...	172
Table A-1: Symbol Key for Table A-2	283
Table A-2: Parameter Information and Camera Support	283

Chapter 1: About this Manual

This manual describes terms and concepts used in PICam and provides descriptions of functions, parameters, and values used to create a user-designed interface to Princeton Instruments cameras. The majority of the manual focuses on the basic PICam functions (`picam.h`). The remainder focuses on the more complex PICam functions (`picam_advanced.h`) and on the EM gain calibration functions (`picam_em_calibration.h`).



NOTE:

Functions that are specific to a particular OEM are included in `picam_special.h` and are not described in this manual.

1.1 Manual Organization

This manual includes the following chapters:

- [Chapter 1, About this Manual](#)
This chapter provides general information about this manual, as well as contact information for Princeton Instruments.
- [Chapter 2, Introduction to PICam™](#)
This chapter provides information about concepts, terms, and data types used in PICam. It also provides information about the general sequence of making functions calls when writing a program.
- [Chapter 3, General Library APIs](#)
Provides programming reference information for each of the basic functions (`picam.h`).
- [Chapter 4, Camera Identification APIs](#)
Provides programming reference information for each of the basic functions (`picam.h`).
- [Chapter 5, Camera Configuration APIs](#)
Provides programming reference pages for each of the basic functions (`picam.h`).
- [Chapter 6, Camera Data Acquisition APIs](#)
Provides programming reference pages for each of the basic functions (`picam.h`).
- [Chapter 7, Advanced Function APIs](#)
Provides programming reference information about advanced functions included in `picam_advanced.h`.
- [Chapter 8, EM Calibration APIs](#)
Provides programming reference information for EM Calibration functions included in `picam_em_calibration.h`.

- [Appendix A, Available Parameters](#)
Provides parameter information and camera support for customer-accessible parameters.
- [Appendix B, EM Gain Calibration Code Sample](#)
Provides information about building and using the `EMGainCalibration.exe` sample file included with PICam.
- [Appendix C, Firmware Upgrade/Restore](#)
Provides information about upgrading GigE camera firmware to be compatible with PICam 4.x. Information is also provided about restoring firmware to PICam 3.x.
- [Appendix D, Debugging GigE Cameras](#)
Provides information about using the Heartbeat Timeout system variable.

Wherever possible, this manual uses the headings in the PICam header files (i.e., `pil_platform.h`, `picam.h`, `picam_advanced.h`, and `picam_em_calibration.h`) when grouping functions.

1.2 Contact Information

Princeton Instruments' manufacturing facility is located at the following address:

Princeton Instruments
3660 Quakerbridge Road
Trenton, NJ 08619 (USA)
TEL: 800-874-9789 / 609-587-9797
TEL: 609-587-1970

Customer Support E-mail: techsupport@princetoninstruments.com

For technical support and service outside the United States, refer to the Princeton Instruments web page at www.princetoninstruments.com. An up-to-date list of addresses, telephone numbers, and e-mail addresses for Princeton Instruments' overseas offices and representatives is maintained on the web page.

Chapter 2: Introduction to PICam™

PICam is an ANSI C library of camera control and data acquisition functions.

2.1 System Overview

To use PICam, a system must include supported camera hardware, and a host computer with the PICam runtime installed.

2.2 Hardware Support

Version 4.x of the PICam library supports the following Princeton Instruments hardware:

- PI-MAX3
- PI-MAX4, PI-MAX4:RF, PI-MAX4:EM
- PI-MTE
- PIoNIR/NIRvana
- NIRvana-LN
- PIXIS, PIXIS-XB, PIXIS-XF, PIXIS-XO
- ProEM
- ProEM+
- ProEM-HS
- PyLoN
- PyLoN-IR
- Quad-RO
- SOPHIA

2.2.1 Camera Firmware [GigE Cameras Only]

For GigE cameras, PICam 4.x is not backwards compatible with prior releases of PICam. Therefore, when using PICam 4.x with any GigE camera, the camera's firmware must be PICam 4.x compatible. Upgrading PICam 3.x camera firmware is easily achieved using the Upgrade Tool supplied by Princeton Instruments.

The key symptom of a firmware mismatch between PICam and a GigE camera is the inability to see the camera from within PICam. When this occurs, the firmware within the camera must be updated to be compatible with the version of PICam being used.

- For information about installing PICam 4.x firmware onto a GigE camera with PICam 3.x firmware, refer to [Section C.1, Firmware Upgrade Procedure](#), on page 291.
- For information about restoring firmware, refer to [Section C.2, Restore Firmware](#), on page 293.

2.3 Supported Operating Systems

PICam currently supports the following 64-bit operating systems:

- Windows Vista®;
- Windows® 7;
- Windows 8/8.1
- Windows 10
- RedHat® Enterprise Linux®, version 6.4 (RHEL6.4)

In the future, the functions described in this manual may work with additional operating systems.

2.3.1 WoW64 Support

PICam supports WoW64 which enables 32-bit programs to work with PICam and operate Princeton Instruments detectors in a 64-bit operating system.



NOTE:

64-bit programs link with `picam.dll`.

32-bit programs link with `picam32.dll`.

2.4 Sample Code

Code samples are provided with PICam. When the PICam Software Development Kit (SDK) is installed, these samples are installed, by default, in the PICam installation directory.



NOTE:

The specific directory in which code samples are installed varies by operating system.

Table 2-1: List of Sample Code Files Provided (Sheet 1 of 2)

Sample Name	Description
AcquisitionState	This sample demonstrates an advanced acquisition scenario where the program can be notified when the camera transitions through important acquisition states (e.g., the beginning of readout.)
Acquire	This is the basic data acquisition sample. It calls <code>Picam_Acquire()</code> and waits for all frames to be completed. The second part of this sample waits in a loop for N frames, acquiring 1 frame at a time.
Advanced	This sample illustrates features of <code>picam_advanced.h</code> .
Configure	This sample illustrates how to change settings during camera setup as well as online while polling for data.
EMGainCalibration	This sample illustrates how to set up EM Gain Calibration. For additional information about incorporating this sample into production code, refer to Appendix B, EM Gain Calibration Code Sample .
Gating	This sample illustrates how to set up repetitive and sequential gating. Aso demonstrates RF features on cameras which support RF functionality.
Kinetics	This sample provides a sequence of API calls used to request acquisition of image data using the kinetics window capture mode. The demo also illustrates how to make calls to utilize external triggering of captures. The captured pixel data are stored to a raw data file.
Metadata	This sample enables metadata (i.e., Time Stamp(s) and Frame Tracking.) It illustrates how to extract metadata from the data stream.
MultiCam	This example opens multiple (i.e., 2,) cameras and collects data from all simultaneously.
ParamInfo	This sample accesses all parameter information for all camera parameters, and then prints them to the screen.

Table 2-1: List of Sample Code Files Provided (Sheet 2 of 2)

Sample Name	Description
Poll	This sample illustrates how to use the polling method for collecting data by using <code>Picam_WaitForAcquisitionUpdate()</code> .
Rois	This sample demonstrates the API for setting a simple single region of interest. It also shows how to set up a camera for multiple regions of interest and then acquires data for the given region(s).
SaveData	This sample acquires data synchronously and writes the returned data buffer to disk.
WaitForTrig	This sample waits for an external trigger to start data acquisition.

2.5 Naming Conventions


The following naming conventions are used in PICam:

- All primitive types have a typedef with a `pi` prefix (e.g., `piint`, `pi64s`.)
- All functions defined by PICam are prefixed with `PICam_` and return an error code of `PicamError` (e.g., `Picam_GetParameterIntegerValue`, `Picam_CloseCamera`.)
- All functions that allocate memory to store the results of a function call return a pointer to a constant allocation of the appropriate type. For example:
 - `Picam_GetEnumerationString` returns a string by taking the address of a pointer to a constant string. In other words an argument to the function is `const pichar**`.
 - `Picam_GetParameterCollectionConstraint` returns a collection constraint by taking the address of a pointer to a constant collection constraint. In other words, an argument to the function is `const PicamCollectionConstraint**`.
- All functions that allocate an array of memory to store the results of the function call return a pointer to a constant array allocation of the appropriate type as well as the number of items in the array.
For example, `Picam_GetParameters` returns an array of parameters by taking the addresses of a pointer to a constant parameter array and a count. In other words, two arguments to the function are `const PicamParameter**` and `piint*`.
- All functions that free memory allocated by PICam have a `Picam_Destroy` prefix (e.g., `Picam_DestroyString`, `Picam_Destroy_CollectionConstraints`, `Picam_DestroyRois`.)
- All types defined by PICam are prefixed with `PICam` and have a typedef to `<TypeName>` (e.g., `PicamParameter`, `PicamRoi`.)
- All enum type members defined by PICam are prefixed with `<EnumName>_` (e.g., `PicamValueType` enum has a `PicamValueType_Integer` constant.)
- All enum types that represent multiple values with bitmasks have a `Mask` suffix (e.g., `PicamCcdCharacteristicsMask`, `PicamTimeStampsMask`.)

2.6 Concepts

Figure 2-1 is a high-level block diagram of the basic PICam structure. Cameras that are powered on and plugged into the host computer are initially represented by camera IDs. The content of the camera ID will be unique for each camera. From the camera ID, basic information can be garnered such as availability and basic camera information. It is also from a camera ID that a camera can be opened. Once opened, the camera can be configured by adjusting the values of its parameters. The permitted values a parameter can take are defined by its constraints. Different cameras not only possess different parameters, but different rules for interacting with those parameters. This information for each parameter may also be queried. Once a camera has been configured, data can be acquired from it.

Figure 2-1: Basic PICam Structure



2.6.1 Camera Handles

Most PICam APIs require handles to identify the specific camera with which they are currently interacting. When a camera is brought online, it is assigned a specific handle that is then used to identify the camera throughout the active session.

The following handle(s) may be passed as an API parameter:

- `camera`

Identifies a specific physical camera within the system.

When `camera` is passed to an API, PICam determines the appropriate actions depending on the API that has been called.

This handle is passed as a Basic API parameter.

- `device`

Identifies a specific PHYSICAL camera within the system.

When `device` is passed to an API, any resulting interaction or configuration performed by the API is done on a physical camera that is attached to the system.

This handle is passed as an Advanced API parameter, and must be used in conjunction with `model`.

- `model`


Identifies a specific VIRTUAL camera within application memory.

When `model` is passed to an API, any parameter configuration changes are temporarily stored in system memory (i.e., within the host computer.) The actual camera configuration remains unchanged until an API is called that commits values to the `device` (i.e., the physical camera.)

This handle is passed as an Advanced API parameter, and must be used in conjunction with `device`.

Figure 2-2 illustrates the hierarchical relationship between PICam handles and with which set of APIs they are valid.

Figure 2-2: Block Diagram of Handle Hierarchy



2.7 Defined Data Types

The typedefs are given in the header file `pic_platform.h`.

Table 2-2: Data Type Definitions

Type	Definition
piint	Integer native to platform
piflt	Floating point native to platform
pibln	Boolean native to platform
pichar	Character native to platform
pibyte	Byte native to platform
pibool	C++ Boolean native to platform

Table 2-3: Sized Data Type Definitions

Type	Definition
pi8s	8-bit signed integer
pi8u	8-bit unsigned integer
pi16s	16-bit signed integer
pi16u	16-bit unsigned integer
pi32s	32-bit signed integer
pi32u	32-bit unsigned integer
pi64s	64-bit signed integer
pi64u	64-bit unsigned integer
pi32f	32-bit floating point
pi64f	64-bit floating point

2.8 Include Files

Any program using PICam must include the following header files:

- `pic_platform.h`
Princeton Instruments' library platform support. This is included indirectly via `picam.h`.
- `picam.h`
Princeton Instruments' camera control Application Programming Interface (API.)

2.8.1 Optional and Advanced Files

The following files are optional and only need to be included when one or more of the functions found within them are required:

- `picam_special.h`
Only include `picam_special.h` when using a special function defined in that file.
- `picam_advanced.h`
This is the Princeton Instruments advanced camera control API.

This header file contains advanced functionality such as camera discovery, change notification, circular buffering, user state, defect map, and data acquisition callbacks.
- `picam_em_calibration.h`
This header EM Gain Calibration file provides the APIs and functionality needed to perform EM gain calibration for a ProEM camera.

This page is intentionally blank.

Chapter 3: General Library APIs

The first section of `picam.h` includes functions to determine if the PICam library has been initialized, to initialize the library, uninitialize it, or get the version. This section also includes error codes that may be returned during library initialization, camera initialization, demo and camera set up, and data acquisition.

The first step in using the PICam library is library initialization. This is typically done at the start of the program. Once the library has been initialized, PICam function can then be called. The success of every function call is determined by the error code that is returned. It is paramount this error code be checked as most results are invalidated if a function fails. To facilitate debugging, PICam can convert an error code into a string. (In fact, any PICam enum can be converted into a string.) Once the program is finished with the library, it should clean up and uninitialized the library. This often occurs during program shutdown.

3.1 Data Type Definitions

Refer to [Table 3-1](#) for information about data definitions.

Table 3-1: Data Enumeration Definitions for General Library APIs

Name	Type	Description
<code>PicamError</code>	enum	The set of error codes returned from all APIs declared as <code>PICAM_API</code> .
<code>PicamEnumeratedType</code>	enum	The set of all PICam enumeration types.

3.2 Programmers' Reference for General Use Library APIs

This section provides a detailed programmers' reference guide for the following APIs, including their syntax and behavior:

- Library Version
 - `Picam_GetVersion()`
- Library Initialization
 - `Picam_IsLibraryInitialized()`
 - `Picam_InitializeLibrary()`
 - `Picam_UninitializeLibrary()`
- General String Handling
 - `Picam_DestroyString()`
 - `Picam_GetEnumerationString()`

3.2.1 `Picam_GetVersion()`

Description

`Picam_GetVersion()` returns PICam version information.

The following version information may be requested:

- Major

This is the Major release version which is incremented with each major feature addition or breaks backward-compatibility.
- Minor

This is the Minor release version which is incremented with minor feature additions.
- Distribution

This is the Distribution version which is incremented with bug fix releases.
- Released

This is the date of the current official release in the format **YYMM**.

When a release is classified as a Beta release, requesting this information returns a zero (0).



NOTE:

`Picam_GetVersion()` may be called prior to initializing the library with `Picam_InitializeLibrary()`.

Syntax

The syntax of `Picam_GetVersion()` is:

```
PICAM_API Picam_GetVersion(
    piint* major,
    piint* minor,
    piint* distribution,
    piint* released);
```

Input Parameters

Input parameters for `Picam_GetVersion()` are:

`major`: Used to request Major version.

Valid values are:

- `&major`
Indicates that the Major version is to be returned.
- `0/null`
Indicates that the Major version is not to be returned.

`minor`: Used to request Minor version.

Valid values are:

- `&minor`
Indicates that the Minor version is to be returned.
- `0/null`
Indicates that the Minor version is not to be returned.

`distribution`: Used to request Distribution version.

Valid values are:

- `&distribution`
Indicates that the Distribution version is to be returned.
- `0/null`
Indicates that the Distribution version is not to be returned.

`released`: Used to request official Release date.

Valid values are:

- `&released`
Indicates that the Release date is to be returned.
- `0/null`
Indicates that the Release date is not to be returned.

Output Parameters

Output Parameters for `Picam_GetVersion()` are:

`major`: Returns the Major version.

`minor`: Returns the Minor version.

`distribution`: Returns the Distribution version.

`released`: Returns the Released version.

Examples

If the PICam version is **4.2.1.1006**, it indicates the following version information:

- Major version: **4**
- Minor version: **2**
- Distribution version: **1**
- Release Date: **1006** [i.e., June, 2010.]

Similarly, if the PICam version is **5.1.2.0**, it indicates the following version information:

- Major version: **5**
- Minor version: **1**
- Distribution version: **2**
- Release Date: **0** indicating a Beta release.

3.2.2 Picam_IsLibraryInitialized()

Description

`Picam_IsLibraryInitialized()` determines if the library has been initialized.



NOTE:

`Picam_IsLibraryInitialized()` may be called prior to initializing the library using `Picam_InitializeLibrary()`.

Syntax

The syntax of `Picam_IsLibraryInitialized()` is:

```
PICAM_API Picam_IsLibraryInitialized (pibln* inited);
```

Input Parameters

There are no input parameters associated with `Picam_IsLibraryInitialized()`.

Output Parameters

Output parameters for `Picam_IsLibraryInitialized()` are:

`inited`: Indicates the initialization status for the library.

Valid values are:

- `True`
Indicates that the library has been initialized.
- `False`
Indicates that the library remains uninitialized.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_InitializeLibrary()`

3.2.3 Picam_InitializeLibrary()

Description

`Picam_InitializeLibrary()` initializes the library and prepares it for use.

Syntax

The syntax of `Picam_InitializeLibrary()` is:

```
PICAM_API Picam_InitializeLibrary (void);
```

Usage

Unless specifically noted otherwise, `Picam_InitializeLibrary()` **MUST** be called prior to calling any additional Library API routine.



NOTE:

`Picam_UninitializeLibrary()` **MUST** be called prior to program termination.

Input Parameters

There are no input parameters associated with `Picam_InitializeLibrary()`.

Output Parameters

There are no output parameters associated with `Picam_InitializeLibrary()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_UninitializeLibrary()`

3.2.4 Picam_UninitializeLibrary()

Description

`Picam_UninitializeLibrary()` frees resources that have been used by the API Library, including open cameras and memory.



NOTE:

`Picam_UninitializeLibrary()` **MUST** be called prior to program termination.

Syntax

The syntax of `Picam_UninitializeLibrary()` is:

```
PICAM_API Picam_UninitializeLibrary (void);
```

Input Parameters

There are no input parameters associated with `Picam_UninitializeLibrary()`.

Output Parameters

There are no output parameters associated with `Picam_UninitializeLibrary()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_IsLibraryInitialized()`
- `Picam_InitializeLibrary()`

3.2.5 Picam_DestroyString()

Description

`Picam_DestroyString()` releases PICam- allotted memory that has been associated with a specified character string, `s`.

**NOTE:**

If the character string, `s`, is null, `Picam_DestroyString()` has no effect.

Syntax

The syntax of `Picam_DestroyString()` is:

```
PICAM_API Picam_DestroyString(  
    const pichar* s);
```

Input Parameters

Input parameters for `Picam_DestroyString()` are:

`s`: Pointer to the character string for which memory is to be released.

Output Parameters

There are no output parameters associated with `Picam_DestroyString()`.

3.2.6 Picam_GetEnumerationString()

Description

`Picam_GetEnumerationString()` determines what enumeration strings have been defined for the specified enumerated type. Returns an allocated string representation of the enumeration type with value in `s`.



NOTE:

`Picam_DestroyString()` must be called to free the allocated memory associated with string `s`.

Syntax

The syntax of `Picam_GetEnumerationString()` is:

```
PICAM_API Picam_GetEnumerationString(  
    PicamEnumeratedType type,  
    piint value,  
    const pichar** s);
```

Input Parameters

Input parameters for `Picam_GetEnumerationString()` are:

- `type`: The type for which enumeration strings are being requested.
- `value`: The numeric value associated with enumeration string being requested.

Output Parameters

Output parameters for `Picam_GetEnumerationString()` are:

- `s`: Pointer to the enumeration string.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_DestroyString()`

Chapter 4: Camera Identification APIs

The APIs in this section of `picam.h` deal with determining what cameras or demo cameras are available or being used in another instance, retrieving information from camera firmware, opening and closing a camera, and connecting/disconnecting a demo camera.

Once the library has been initialized, all cameras that are powered on and connected to the host computer will have a corresponding camera ID. Accessing a camera is as simple as opening an available camera with its corresponding ID.

A demo camera is a software-simulated camera. This allows program development without a camera connected. A demo camera can be instantiated by choosing a particular camera model and connecting it. Once connected, it can be interacted with as any other camera.

Once a camera (possibly a demo) is no longer used, it should be closed.

The following factors affect camera availability to the program:

- Connectivity

In order for a camera to be detected by the program it must be:

- Connected to the host computer;
- The camera must be powered on.

- Open Elsewhere

A camera can only be controlled by a single instance of a program. If a camera has already been opened by another program (i.e., it is open elsewhere,) it is unavailable and cannot be used until it is closed.

Information contained in the camera's firmware can be read if the specified camera is connected. The information identifies the sensor and camera model and provides the logic program IDs and revision levels. This information may not be available for cameras that have been opened elsewhere (in another process).

4.1 Data Type Definitions

Refer to [Table 4-1](#) for information about data definitions specific to Camera APIs.

Table 4-1: Data Definitions for Camera APIs

Name	Type	Description
PicamModel	enum	The camera model. Series models represent a model family and may be used to represent older cameras whose exact model is not known.
PicamComputerInterface	enum	The interface used to communicate with the camera.
PicamStringSize	enum	Fixed sizes limiting the maximum size of some <code>picam</code> strings.
PicamHandle	void*	A <code>PICam</code> allocated resource.

4.2 Structure Definitions

This section provides information about structures required by the camera APIs.

4.2.1 PicamCameraID

Structure Definition

The structure definition for `PicamCameraID` is:

```
typedef struct PicamCameraID
{
    PicamModel    model;
    PicamComputerInterface  computer_interface;
    pichar    sensor_name [ ];
    pichar    serial_number [ ];
} PicamCameraID;
```

Variable Definitions

The variables required for `PicamCameraID` are:

model: This is the camera model.

computer_interface: This is the method by which the camera communicates with the host computer.

sensor_name: This is the name of the sensor in the camera.

serial_number: This is the unique serial number that corresponds with the camera.

4.2.2 PicamFirmwareDetail

Structure Definition

The structure definition for `PicamFirmwareDetail` is:

```
typedef struct PicamFirmwareDetail
{
    pichar    name [ ];
    pichar    detail [ ];
} PicamFirmwareDetail;
```

Variable Definitions

The variables required for `PicamFirmwareDetail` are:

name: This is the name of a hardware device containing firmware.

detail: This stores information about the hardware device, such as version number.

4.3 Programmers' Reference for Camera APIs

This section provides a detailed programmers' reference guide for the following APIs:

- Camera Identification APIs
 - `Picam_DestroyCameraIDs()`
 - `Picam_GetAvailableCameraIDs()`
 - `Picam_GetUnavailableCameraIDs()`
 - `Picam_IsCameraIDConnected()`
 - `Picam_IsCameraIDOpenElsewhere()`
- Camera Access APIs
 - `Picam_DestroyHandles()`
 - `Picam_OpenFirstCamera()`
 - `Picam_OpenCamera()`
 - `Picam_CloseCamera()`
 - `Picam_GetOpenCameras()`
 - `Picam_IsCameraConnected()`
 - `Picam_IsCameraFaulted()`
 - `Picam_GetCameraID()`
- Camera Information APIs
 - `Picam_DestroyFirmwareDetails()`
 - `Picam_GetFirmwareDetails()`
- Demo Camera Identification APIs
 - `Picam_DestroyModels()`
 - `Picam_GetAvailableDemoCameraModels()`
 - `Picam_ConnectDemoCamera()`
 - `Picam_DisconnectDemoCamera()`
 - `Picam_IsDemoCamera()`

4.3.1 Camera Identification APIs

This section provides programming information for Camera Identification APIs.

4.3.1.1 *Picam_DestroyCameraIDs()*

Description

`Picam_DestroyCameraIDs()` releases PICam-allotted memory associated with `id_array`.



NOTE:

`id_array` may be a single `PicamCameraID` allocated by PICam.

If `id_array` is a null array, calling `Picam_DestroyCameraIDs()` has no effect.

Syntax

The syntax for `Picam_DestroyCameraIDs()` is:

```
PICAM_API Picam_DestroyCameraIDs (  
    const PicamCameraID* id_array);
```

Input Parameters

Input parameters for `Picam_DestroyCameraIDs()` are:

`id_array`: Pointer to the `id_array` for which memory is to be released.

Output Parameters

There are no output parameters associated with `Picam_DestroyCameraIDs()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_GetAvailableCameraIDs()`;
- `Picam_GetUnavailableCameraIDs()`.

4.3.1.2 *Picam_GetAvailableCameraIDs()*

Description

`Picam_GetAvailableCameraIDs()` dynamically creates an array of length N. This array stores camera IDs for all available cameras.

**NOTE:**

Cameras that have been disconnected or are currently open in another process are not available.

**NOTE:**

Prior to program termination, memory that has been dynamically allocated to `id_array` must be released by calling `Picam_DestroyCameraIDs()`.

Syntax

The syntax for `Picam_GetAvailableCameraIDs()` is:

```
PICAM_API Picam_GetAvailableCameraIDs(  
    const PicamCameraID** id_array,  
    piint* id_count);
```

Input Parameters

There are no input parameters associated with `Picam_GetAvailableCameraIDs()`.

Output Parameters

Output parameters for `Picam_GetAvailableCameraIDs()` are:

- `id_array`: Pointer to the memory address for the array in which the list of available camera IDs is stored.
When there are no available camera IDs, a null value is returned.
- `id_count`: The total number of available camera IDs stored in `id_array`. This equals the length of the array that has been created.
When there are no available camera IDs, a value of 0 [zero] is returned.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_DestroyCameraIDs()`.

4.3.1.3 *Picam_GetUnavailableCameraIDs()*

Description

Picam_GetUnavailableCameraIDs() dynamically creates an array of length N. This array stores camera IDs for all unavailable cameras.

**NOTE:**

Cameras that have been disconnected or are currently open in another process are not available.

**NOTE:**

Prior to program termination, memory that has been dynamically allocated to *id_array* must be released by calling *Picam_DestroyCameraIDs()*.

Syntax

The syntax for *Picam_GetAvailableCameraIDs()* is:

```
PICAM_API Picam_GetUnavailableCameraIDs(  
    const PicamCameraID** id_array,  
    piint* id_count);
```

Input Parameters

There are no input parameters associated with *Picam_GetUnavailableCameraIDs()*.

Output Parameters

Output parameters for *Picam_GetUnavailableCameraIDs()* are:

- id_array*: Pointer to the memory address for the array in which the list of unavailable camera IDs is stored.
When there are no unavailable camera IDs, a null value is returned.
- id_count*: The total number of unavailable camera IDs stored in *id_array*. This equals the length of the array that has been created.
When there are no unavailable camera IDs, a value of 0 [zero] is returned.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_DestroyCameraIDs()*.

4.3.1.4 *Picam_IsCameraIDConnected()*

Description

`Picam_IsCameraIDConnected()` determines if a specified camera ID is plugged into the host computer and turned on.

Syntax

The syntax for `Picam_IsCameraIDConnected()` is:

```
PICAM_API Picam_IsCameraIDConnected(  
    const PicamCameraID* id,  
    pibln* connected);
```

Input Parameters

Input parameters for `Picam_IsCameraIDConnected()` are:

`id`: Specifies the ID of the camera for which the connection status is being tested.

Output Parameters

Output parameters for `Picam_IsCameraIDConnected()` are:

`connected`: Returns the connection status for the specified camera ID.

Valid values are:

- `True`
Indicates that the specified camera ID is connected to the host computer and is turned on;
- `False`
Indicates that the specified camera ID is not connected to the host computer or is not turned on.

4.3.1.5 *Picam_IsCameraIDOpenElsewhere()*

Description

`Picam_IsCameraIDOpenElsewhere()` determines if a specified camera ID has been opened by another process.

Syntax

The syntax for `Picam_IsCameraIDOpenElsewhere()` is:

```
PICAM_API Picam_IsCameraIDOpenElsewhere(  
    const PicamCameraID* id,  
    pibln* open_elsewhere);
```

Input Parameters

Input parameters for `Picam_IsCameraIDOpenElsewhere()` are:

`id`: Specifies the ID of the camera for which the connection status is being tested.

Output Parameters

Output parameters for `Picam_IsCameraIDOpenElsewhere()` are:

`open_elsewhere`: Returns the connection status for the specified camera ID.

Valid values are:

- `True`
Indicates that the specified camera ID is currently open in another process;
- `False`
Indicates that the specified camera ID is not currently open in another process.

4.3.2 Camera Access APIs

This section provides programming information for Camera Access APIs.

4.3.2.1 *Picam_DestroyHandles()*

Description

`Picam_DestroyHandles()` releases memory that has been allocated by PICam for use by `handle_array`.

**NOTE:**

`handle_array` may be a single `PicamHandle` allocated by PICam.

If `handle_array` is a null array, calling `Picam_DestroyHandles()` has no effect.

**NOTE:**

`Picam_DestroyHandles()` releases the memory used to store the handles. It does NOT free the resources to which the handles refer.

Syntax

The syntax for `Picam_DestroyHandles()` is:

```
PICAM_API Picam_DestroyHandles(  
    const PicamHandle* handle_array);
```

Input Parameters

Input parameters for `Picam_DestroyHandles()` are:

`handle_array`: Pointer to array memory that is to be released.

Output Parameters

There are no output parameters associated with `Picam_DestroyHandles()`.

4.3.2.2 *Picam_OpenFirstCamera()*

Description

Picam_OpenFirstCamera() opens the first available camera, and returns a handle to the camera.



NOTE:

Prior to program termination, all open cameras must be closed by calling *Picam_CloseCamera()*.

Syntax

The syntax for *Picam_OpenFirstCamera()* is:

```
PICAM_API Picam_OpenFirstCamera(  
                                PicamHandle* camera);
```

Input Parameters

There are no input parameters associated with *Picam_OpenFirstCamera()*.

Output Parameters

Output parameters for *Picam_OpenFirstCamera()* are:

camera: The handle corresponding to the camera that has been opened.

Advanced API Usage

When used in conjunction with Advanced APIs, the handle returned is for the model.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_CloseCamera()*.

4.3.2.3 *Picam_OpenCamera()*

Description

`Picam_OpenCamera()` opens a specified camera, and returns a handle to the camera.



NOTE:

Prior to program termination, all open cameras must be closed by calling `Picam_CloseCamera()`.

Syntax

The syntax for `Picam_OpenCamera()` is:

```
PICAM_API Picam_OpenCamera(  
    const PicamCameraID* id,  
    PicamHandle* camera);
```

Input Parameters

Input parameters for `Picam_OpenCamera()` are:

`id`: The `id` for camera to be opened.

Output Parameters

Output parameters for `Picam_OpenCamera()` are:

`camera`: The handle corresponding to the open camera.

Advanced API Usage

When used in conjunction with Advanced APIs, the handle returned is for the camera model.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_CloseCamera()`.

4.3.2.4 *Picam_CloseCamera()*

Description

`Picam_CloseCamera()` releases all resources that have been associated with a specified camera.

Syntax

The syntax for `Picam_CloseCamera()` is:

```
PICAM_API Picam_CloseCamera(  
    PicamHandle camera);
```

Input Parameters

Input parameters for `Picam_CloseCamera()` are:

`camera`: The handle associated with the camera that is to be closed.

Output Parameters

There are no output parameters associated with `Picam_CloseCamera()`.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` can be a handle to either the:

- `device`, or
- `model`.

In either case, when `Picam_CloseCamera()` is called, it always closes both the specified device and model.

4.3.2.5 *Picam_GetOpenCameras()*

Description

`Picam_GetOpenCameras()` dynamically creates an array of length N. This array stores camera handles for all open cameras in the current process.



NOTE:

Prior to program termination, memory that has been dynamically allocated to `camera_array` must be released by calling `Picam_DestroyHandles()`.

Syntax

The syntax for `Picam_GetOpenCameras()` is:

```
PICAM_API Picam_GetOpenCameras(  
    const PicamHandle** camera_array,  
    piint* camera_count);
```

Input Parameters

There are no input parameters associated with `Picam_GetOpenCameras()`.

Output Parameters

Output parameters for `Picam_GetOpenCameras()` are:

- `camera_array`: Pointer to the memory address for the array in which the list of camera handles is stored.
When there are no available camera handles, a null value is returned.
- `camera_count`: The total number of camera handles stored in `camera_array`. This equals the length of the array that has been created.
When there are no available camera handles, a value of 0 [zero] is returned.

Advanced API Usage

When used in conjunction with Advanced APIs, this array (`camera_array`) stores a list of model handles.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_DestroyHandles()`.

4.3.2.6 *Picam_IsCameraConnected()*

Description

`Picam_IsCameraConnected()` determines if the specified camera is plugged into the host computer and is turned on.

Syntax

The syntax for `Picam_IsCameraConnected()` is:

```
PICAM_API Picam_IsCameraConnected(  
                                PicamHandle camera,  
                                pibln* connected);
```

Input Parameters

Input parameters for `Picam_IsCameraConnected()` are:

`camera`: The handle for the camera for which the status is being determined.

Output Parameters

Output parameters for `Picam_IsCameraConnected()` are:

`connected`: Returns the connection status for the specified camera.

Valid values are:

- `True`
Indicates that the specified camera is connected to the host computer and is turned on.
- `False`
Indicates that the specified camera is not connected to the host computer and/or not turned on.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` can be a handle to either the:

- `device`, or
- `model`.

Both `device` and `model` share the same connected state.

4.3.2.7 *Picam_IsCameraFaulted()*

Description

`Picam_IsCameraFaulted()` determines if the specified camera has experienced a critical malfunction and is in need of service. Any acquisition in progress will be stopped and further acquisitions are not possible until the camera has been serviced.

Syntax

The syntax for `Picam_IsCameraFaulted()` is:

```
PICAM_API Picam_IsCameraFaulted(  
                                PicamHandle camera,  
                                pibln*   faulted);
```

Input Parameters

Input parameters for `Picam_IsCameraFaulted()` are:

`camera`: The handle for the camera for which the status is being determined.

Output Parameters

Output parameters for `Picam_IsCameraFaulted()` are:

`faulted`: Returns the faulted status for the specified camera.

Valid values are:

- `True`
Indicates that the specified camera has experienced a critical malfunction.
- `False`
Indicates that the specified camera is working properly.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` can be a handle to either the:

- `device`, or
- `model`.

Both `device` and `model` share the same faulted state.

4.3.2.8 *Picam_GetCameraID()*

Description

`Picam_GetCameraID()` returns the ID associates with a specified camera handle.

Syntax

The syntax for `Picam_GetCameraID()` is:

```
PICAM_API Picam_GetCameraID(  
    PicamHandle camera,  
    PicamCameraID* id);
```

Input Parameters

Input parameters for `Picam_GetCameraID()` are:

`camera`: The handle associated with the camera for which the ID is to be determined.

Output Parameters

Output parameters for `Picam_GetCameraID()` are:

`id`: The camera ID associated with the specified handle.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` can be a handle to either the:

- `device`, or
- `model`.

Both `device` and `model` share the same camera ID.

4.3.3 Camera Information APIs

This section provides programming information for Camera Information APIs.

4.3.3.1 *Picam_DestroyFirmwareDetails()*

Description

`Picam_DestroyFirmwareDetails()` releases memory that has been allocated for use by the `firmware_array`.



NOTE:

`firmware_array` may be a single `PicamFirmwareDetail` allocated by `PICam`.

If `firmware_array` is a null array, calling `Picam_DestroyFirmwareDetails()` has no effect.

Syntax

The syntax for `Picam_DestroyFirmwareDetails()` is:

```
PICAM_API Picam_DestroyFirmwareDetails(  
    const PicamFirmwareDetail* firmware_array);
```

Input Parameters

Input parameters for `Picam_DestroyFirmwareDetails()` are:

`firmware_array`: Pointer to the memory location where the array is stored.

Output Parameters

There are no output parameters associated with `Picam_DestroyFirmwareDetails()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_GetFirmwareDetails()`

4.3.3.2 *Picam_GetFirmwareDetails()*

Description

`Picam_GetFirmwareDetails()` dynamically creates an array of length N. This array stores firmware details associated with a specified camera ID.



NOTE:

Prior to program termination, memory that has been dynamically allocated to `firmware_array` must be released by calling `Picam_DestroyFirmwareDetails()`.

Syntax

The syntax for `Picam_GetFirmwareDetails()` is:

```
PICAM_API Picam_GetFirmwareDetails(  
    const PicamCameraID* id,  
    const PicamFirmwareDetail** firmware_array,  
    piint* firmware_count);
```

Input Parameters

Input parameters for `Picam_GetFirmwareDetails()` are:

`id`: Camera id for which firmware details are to be retrieved.

Output Parameters

Output parameters for `Picam_GetFirmwareDetails()` are:

`firmware_array`: Pointer to the memory address for the array in which firmware information is stored.
When no information is stored, a null value is returned.

`firmware_count`: The total number of firmware details stored in `firmware_array`. This equals the length of the array that has been created.
When no information is available, a value of 0 [zero] is returned.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_DestroyFirmwareDetails()`.

4.3.4 Demo Camera Identification APIs

This section provides programming information for Demo Camera Identification APIs.

4.3.4.1 *Picam_DestroyModels()*

`Picam_DestroyModels()` releases memory that has been allocated for use by the `model_array`.

**NOTE:**

`model_array` may be a single `PicamModel` allocated by PICam.

If `model_array` is a null array, calling `Picam_DestroyModels()` has no effect.

Syntax

The syntax for `Picam_DestroyModels()` is:

```
PICAM_API Picam_DestroyModels(  
    const PicamModel* model_array);
```

Input Parameters

Input parameters for `Picam_DestroyModels()` are:

`model_array`: Pointer to the memory location where the array is stored.

Output Parameters

There are no output parameters associated with `Picam_DestroyModels()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_GetAvailableDemoCameraModels()`

4.3.4.2 *Picam_GetAvailableDemoCameraModels()*

Description

`Picam_GetAvailableDemoCameraModels()` dynamically creates an array of length N. This array stores a list of virtual camera models which are available for use in Demo Mode.



NOTE:

Prior to program termination, memory that has been dynamically allocated to `model_array` must be released by calling `Picam_DestroyModels()`.

Syntax

The syntax for `Picam_GetAvailableDemoCameraModels()` is:

```
PICAM_API Picam_GetAvailableDemoCameraModels(  
    const PicamModel** model_array,  
    piint* model_count);
```

Input Parameters

There are no input parameters associated with `Picam_GetAvailableDemoCameraModels()`.

Output Parameters

Output parameters for `Picam_GetAvailableDemoCameraModels()` are:

- `model_array`: Pointer to the memory address for the array in which the list of virtual camera models is stored.
When there are no virtual camera models available, a null value is returned.
- `model_count`: The total number of virtual models being stored in `model_array`. This equals the length of the array that has been created.
When there are no virtual models available, a value of 0 [zero] is returned.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_DestroyModels()`.

4.3.4.3 *Picam_ConnectDemoCamera()*

Description

`Picam_ConnectDemoCamera()` establishes a connection with the specified virtual camera.

Syntax

The syntax for `Picam_ConnectDemoCamera()` is:

```
PICAM_API Picam_ConnectDemoCamera(  
    PicamModel model,  
    const pichar* serial_number,  
    PicamCameraID* id);
```

Input Parameters

Input parameters for `Picam_ConnectDemoCamera()` are:

`model`: Model for the virtual camera for which a connection is to be established.

`serial_number`: Serial number of the virtual camera for which a connection is to be established.

Output Parameters

Output parameters for `Picam_ConnectDemoCamera()` are:

`id`: ID of the virtual camera for which a connection is to be established



NOTE:

`id` is an optional parameter and may be null.

4.3.4.4 *Picam_DisconnectDemoCamera()*

Description

`Picam_DisconnectDemoCamera()` breaks an established connection with the specified virtual camera.

Syntax

The syntax for `Picam_DisconnectDemoCamera()` is:

```
PICAM_API Picam_DisconnectDemoCamera(  
    const PicamCameraID* id);
```

Input Parameters

Input parameters for `Picam_DisconnectDemoCamera()` are:

`id`: ID of the virtual camera for which the connection is to be broken.

Output Parameters

There are no output parameters associated with `Picam_DisconnectDemoCamera()`.

4.3.4.5 *Picam_IsDemoCamera()*

Description

`Picam_IsDemoCamera()` determines if the specified camera is a virtual camera.

Syntax

The syntax for `Picam_IsDemoCamera()` is:

```
PICAM_API Picam_IsDemoCamera(  
    const PicamCameraID* id,  
    pibln* demo);
```

Input Parameters

Input parameters for `Picam_IsDemoCamera()` are:

`id`: ID of the camera being identified.

Output Parameters

Output parameters for `Picam_IsDemoCamera()` are:

`demo`: Indicates if the specified camera is a software-simulated camera.

Valid values are:

- `True`
Indicates that the specified camera is a virtual camera.
- `False`
Indicates that the specified camera is an actual physical camera.

This page is intentionally blank.

Chapter 5: Camera Configuration APIs

The functions in this grouping set or query parameter values, parameter information, and parameter constraints that characterize a camera. A parameter is a camera setting. Parameters have varying qualities as well as values and constraints. A parameter may have several different values but constraints determine which kinds of values a parameter can have based on camera type, read/write capability, or other parameters used or to be used in describing and setting up a particular camera. Once a camera has been opened, it can be configured by changing its parameters through software and applying them to the camera hardware. Each camera model has a different set of parameters. Parameters contain different attributes.

The most important parameter attribute is its value. Values are represented by different types (i.e., integer, floating point, enumeration, etc.)

All parameter values can be read, but not all can be written. This is determined by the parameter's value access:

- Read/Write
- Read Only.



NOTE:

Note a special case of value access is when a parameter value can be written, but only one particular value is permitted. This is called read/write trivial.

Parameter values that can be written have constraints. Constraints describe the set of values a parameter value can take. The nature of this set determines the constraint type (e.g., a numeric range, a set of options, etc.) It is useful to describe a different constraint based on purpose. This is where constraint categories come into play. These categories differentiate “is this parameter capable of x ?” from “based on my current configuration, is it valid to set parameter to x ?”

Due to the complex nature of configuration, some parameters override others when certain values are set. A parameter is relevant if it has an effect on the current configuration.

Most parameters are only used for acquisition setup. However, this is not always the case. Some can be modified while the camera is acquiring. These parameters are deemed onlineable. Others reflect the current state of the hardware. These parameters only have meaning when read directly from hardware. These are termed readable.

Another parameter may reflect the status of hardware that is not directly controllable by the software (e.g., may be changed due to external influences,) yet its value may impact the decisions and/or further progress of the software. Such a parameter is a waitable status.

Once the parameters values are adjusted as desired they must be committed to the camera hardware before the camera can acquire data.

5.1 Data Type Definitions

This section provides programming information about PICam data definitions.

5.1.1 Camera Parameter Enumerations

This section provides detailed information about the following Camera Parameter enumerations:

- [PicamValueType](#)
- [PicamConstraintType](#)
- [PicamParameter](#)
- [PicamActiveShutter](#)
- [PicamAdcAnalogGain](#)
- [PicamAdcQuality](#)
- [PicamCcdCharacteristicsMask](#)
- [PicamCoolingFanStatus](#)
- [PicamEMICcdGainControlMode](#)
- [PicamGateTrackingMask](#)
- [PicamGatingMode](#)
- [PicamGatingSpeed](#)
- [PicamIntensifierOptionsMask](#)
- [PicamIntensifierStatus](#)
- [PicamModulationTrackingMask](#)
- [PicamOrientationMask](#)
- [PicamOutputSignal](#)
- [PicamPhosphorType](#)
- [PicamPhotocathodeSensitivity](#)
- [PicamPhotonDetectionMode](#)
- [PicamPixelFormat](#)
- [PicamReadoutControlMode](#)
- [PicamSensorTemperatureStatus](#)
- [PicamSensorType](#)
- [PicamShutterStatus](#)
- [PicamShutterTimingMode](#)
- [PicamShutterType](#)
- [PicamTimeStampsMask](#)
- [PicamTriggerCoupling](#)
- [PicamTriggerDetermination](#)
- [PicamTriggerResponse](#)
- [PicamTriggerSource](#)
- [PicamTriggerTermination](#)

5.1.1.1 *PicamValueType*

Data Type

PicamValueType is defined as enum.

Description

PicamValueType is the set of parameter value data types.

Enumerator Definitions

Refer to [Table 5-1](#) for enumerator definitions.

Table 5-1: *PicamValueType* Enumerator Definitions

Enumerator	Description
<i>PicamValueType_Boolean</i>	Accessed as <i>piint</i> . <ul style="list-style-type: none">• FALSE = 0• TRUE = 1
<i>PicamValueType_Enumeration</i>	Any enum accessed as <i>piint</i> .
<i>PicamValueType_FloatingPoint</i>	Accessed as <i>piflt</i> .
<i>PicamValueType_Integer</i>	Accessed as <i>piint</i> .
<i>PicamValueType_LargeInteger</i>	Accessed as <i>pi64s</i> .
<i>PicamValueType_Modulations</i>	Accessed as <i>PicamModulations</i>
<i>PicamValueType_Pulse</i>	Accessed as <i>PicamPulse</i> .
<i>PicamValueType_Rois</i>	Accessed as <i>PicamRois</i> .

5.1.1.2 *PicamConstraintType*

Data Type

PicamConstraintType is defined as enum.

Description

PicamConstraintType is the set of constraints that may be placed on a parameter's value.

Enumerator Definitions

Refer to [Table 5-2](#) for enumerator definitions.

Table 5-2: *PicamConstraintType* Enumerator Definitions

Enumerator	Description
<i>PicamConstraintType_Collection</i>	The value can be one in a collection of choices.
<i>PicamConstraintType_Modulations</i>	The value is a custom modulation sequence.
<i>PicamConstraintType_None</i>	The value is read-only and not constrained.
<i>PicamConstraintType_Pulse</i>	The value is a gate pulse.
<i>PicamConstraintType_Range</i>	The value is numeric and naturally constrained within a linear range.
<i>PicamConstraintType_Rois</i>	The value is a set of regions of interests.

5.1.1.3 *PicamParameter*

Data Type

PicamParameter is defined as enum.

Description

PicamParameter is the set of user-accessible camera parameters.

Enumerator Definitions

Refer to [Table 5-3](#) for enumerator definitions.

Table 5-3: *PicamParameter* Enumerator Definitions (Sheet 1 of 9)

Enumerator	Description
<i>PicamParameter_Accumulations</i>	Controls the number of on-sensor accumulations
<i>PicamParameter_ActiveBottomMargin</i>	Controls the inactive number of rows on the bottom.
<i>PicamParameter_ActiveExtendedHeight</i>	Controls the number of additional active rows that can be used for storage. NOTE: These rows cannot be imaged directly.
<i>PicamParameter_ActiveHeight</i>	Controls the active number of rows.
<i>PicamParameter_ActiveLeftMargin</i>	Controls the inactive number of columns on the left.
<i>PicamParameter_ActiveRightMargin</i>	Controls the inactive number of columns on the right.
<i>PicamParameter_ActiveShutter</i>	Selects the shutter via the <i>PicamActiveShutter</i> data enumeration. Refer to Section 5.1.1.4, <i>PicamActiveShutter</i> , on page 62 for additional information.
<i>PicamParameter_ActiveTopMargin</i>	Controls the inactive number of rows on the top.
<i>PicamParameter_ActiveWidth</i>	Controls the active number of columns.
<i>PicamParameter_AdcAnalogGain</i>	Controls the electronic gain of the pixel digitization via the <i>PicamAdcAnalogGain</i> data enumeration. Refer to Section 5.1.1.5, <i>PicamAdcAnalogGain</i> , on page 62 for additional information.
<i>PicamParameter_AdcBitDepth</i>	Controls the resolution of the pixel digitization in bits-per-pixel.
<i>PicamParameter_AdcEMGain</i>	Controls the electromagnetic gain in terms of multiples.
<i>PicamParameter_AdcQuality</i>	Controls the nature of pixel digitization via the <i>PicamAdcQuality</i> data enumeration. Refer to Section 5.1.1.6, <i>PicamAdcQuality</i> , on page 63 for additional information.
<i>PicamParameter_AdcSpeed</i>	Controls the rate pixels are digitized in MHz.
<i>PicamParameter_AnticiapteTrigger</i>	Uses an external pre-trigger to anticipate an external trigger.
<i>PicamParameter_AuxOutput</i>	Controls the auxiliary output gate pulse.
<i>PicamParameter_BracketGating</i>	Enables bracket pulsing for an intensified camera.

Table 5-3: PicamParameter Enumerator Definitions (Sheet 2 of 9)

Enumerator	Description
PicamParameter_CcdCharacteristics	Reports characteristics of a CCD sensor via the PicamCcdCharacteristicsMask data enumeration. Refer to Section 5.1.1.7, PicamCcdCharacteristicsMask , on page 64 for additional information.
PicamParameter_CleanBeforeExposure	Controls cleaning before each exposure.
PicamParameter_CleanCycleCount	Controls the number of clean cycles to run before acquisition begins.
PicamParameter_CleanCycleHeight	Controls the number of rows in a clean cycle.
PicamParameter_CleanSectionFinalHeight	Controls the final height rows for exponential decomposition cleaning.
PicamParameter_CleanSectionFinalHeightCount	Controls the final height iterations for exponential decomposition cleaning.
PicamParameter_CleanSerialRegister	Controls the cleaning of the serial register itself.
PicamParameter_CleanUntilTrigger	Controls the nature of cleaning while waiting for an external trigger.
PicamParameter_CoolingFanStatus	Reports the status of the cooling fan via the PicamCoolingFanStatus data enumeration. Refer to Section 5.1.1.8, PicamCoolingFanStatus , on page 65 for additional information.
PicamParameter_CorrectPixelBias	Enables pixel bias correction.
PicamParameter_CustomModulationSequence	Customizes a modulation sequence.
PicamParameter_DelayFromPreTrigger	Specifies the delay from pre-trigger to trigger in microseconds (μ S).
PicamParameter_DifEndingGate	Controls the second gate pulse in DIF readout in nanoseconds (nS).
PicamParameter_DifStartingGate	Controls the initial gate pulse in DIF readout in nanoseconds (nS).
PicamParameter_DisableCoolingFan	Enables/disables the thermoelectric cooling fan.
PicamParameter_DisableDataFormatting	Controls the basic processing necessary to receive data in the expected format.
PicamParameter_EMIccdGain	Optimally controls the intensifier gain and electron multiplication gain in an eMCCD camera in terms of multiples.
PicamParameter_EMIccdGainControlMode	Determines how the intensifier gain and electron multiplication gain are controlled in an eMCCD camera via the PicamEMIccdGainControlMode data enumeration. Refer to Section 5.1.1.9, PicamEMIccdGainControlMode , on page 65, for additional information.
PicamParameter_EnableIntensifier	Enables the intensifier. NOTE: The intensifier must be enabled and powered on for it to function.
PicamParameter_EnableModulation	Enables RF modulation for an intensified camera.

Table 5-3: PicamParameter Enumerator Definitions (Sheet 3 of 9)

Enumerator	Description
PicamParameter_EnableModulationOutputSignal	Enables an RF output signal from the intensified camera to be used as the user sees fit.
PicamParameter_EnableNondestructiveReadout	Allows the camera to periodically readout while exposing.
PicamParameter_EnableSensorWindowHeater	Enables the sensor window to heat up in an effort to prevent condensation.
PicamParameter_EnableSyncMaster	Enables SyncMASTER1 and SyncMASTER2 gate pulses.
PicamParameter_ExactReadoutCountMaximum	Reports the maximum number of readouts the camera can acquire. NOTE: This does not include non-destructive readouts from cameras that support such a feature.
PicamParameter_ExposureTime	Controls the time the sensor is exposed in milliseconds (mS).
PicamParameter_ExternalShutterStatus	Reports the status of the shutter that is external to the camera via the PicamShutterStatus data enumeration. Refer to Section 5.1.1.25, PicamShutterStatus , on page 74 for additional information.
PicamParameter_ExternalShutterType	Reports the type of shutter that is external to and can be driven by the camera via the PicamShutterType data enumeration. Refer to Section 5.1.1.27, PicamShutterType , on page 75 for additional information.
PicamParameter_FrameRateCalculation	Reports the estimated frame rate in frames-per-second. NOTE: If there is more than one frame-per-readout, this represents the burst frame rate within the readout. NOTE: If the camera is being externally triggered, this represents the fastest possible rate.
PicamParameter_FrameSize	Reports the size, in bytes, of a data frame.
PicamParameter_FramesPerReadout	Reports the number of frames contained in one readout.
PicamParameter_FrameStride	Reports the length, in bytes, necessary to traverse to the next frame.
PicamParameter_FrameTrackingBitDepth	Controls the frame tracking number size in bits-per-pixel.
PicamParameter_GateTracking	Controls the tracking of a sequential gate pulse in metadata via the PicamGateTrackingMask data enumeration. Refer to Section 5.1.1.10, PicamGateTrackingMask , on page 66 for additional information.
PicamParameter_GateTrackingBitDepth	Controls the size of one component in a varying sequential gate pulse. NOTE: This metadata is floating point.
PicamParameter_GatingMode	Controls the nature of gate pulse timing via the PicamGatingMode data enumeration. Refer to Section 5.1.1.11, PicamGatingMode , on page 66 for additional information.
PicamParameter_GatingSpeed	Classifies the narrowest gate pulse.

Table 5-3: PicamParameter Enumerator Definitions (Sheet 4 of 9)

Enumerator	Description
PicamParameter_InactiveShutterTimingMode Result	Reports the state of the inactive shutter via the PicamShutterTimingMode data enumeration. Refer to Section 5.1.1.26, PicamShutterTimingMode , on page 75 for additional information.
PicamParameter_IntensifierDiameter	Reports the diameter of the intensifier in millimeters (mm).
PicamParameter_IntensifierGain	Controls the gain of the intensifier in terms of multiples.
PicamParameter_IntensifierOptions	Reports additional features of the intensified camera via the PicamIntensifierOptionsMask data enumeration. Refer to Section 5.1.1.13, PicamIntensifierOptionsMask , on page 67 for additional information.
PicamParameter_IntensifierStatus	Reports the status of the intensifier power via the PicamIntensifierStatus data enumeration. Refer to Section 5.1.1.14, PicamIntensifierStatus , on page 68 for additional information.
PicamParameter_InternalShutterStatus	Reports the status of the shutter that is internal to the camera via the PicamShutterStatus data enumeration. Refer to Section 5.1.1.25, PicamShutterStatus , on page 74 for additional information.
PicamParameter_InternalShutterType	Reports the type of shutter that is internal to and can be driven by the camera via the PicamShutterType data enumeration. Refer to Section 5.1.1.27, PicamShutterType , on page 75 for additional information.
PicamParameter_InvertOutputSignal	Controls if the timing signal is inverted when viewed from the camera monitor.
PicamParameter_InvertOutputSignal2	Controls if the timing signal is inverted when viewed from the second camera monitor.
PicamParameter_KineticsWindowHeight	Controls the number of rows used for the sensing window in a kinetics readout.
PicamParameter_MaskedBottomMargin	Controls the number of masked rows akin to active bottom margin.
PicamParameter_MaskedHeight	Controls the number of masked rows akin to active height.
PicamParameter_MaskedTopMargin	Controls the number of masked rows akin to active top margin.
PicamParameter_ModulationDuration	Controls the time the intensifier is modulating in milliseconds (mS).
PicamParameter_ModulationFrequency	Controls the frequency of the intensifier modulation in MHz.
PicamParameter_ModulationOutputSignal Amplitude	Controls the peak-to-peak amplitude of the user RF output signal in volts (V).
PicamParameter_ModulationOutputSignal Frequency	Controls the frequency of the user RF output signal in MHz

Table 5-3: `PicamParameter` Enumerator Definitions (Sheet 5 of 9)

Enumerator	Description
<code>PicamParameter_ModulationTracking</code>	Controls the tracking of a sequential phase or custom modulation sequence in metadata via the PicamModulationTrackingMask data enumeration. Refer to Section 5.1.1.15, PicamModulationTrackingMask , on page 68 for additional information.
<code>PicamParameter_ModulationTrackingBitDepth</code>	Controls the size of one component in a varying sequential modulation phase or custom modulation sequence. NOTE: This metadata is floating point.
<code>PicamParameter_NondestructiveReadoutPeriod</code>	Controls the rate at which the camera will non-destructively readout during exposure in seconds (S). NOTE: This duration must be less than exposure time for any non-destructive readouts to occur.
<code>PicamParameter_NormalizeOrientation</code>	Controls automatic orientation correction for data due to readout ports used.
<code>PicamParameter_OnlineReadoutRateCalculation</code>	Reports the fastest possible readout rate that could occur given the current setup and accounting for possible changes to online camera parameters while acquiring
<code>PicamParameter_Orientation</code>	Reports the orientation of the data via the PicamOrientationMask data enumeration.. Refer to Section 5.1.1.16, PicamOrientationMask , on page 69 for additional information.
<code>PicamParameter_OutputSignal</code>	Controls what timing signal is issued from the camera monitor via the PicamOutputSignal data enumeration. Refer to Section 5.1.1.17, PicamOutputSignal , on page 69 for additional information.
<code>PicamParameter_OutputSignal2</code>	Controls what timing signal is issued from the second camera monitor via the PicamOutputSignal data enumeration. Refer to Section 5.1.1.17, PicamOutputSignal , on page 69 for additional information.
<code>PicamParameter_PhosphorDecayDelay</code>	Controls the length of time a camera waits for the phosphor to decay before reading out. The time unit depends on PicamParameter_PhosphorDecayDelayResolution .
<code>PicamParameter_PhosphorDecayDelayResolution</code>	Controls the time unit used for phosphor decay delay. This value is in microseconds (μ S). Example: A resolution of: <ul style="list-style-type: none"> • 1 signifies delay is in microseconds (μS); • 1000 for milliseconds (mS); • 0.01 for tens-of-nanoseconds (nS).
<code>PicamParameter_PhosphorType</code>	Reports the type of phosphor used in the intensified camera.
<code>PicamParameter_PhotocathodeSensitivity</code>	Classifies the wavelength sensitivity of the photocathode.

Table 5-3: PicamParameter Enumerator Definitions (Sheet 6 of 9)

Enumerator	Description
PicamParameter_PhotonDetectionMode	Enables/disables photon detection and controls how it is done via the PicamPhotonDetectionMode data enumeration. Refer to Section 5.1.1.20, PicamPhotonDetectionMode , on page 72 for additional information.
PicamParameter_PhotonDetectionThreshold	The threshold, in counts, used to distinguish photons from background for each pixel.
PicamParameter_PixelBitDepth	Reports the size of a data pixel in bits-per-pixel.
PicamParameter_PixelFormat	Controls the format of a data pixel via the PicamPixelFormat data enumeration. Refer to Section 5.1.1.21, PicamPixelFormat , on page 72 for additional information.
PicamParameter_PixelGapHeight	Reports the vertical distance between pixels, in microns.
PicamParameter_PixelGapWidth	Reports the horizontal distance between pixels, in microns.
PicamParameter_PixelHeight	Reports the pixel height, in microns.
PicamParameter_PixelWidth	Reports the pixel width, in microns.
PicamParameter_ReadoutControlMode	Controls how the sensor is read out via the PicamReadoutControlMode data enumeration. Refer to Section 5.1.1.22, PicamReadoutControlMode , on page 73 for additional information.
PicamParameter_ReadoutCount	Controls the number of readouts to acquire before stopping the camera. NOTE: The camera may acquire more than the readouts requested for large requests (i.e., more readouts than those specified by PicamParameter_ExactReadoutCountMaximum). NOTE: This does not include non-destructive readouts from cameras that support such a feature. NOTE: [Advanced-API Usage Only] The value 0 indicates the camera will run forever until explicitly stopped or an error occurs.
PicamParameter_ReadoutOrientation	Reports the orientation of the data due to readout port location via the PicamOrientationMask data enumeration. Refer to Section 5.1.1.16, PicamOrientationMask , on page 69 for additional information.
PicamParameter_ReadoutPortCount	Controls the number of readout ports from which the camera should simultaneously read data.
PicamParameter_ReadoutRateCalculation	Reports the estimated rate of data in readouts-per-second. NOTE: If the camera is being externally triggered, this represents the fastest possible rate.
PicamParameter_ReadoutStride	Reports the length, in bytes, necessary to traverse to the next readout.
PicamParameter_ReadoutTimeCalculation	Reports the duration of time it takes for the camera to read out one readout in milliseconds (mS).

Table 5-3: `PicamParameter` Enumerator Definitions (Sheet 7 of 9)

Enumerator	Description
<code>PicamParameter_RepetitiveGate</code>	Controls the constant gate pulse in nanoseconds (nS).
<code>PicamParameter_RepetitiveModulationPhase</code>	Controls the constant phase of the intensifier with respect to the modulation output signal in degrees.
<code>PicamParameter_Rois</code>	Controls the area of the sensor to be digitized via the PicamRois structure. Refer to Section 5.2.1.2, PicamRois , on page 86 for additional information.
<code>PicamParameter_SecondaryActiveHeight</code>	Controls the number of secondary active rows.
<code>PicamParameter_SecondaryMaskedHeight</code>	Controls the number of secondary masked rows.
<code>PicamParameter_SensorActiveBottomMargin</code>	Reports the inactive rows on the bottom.
<code>PicamParameter_SensorActiveExtendedHeight</code>	Reports the number of additional active rows that can be used for storage. NOTE: These rows cannot be imaged directly.
<code>PicamParameter_SensorActiveHeight</code>	Reports the active number of rows.
<code>PicamParameter_SensorActiveLeftMargin</code>	Reports the inactive columns on the left.
<code>PicamParameter_SensorActiveRightMargin</code>	Reports the inactive columns on the right.
<code>PicamParameter_SensorActiveTopMargin</code>	Reports the inactive rows on the top.
<code>PicamParameter_SensorActiveWidth</code>	Reports the active number of columns.
<code>PicamParameter_SensorMaskedBottomMargin</code>	Reports the number of masked rows akin to active bottom margin.
<code>PicamParameter_SensorMaskedHeight</code>	Reports the number of masked rows akin to active height.
<code>PicamParameter_SensorMaskedTopMargin</code>	Reports the number of masked rows akin to active top margin.
<code>PicamParameter_SensorSecondaryActiveHeight</code>	Reports the number of secondary active rows.
<code>PicamParameter_SensorSecondaryMaskedHeight</code>	Reports the number of secondary masked rows.
<code>PicamParameter_SensorTemperatureReading</code>	Reports the temperature of the sensor.
<code>PicamParameter_SensorTemperatureSetPoint</code>	Controls the target temperature for the sensor.
<code>PicamParameter_SensorTemperatureStatus</code>	Reports the status of the sensor temperature via the PicamSensorTemperatureStatus data enumeration. Refer to Section 5.1.1.23, PicamSensorTemperatureStatus , on page 73 for additional information.
<code>PicamParameter_SensorType</code>	Reports the kind of sensor being used via the PicamSensorType data enumeration. Refer to Section 5.1.1.24, PicamSensorType , on page 74 for additional information.
<code>PicamParameter_SequentialEndingGate</code>	Controls the last gate pulse in a sequence in nanoseconds (nS).
<code>PicamParameter_SequentialEndingModulationPhase</code>	Controls the last modulation phase of the intensifier with respect to the modulation output signal in a sequence in degrees.
<code>PicamParameter_SequentialGateStepCount</code>	Controls the number of gate pulse steps in a sequence.

Table 5-3: PicamParameter Enumerator Definitions (Sheet 8 of 9)

Enumerator	Description
PicamParameter_SequentialGateStepIterations	Controls the number of gate pulses at each step in a sequence.
PicamParameter_SequentialStartingGate	Controls the first gate pulse in a sequence in nanoseconds (nS).
PicamParameter_SequentialStartingModulation Phase	Controls the first modulation phase of the intensifier with respect to the modulation output signal in a sequence in degrees.
PicamParameter_ShutterClosingDelay	Controls the duration of time the camera waits for the shutter to close before reading out. The time unit depends on PicamParameter_ShutterDelayResolution .
PicamParameter_ShutterDelayResolution	Controls the time unit used for shutter opening/closing delay. This value is in microseconds. Example: A resolution of: <ul style="list-style-type: none"> • 1 signifies delay is in microseconds (μS); • 1000 for milliseconds (mS); • 0.01 for tens-of-nanoseconds (nS).
PicamParameter_ShutterOpeningDelay	Controls the duration of time the camera waits for the shutter to open before exposing. The time unit depends on PicamParameter_ShutterDelayResolution .
PicamParameter_ShutterTimingMode	Controls the behavior of the shutter during acquisition.
PicamParameter_StopCleaningOnPreTrigger	Stops sensor cleaning when an external pre-trigger is acknowledged.
PicamParameter_SyncMaster2Delay	Controls the delay of SyncMASTER2 relative to SyncMASTER1 in microseconds (μS).
PicamParameter_TimeStampBitDepth	Controls the time stamp size in bits-per-pixel NOTE: Because time stamps may be negative one bit is reserved for sign.
PicamParameter_TimeStampResolution	Controls the time stamp resolution in ticks-per-second. NOTE: This value is computer-dependent when time stamps are software generated.
PicamParameter_TimeStamps	Controls time stamp metadata via the PicamTimeStampsMask data enumeration. Refer to Section 5.1.1.28, PicamTimeStampsMask , on page 76 for additional information.
PicamParameter_TrackFrames	Controls frame tracking metadata.
PicamParameter_TriggerCoupling	Controls the coupling between an external trigger source and the camera input via the PicamTriggerCoupling data enumeration. Refer to Section 5.1.1.29, PicamTriggerCoupling , on page 76 for additional information.

Table 5-3: `PicamParameter` Enumerator Definitions (Sheet 9 of 9)

Enumerator	Description
<code>PicamParameter_TriggerDetermination</code>	Controls what the camera recognizes as an external trigger via the PicamTriggerDetermination data enumeration. Refer to Section 5.1.1.30, PicamTriggerDetermination , on page 78 for additional information.
<code>PicamParameter_TriggerFrequency</code>	Controls the internal trigger and SyncMASTER frequency in Hz.
<code>PicamParameter_TriggerResponse</code>	Controls the camera's behavior in response to a trigger via the PicamTriggerResponse data enumeration. Refer to Section 5.1.1.31, PicamTriggerResponse , on page 79 for additional information.
<code>PicamParameter_TriggerResponse</code>	Controls the camera's behavior in response to an external trigger via the PicamTriggerResponse data enumeration. Refer to Section 5.1.1.31, PicamTriggerResponse , on page 79 for additional information.
<code>PicamParameter_TriggerSource</code>	Controls the source of a trigger via the PicamTriggerSource data enumeration. Refer to Section 5.1.1.32, PicamTriggerSource , on page 79 for additional information.
<code>PicamParameter_TriggerTermination</code>	Controls the termination of an external trigger source at the camera input via the PicamTriggerTermination data enumeration. Refer to Section 5.1.1.33, PicamTriggerTermination , on page 80 for additional information.
<code>PicamParameter_TriggerThreshold</code>	Controls the voltage threshold necessary for the camera to recognize a trigger in volts (V).
<code>PicamParameter_VerticalShiftRate</code>	Controls the rate to shift one row towards the serial register in a CCD in microseconds. (μ S)

5.1.1.4 *PicamActiveShutter*

Data Type

`PicamActiveShutter` is defined as enum.

Description

`PicamActiveShutter` is the shutter that will be controlled during an acquisition.

Enumerator Definitions

Refer to [Table 5-5](#) for enumerator definitions.

Table 5-4: `PicamActiveShutter` Enumerator Definitions

Enumerator	Description
<code>PicamActiveShutter_External</code>	The shutter external to the camera.
<code>PicamActiveShutter_Internal</code>	The shutter internal to the camera.
<code>PicamActiveShutter_None</code>	There is no shutter installed.

5.1.1.5 *PicamAdcAnalogGain*

Data Type

`PicamAdcAnalogGain` is defined as enum.

Description

`PicamAdcAnalogGain` is the set of electronic gain settings for pixel digitization.

Enumerator Definitions

Refer to [Table 5-5](#) for enumerator definitions.

Table 5-5: `PicamAdcAnalogGain` Enumerator Definitions

Enumerator	Description
<code>PicamAdcAnalogGain_High</code>	Large amplification. Refer to the user manual for the specific camera being used for complete information.
<code>PicamAdcAnalogGain_Low</code>	Small amplification. Refer to the user manual for the specific camera being used for complete information.
<code>PicamAdcAnalogGain_Medium</code>	Average amplification. Refer to the user manual for the specific camera being used for complete information.

5.1.1.6 *PicamAdcQuality*

Data Type

`PicamAdcQuality` is defined as enum.

Description

`PicamAdcQuality` is the set of Analog-to-Digital conversion techniques and quality settings for pixel digitization.

Enumerator Definitions

Refer to [Table 5-6](#) for enumerator definitions.

Table 5-6: `PicamAdcQuality` Enumerator Definitions

Enumerator	Description
<code>PicamAdcQuality_ElectronMultiplied</code>	Provides electron multiplication.
<code>PicamAdcQuality_HighCapacity</code>	Optimized for sensing high levels of radiation.
<code>PicamAdcQuality_HighSpeed</code>	Provides faster readout speeds.
<code>PicamAdcQuality_LowNoise</code>	Optimized for the lowest noise.

5.1.1.7 *PicamCcdCharacteristicsMask*

Data Type

`PicamCcdCharacteristicsMask` is defined as enum.

Description

`PicamCcdCharacteristicsMask` is the set of CCD sensor characteristics.

Enumerator Definitions

Refer to [Table 5-7](#) for enumerator definitions.

Table 5-7: `PicamCcdCharacteristicsMask` Enumerator Definitions

Enumerator	Description
<code>PicamCcdCharacteristicsMask_AdvancedInvertedMode</code>	The CCD has reduced dark current.
<code>PicamCcdCharacteristicsMask_BackIlluminated</code>	Indicates the type of illumination used. Valid values are: <ul style="list-style-type: none"> 1 CCD is back-illuminated 0 CCD is front-illuminated
<code>PicamCcdCharacteristicsMask_DeepDepleted</code>	The CCD is deep depleted.
<code>PicamCcdCharacteristicsMask_ExcelonEnabled</code>	The CCD is enhanced with eXcelon technology.
<code>PicamCcdCharacteristicsMask_HighResistivity</code>	The CCD is enhanced for sensing infrared radiation.
<code>PicamCcdCharacteristicsMask_Multiport</code>	The CCD has multiple readout ports that can be used simultaneously.
<code>PicamCcdCharacteristicsMask_None</code>	No additional characteristics.
<code>PicamCcdCharacteristicsMask_OpenElectrode</code>	The CCD is open electrode.
<code>PicamCcdCharacteristicsMask_SecondaryMask</code>	The CCD has an additional masked area.
<code>PicamCcdCharacteristicsMask_UVEnhanced</code>	The CCD is enhanced for sensing ultraviolet radiation.

5.1.1.8 *PicamCoolingFanStatus*

Data Type

`PicamCoolingFanStatus` is defined as enum.

Description

`PicamCoolingFanStatus` is the set of cooling fan statuses.

Enumerator Definitions

Refer to [Table 5-8](#) for enumerator definitions.

Table 5-8: `PicamCoolingFanStatus` Enumerator Definitions

Enumerator	Description
<code>PicamCoolingFanStatus_ForcedOn</code>	The cooling fan has been forced on to prevent overheating.
<code>PicamCoolingFanStatus_Off</code>	The cooling fan is off.
<code>PicamCoolingFanStatus_On</code>	The cooling fan is on.

5.1.1.9 *PicamEMICcdGainControlMode*

Data Type

`PicamEMICcdGainControlMode` is defined as enum.

Description

`PicamEMICcdGainControlMode` is the set of Control Modes which control intensifier gain and electron multiplication gain for an emICCD camera.

Enumerator Definitions

Refer to [Table 5-9](#) for enumerator definitions.

Table 5-9: `PicamEMICcdGainControlMode` Enumerator Definitions

Enumerator	Description
<code>PicamEMICcdGainControlMode_Manual</code>	Allows each gain to be controlled independently.
<code>PicamEMICcdGainControlMode_Optimal</code>	Controls both gains simultaneously as a single emICCD gain.

5.1.1.10 *PicamGateTrackingMask*

Data Type

PicamGateTrackingMask is defined as enum.

Description

PicamGateTrackingMask is the set of sequential gate pulse components that are to be tracked.

Enumerator Definitions

Refer to [Table 5-10](#) for enumerator definitions.

Table 5-10: *PicamGateTrackingMask* Enumerator Definitions

Enumerator	Description
<i>PicamGateTrackingMask_Delay</i>	The delay of the gate pulse is tracked.
<i>PicamGateTrackingMask_None</i>	No components are tracked.
<i>PicamGateTrackingMask_Width</i>	The width of the gate pulse is tracked.

5.1.1.11 *PicamGatingMode*

Data Type

PicamGatingMode is defined as enum.

Description

PicamGatingMode is the set of supported gate pulse timing modes.

Enumerator Definitions

Refer to [Table 5-11](#) for enumerator definitions.

Table 5-11: *PicamGatingMode* Enumerator Definitions

Enumerator	Description
<i>PicamGatingMode_Custom</i>	A customized modulation per readout.
<i>PicamGatingMode_Repetitive</i>	A constant gate pulse is used per readout.
<i>PicamGatingMode_Sequential</i>	A varying gate pulse is used per readout.

5.1.1.12 *PicamGatingSpeed*

Data Type

[PicamGatingSpeed](#) is defined as enum.

Description

[PicamGatingSpeed](#) is the set of classifications of the narrowest gate pulse.

Enumerator Definitions

Refer to [Table 5-12](#) for enumerator definitions.

Table 5-12: [PicamGatingSpeed](#) Enumerator Definitions

Enumerator	Description
PicamGatingSpeed_Fast	The gate pulse can be very narrow.
PicamGatingSpeed_Slow	The gate pulse width is limited by the intensifier.

5.1.1.13 *PicamIntensifierOptionsMask*

Data Type

[PicamIntensifierOptionsMask](#) is defined as enum.

Description

[PicamIntensifierOptionsMask](#) is the set of intensifier characteristics.

Enumerator Definitions

Refer to [Table 5-13](#) for enumerator definitions.

Table 5-13: [PicamIntensifierOptionsMask](#) Enumerator Definitions

Enumerator	Description
PicamIntensifierOptionsMask_Modulation	The intensifier can be modulated.
PicamIntensifierOptionsMask_SubNanosecondGating	The pulse can be gated narrower than a nanosecond.
PicamIntensifierOptionsMask_McpGating	The microchannel plate is gated instead of the photocathode.
PicamIntensifierOptionsMask_None	No additional options.

5.1.1.14 *PicamIntensifierStatus*

Data Type

`PicamIntensifierStatus` is defined as enum.

Description

`PicamIntensifierStatus` is the set of intensifier power statuses.

Enumerator Definitions

Refer to [Table 5-14](#) for enumerator definitions.

Table 5-14: `PicamIntensifierStatus` Enumerator Definitions

Enumerator	Description
<code>PicamIntensifierStatus_PoweredOff</code>	The physical switch is in the off position.
<code>PicamIntensifierStatus_PoweredOn</code>	The physical switch is in the on position.

5.1.1.15 *PicamModulationTrackingMask*

Data Type

`PicamModulationTrackingMask` is defined as enum.

Description

`PicamModulationTrackingMask` is the set of modulation parameters that are to be tracked.

Enumerator Definitions

Refer to [Table 5-15](#) for enumerator definitions.

Table 5-15: `PicamModulationTrackingMask` Enumerator Definitions

Enumerator	Description
<code>PicamModulationTrackingMask_Duration</code>	The modulation duration is tracked.
<code>PicamModulationTrackingMask_Frequency</code>	The modulation frequency is tracked.
<code>PicamModulationTrackingMask_None</code>	No components are tracked.
<code>PicamModulationTrackingMask_Output SignalFrequency</code>	The modulation output signal frequency is tracked.
<code>PicamModulationTrackingMask_Phase</code>	The modulation phase is tracked.

5.1.1.16 *PicamOrientationMask*

Data Type

`PicamOrientationMask` is defined as enum.

Description

`PicamOrientationMask` is the set of image orientation descriptors.

Enumerator Definitions

Refer to [Table 5-16](#) for enumerator definitions.

Table 5-16: `PicamOrientationMask` Enumerator Definitions

Enumerator	Description
<code>PicamOrientationMask_Flipped Horizontally</code>	The data is flipped about the centered, vertical axis relative to normal.
<code>PicamOrientationMask_Flipped Vertically</code>	The data is flipped about the centered, horizontal axis relative to normal.
<code>PicamOrientationMask_Normal</code>	This defines a standard orientation.

5.1.1.17 *PicamOutputSignal*

Data Type

`PicamOutputSignal` is defined as enum.

Description

`PicamOutputSignal` is the set of parameters defining a camera's MONITOR OUTPUT signal.

Enumerator Definitions

Refer to [Table 5-17](#) for enumerator definitions.

Table 5-17: `PicamOutputSignal` Enumerator Definitions (Sheet 1 of 2)

Enumerator	Description
<code>PicamOutputSignal_AlwaysHigh</code>	The signal is always high.
<code>PicamOutputSignal_AlwaysLow</code>	The signal is always low.
<code>PicamOutputSignal_Busy</code>	The signal is high when the camera is busy and cannot react to triggers.
<code>PicamOutputSignal_EffectivelyExposing</code>	The signal is high while the sensor is exposed to radiation.
<code>PicamOutputSignal_Exposing</code>	The signal is high for the duration of the programmed exposure time.
<code>PicamOutputSignal_NotReadingOut</code>	The signal is low when the camera is reading out.

Table 5-17: PicamOutputSignal Enumerator Definitions (Sheet 2 of 2)

Enumerator	Description
PicamOutputSignal_ReadingOut	The signal is high while reading data out of the sensor.
PicamOutputSignal_ReadyToStart	The signal is high when the camera is ready to begin acquisition.
PicamOutputSignal_ShiftingUnderMask	The signal is high while image data is shifted under the frame transfer mask.
PicamOutputSignal_ShutterOpen	The signal is high when the shutter is driven open.
PicamOutputSignal_WaitingForTrigger	The signal is high while the camera is waiting to react to a trigger.

5.1.1.18 PicamPhosphorType

Data Type

`PicamPhosphorType` is defined as enum.

Description

`PicamPhosphorType` is the set of phosphor types within an intensified camera.

Enumerator Definitions

Refer to [Table 5-18](#) for enumerator definitions.

Table 5-18: PicamPhosphorType Enumerator Definitions

Enumerator	Description
PicamPhosphorType_P43	The phosphor is P43.
PicamPhosphorType_P46	The phosphor is P46.

5.1.1.19 *PicamPhotocathodeSensitivity*

Data Type

`PicamPhotocathodeSensitivity` is defined as `enum`.

Description

`PicamPhotocathodeSensitivity` is the set of parameters used to define the photocathode's wavelength range.

Enumerator Definitions

Refer to [Table 5-19](#) for enumerator definitions.

Table 5-19: `PicamPhotocathodeSensitivity` Enumerator Definitions

Enumerator	Description
<code>PicamPhotocathodeSensitivity_HighBlueFilmless</code>	Improved quantum efficiency and optimized for blue wavelengths.
<code>PicamPhotocathodeSensitivity_HighQEFilmless</code>	Improved quantum efficiency.
<code>PicamPhotocathodeSensitivity_HighRedFilmless</code>	Improved quantum efficiency and optimized for red wavelengths.
<code>PicamPhotocathodeSensitivity_InGaAsFilmless</code>	Extends into near-infrared wavelengths.
<code>PicamPhotocathodeSensitivity_RedBlue</code>	Spans red and blue wavelengths.
<code>PicamPhotocathodeSensitivity_SolarBlind</code>	Optimized only for ultraviolet wavelengths.
<code>PicamPhotocathodeSensitivity_SuperBlue</code>	Optimized for blue wavelengths.
<code>PicamPhotocathodeSensitivity_SuperRed</code>	Optimized for red wavelengths.
<code>PicamPhotocathodeSensitivity_Unigen2Filmless</code>	Coated with UNIGEN2.
<code>PicamPhotocathodeSensitivity_UV</code>	Optimized for ultraviolet wavelengths.

5.1.1.20 *PicamPhotonDetectionMode*

Data Type

`PicamPhotonDetectionMode` is defined as enum.

Description

`PicamPhotonDetectionMode` is the set of photon detection modes.

Enumerator Definitions

Refer to [Table 5-20](#) for enumerator definitions.

Table 5-20: `PicamPhotonDetectionMode` Enumerator Definitions

Enumerator	Description
<code>PicamPhotonDetectionMode_Clipping</code>	Each pixel whose intensity is greater than or equal to the threshold is a photon and retains its original value. Otherwise the value is 0.
<code>PicamPhotonDetectionMode_Disabled</code>	Photon detection is disabled.
<code>PicamPhotonDetectionMode_Thresholding</code>	Each pixel whose intensity is greater than or equal to the threshold is a photon and replaced with a count of 1. Otherwise the value is 0.

5.1.1.21 *PicamPixelFormat*

Data Type

`PicamPixelFormat` is defined as enum.

Description

`PicamPixelFormat` is the set of characteristics that defines the format of a data pixel.

Enumerator Definitions

Refer to [Table 5-21](#) for enumerator definitions.

Table 5-21: `PicamPixelFormat` Enumerator Definitions

Enumerator	Description
<code>PicamPixelFormat_Monochrome16Bit</code>	16 bits of monochrome data

5.1.1.22 *PicamReadoutControlMode*

Data Type

[PicamReadoutControlMode](#) is defined as enum.

Description

[PicamReadoutControlMode](#) is the set of sensor readout modes.

Enumerator Definitions

Refer to [Table 5-22](#) for enumerator definitions.

Table 5-22: [PicamReadoutControlMode](#) Enumerator Definitions

Enumerator	Description
PicamReadoutControlMode_Dif	The Dual Imaging Feature where the sensor acquires two frames rapidly and then reads them both out.
PicamReadoutControlMode_FrameTransfer	The sensor is reading out a frame while exposing the next frame.
PicamReadoutControlMode_FullFrame	The sensor is read one frame at a time.
PicamReadoutControlMode_Kinetics	The sensor rapidly stores multiple frames and then reads those out.
PicamReadoutControlMode_SpectraKinetics	Same as kinetics, but optimized to capture a larger burst of spectral frames.

5.1.1.23 *PicamSensorTemperatureStatus*

Data Type

[PicamSensorTemperatureStatus](#) is defined as enum.

Description

[PicamSensorTemperatureStatus](#) is the set of sensor temperature statuses.

Enumerator Definitions

Refer to [Table 5-23](#) for enumerator definitions.

Table 5-23: [PicamSensorTemperatureStatus](#) Enumerator Definitions

Enumerator	Description
PicamSensorTemperatureStatus_Faulted	Sensor cooling has malfunctioned.
PicamSensorTemperatureStatus_Locked	The temperature has stabilized at the set point.
PicamSensorTemperatureStatus_Unlocked	The temperature has not stabilized at the set point.

5.1.1.24 *PicamSensorType*

Data Type

PicamSensorType is defined as enum.

Description

PicamSensorType is the set of sensor types.

Enumerator Definitions

Refer to [Table 5-24](#) for enumerator definitions.

Table 5-24: *PicamSensorType* Enumerator Definitions

Enumerator	Description
<i>PicamSensorType_Ccd</i>	The sensor is a CCD.
<i>PicamSensorType_InGaAs</i>	The sensor is an InGaAs.

5.1.1.25 *PicamShutterStatus*

Data Type

PicamShutterStatus is defined as enum.

Description

PicamShutterStatus is the set of shutter statuses.

Enumerator Definitions

Refer to [Table 5-25](#) for enumerator definitions.

Table 5-25: *PicamShutterStatus* Enumerator Definitions

Enumerator	Description
<i>PicamShutterStatus_Connected</i>	A shutter is connected.
<i>PicamShutterStatus_NotConnected</i>	No shutter is connected.
<i>PicamShutterStatus_Overheated</i>	A connected shutter has overheated and is temporarily disabled. NOTE: If a shutter becomes overheated, data acquisition will stop and cannot be started again until the shutter is no longer overheated.

5.1.1.26 *PicamShutterTimingMode*

Data Type

[PicamShutterTimingMode](#) is defined as enum.

Description

[PicamShutterTimingMode](#) is the set of shutter behaviors during data acquisition.

Enumerator Definitions

Refer to [Table 5-26](#) for enumerator definitions.

Table 5-26: [PicamShutterTimingMode](#) Enumerator Definitions

Enumerator	Description
PicamShutterTimingMode_Normal	The shutter only opens during exposure time. NOTE: During PicamReadoutControlMode_Kinetics readout, the shutter stays open while storing frames.
PicamShutterTimingMode_AlwaysClosed	The shutter is always closed. NOTE: This mode is also valid when not acquiring data.
PicamShutterTimingMode_AlwaysOpen	The shutter is always open. NOTE: This mode is also valid when not acquiring data.
PicamShutterTimingMode_OpenBeforeTrigger	The shutter opens ahead of time while waiting for a trigger. This is different from PicamShutterTimingMode_Normal where the shutter opens in reaction to a trigger.

5.1.1.27 *PicamShutterType*

Data Type

[PicamShutterType](#) is defined as enum.

Description

[PicamShutterType](#) is the set of shutter types.



NOTE:

This does not indicate the presence of a shutter, only the kind of shutter that could be driven. [PicamShutterStatus](#) indicates the presence of a shutter.

Enumerator Definitions

Refer to [Table 5-27](#) for enumerator definitions.

Table 5-27: [PicamShutterType](#) Enumerator Definitions (Sheet 1 of 2)

Enumerator	Description
PicamShutterType_None	No shutter.
PicamShutterType_ProntorMagnetic0	PRONTOR magnetic 0 shutter.

Table 5-27: PicamShutterType Enumerator Definitions (Sheet 2 of 2)

Enumerator	Description
PicamShutterType_ProntorMagneticE40	PRONTOR magnetic E/40 shutter.
PicamShutterType_VincentCS25	Vincent CS25 shutter
PicamShutterType_VincentCS45	Vincent CS45 shutter
PicamShutterType_VincentDSS10	Vincent DSS10 shutter
PicamShutterType_VincentVS25	Vincent VS25 shutter
PicamShutterType_VincentVS35	Vincent VS35 shutter

5.1.1.28 *PicamTimeStampsMask*

Data Type

[PicamTimeStampsMask](#) is defined as enum.

Description

[PicamTimeStampsMask](#) is the set of timestamp metadata.

Enumerator Definitions

Refer to [Table 5-28](#) for enumerator definitions.

Table 5-28: PicamTimeStampsMask Enumerator Definitions

Enumerator	Description
PicamTimeStampsMask_None	No time stamps are generated during acquisition.
PicamTimeStampsMask_ExposureStarted	The time will be stamped when exposure starts.
PicamTimeStampsMask_ExposureEnded	The time will be stamped when exposure ends.

5.1.1.29 *PicamTriggerCoupling*

Data Type

[PicamTriggerCoupling](#) is defined as enum.

Description

[PicamTriggerCoupling](#) is the set of coupling modes between an external trigger and the camera's input.

Enumerator Definitions

Refer to [Table 5-29](#) for enumerator definitions.

Table 5-29: `PicamTriggerCoupling` Enumerator Definitions

Enumerator	Description
<code>PicamTriggerCoupling_AC</code>	The components are AC-coupled.
<code>PicamTriggerCoupling_DC</code>	The components are DC-coupled.

5.1.1.30 *PicamTriggerDetermination*

Data Type

`PicamTriggerDetermination` is defined as enum.

Description

`PicamTriggerDetermination` is the set external trigger styles that are recognized by a camera.

Enumerator Definitions

Refer to [Table 5-30](#) for enumerator definitions.

Table 5-30: `PicamTriggerDetermination` Enumerator Definitions

Enumerator	Description
<code>PicamTriggerDetermination_PositivePolarity</code>	The trigger is initially a signal's rising edge and then level-sensitive to a high signal for the rest of the acquisition.
<code>PicamTriggerDetermination_NegativePolarity</code>	The trigger is initially a signal's falling edge and then level-sensitive to a low signal for the rest of the acquisition.
<code>PicamTriggerDetermination_RisingEdge</code>	The trigger is a signal's rising edge.
<code>PicamTriggerDetermination_FallingEdge</code>	The trigger is a signal's falling edge.

5.1.1.31 *PicamTriggerResponse*

Data Type

`PicamTriggerResponse` is defined as enum.

Description

`PicamTriggerResponse` is the set of a camera's responses to an external trigger.

Enumerator Definitions

Refer to [Table 5-31](#) for enumerator definitions.

Table 5-31: `PicamTriggerResponse` Enumerator Definitions

Enumerator	Description
<code>PicamTriggerResponse_NoResponse</code>	The camera acquires on its own and does not react to triggers.
<code>PicamTriggerResponse_ExposeDuringTriggerPulse</code>	Each trigger begins an exposure that lasts for the duration of the signal level (i.e., until the signal transitions back.)
<code>PicamTriggerResponse_ReadoutPerTrigger</code>	Each trigger leads to one readout. NOTE: For cameras that can non-destructively readout, all non-destructive readouts associated with the normal readout will occur on the same trigger as the normal readout.
<code>PicamTriggerResponse_ShiftPerTrigger</code>	Each trigger acquires another frame and reads out after the last frame is stored.
<code>PicamTriggerResponse_StartOnSingleTrigger</code>	The camera begins acquisition after a single trigger and continues without the need for further triggers.

5.1.1.32 *PicamTriggerSource*

Data Type

`PicamTriggerSource` is defined as enum.

Description

`PicamTriggerSource` is the set of trigger sources.

Enumerator Definitions

Refer to [Table 5-32](#) for enumerator definitions.

Table 5-32: `PicamTriggerSource` Enumerator Definitions

Enumerator	Description
<code>PicamTriggerSource_External</code>	The trigger source is external to the camera.
<code>PicamTriggerSource_Internal</code>	The trigger source is the camera.

5.1.1.33 *PicamTriggerTermination*

Data Type

`PicamTriggerTermination` is defined as enum.

Description

`PicamTriggerTermination` is the set of input terminations provided by the camera for an external trigger source.

Enumerator Definitions

Refer to [Table 5-33](#) for enumerator definitions.

Table 5-33: `PicamTriggerTermination` Enumerator Definitions

Enumerator	Description
<code>PicamTriggerTermination_FiftyOhms</code>	The trigger terminates into 50 ohms.
<code>PicamTriggerTermination_HighImpedance</code>	The trigger terminates into very high impedance.

5.1.2 Camera Parameter Access Enumerations

This section provides detailed information about the following Camera Parameter Access enumerations:

- `PicamValueAccess`

5.1.2.1 *PicamValueAccess*

Data Type

`PicamValueAccess` is defined as enum.

Description

`PicamValueAccess` is the set of permitted parameter access.

Enumerator Definitions

Refer to [Table 5-34](#) for enumerator definitions.

Table 5-34: `PicamValueAccess` Enumerator Definitions

Enumerator	Description
<code>PicamValueAccess_ReadOnly</code>	The stored parameter value can only be read.
<code>PicamValueAccess_ReadWriteTrivial</code>	The stored parameter value can be read and/or overwritten, but there is only one value for this parameter.
<code>PicamValueAccess_ReadWrite</code>	The stored parameter value can be read and/or overwritten.

5.1.3 Camera Parameter Constraint Enumerations

This section provides detailed information about the following Camera Parameter Constraint enumerations:

- [PicamConstraintScope](#)
- [PicamConstraintSeverity](#)
- [PicamConstraintCategory](#)
- [PicamRoisConstraintRulesMask](#)

5.1.3.1 *PicamConstraintScope*

Data Type

[PicamConstraintScope](#) is defined as enum.

Description

[PicamConstraintScope](#) is the set of constraint dependencies.

Enumerator Definitions

Refer to [Table 5-35](#) for enumerator definitions.

Table 5-35: [PicamConstraintScope](#) Enumerator Definitions

Enumerator	Description
PicamConstraintScope_Independent	The constraint has no dependencies and is therefore constant.
PicamConstraintScope_Dependent	The constraint has dependencies and therefore is variable.

5.1.3.2 *PicamConstraintSeverity*

Data Type

[PicamConstraintSeverity](#) is defined as enum

Description

[PicamConstraintSeverity](#) is the set of severities when failing a constraint.

Enumerator Definitions

Refer to [Table 5-36](#) for enumerator definitions.

Table 5-36: [PicamConstraintSeverity](#) Enumerator Definitions

Enumerator	Description
PicamConstraintSeverity_Error	Failure indicates the value is in error.
PicamConstraintSeverity_Warning	Failure indicates the value is in warning and notice should be taken.

5.1.3.3 *PicamConstraintCategory*

Data Type

`PicamConstraintCategory` is defined as enum.

Description

`PicamConstraintCategory` is the set of constraint categories.

Enumerator Definitions

Refer to [Table 5-37](#) for enumerator definitions.

Table 5-37: `PicamConstraintCategory` Enumerator Definitions

Enumerator	Description
<code>PicamConstraintCategory_Capable</code>	Which set of values are ultimately possible.
<code>PicamConstraintCategory_Required</code>	Which set of values are currently permissible.
<code>PicamConstraintCategory_Recommended</code>	Which set of values fall within a recommended range for most scenarios.

5.1.3.4 *PicamRoisConstraintRulesMask*

Data Type

`PicamRoisConstraintRulesMask` is defined as enum.

Description

`PicamRoisConstraintRulesMask` is the set of complex rules that defines a valid set of regions of interest.

Enumerator Definitions

Refer to [Table 5-38](#) for enumerator definitions.

Table 5-38: `PicamRoisConstraintRulesMask` Enumerator Definitions

Enumerator	Description
<code>PicamRoisConstraintRulesMask_None</code>	No additional rules.
<code>PicamRoisConstraintRulesMask_XBinningAlignment</code>	Regions sharing columns must bin those columns equally. This means not only must they contain equal x-binning values, the regions must also begin on x-binning boundaries.
<code>PicamRoisConstraintRulesMask_YBinningAlignment</code>	Regions sharing rows must bin those rows equally. This means not only must they contain equal y-binning values, the regions must also begin on y-binning boundaries.
<code>PicamRoisConstraintRulesMask_HorizontalSymmetry</code>	Regions must be symmetrical about the line between the two center-most columns. Either one region must bisect this line or two regions must be reflective to each other about this line.
<code>PicamRoisConstraintRulesMask_VerticalSymmetry</code>	Regions must be symmetrical about the line between the two center-most rows. Either one region must bisect this line or two regions must be reflective to each other about this line.
<code>PicamRoisConstraintRulesMask_SymmetryBoundsBinning</code>	A region required to bisect a line of symmetry may not bin pixels together that fall on both sides of the line.

5.2 Structure Definitions

This section provides programming information about PICam data structure definitions.

5.2.1 Camera Parameter Structures

This section provides detailed programming information about the following Camera Parameter data structures:

- [PicamRoi](#)
- [PicamRoIs](#)
- [PicamPulse](#)
- [PicamModulation](#)
- [PicamModulations](#)
- [PicamStatusPurview](#)

5.2.1.1 *PicamRoi*

Description

[PicamRoi](#) defines a single Region of Interest (ROI.)

Structure Definition

The structure definition for [PicamRoi](#) is:

```
typedef struct PicamRoi
{
    piint  x;
    piint  width;
    piint  x_binning;
    piint  y;
    piint  height;
    piint  y_binning;
} PicamRoi;
```

Variable Definitions

The variables required by [PicamRoi](#) are:

- `x`: The left-most column coordinate (zero-based).
- `width`: The number of columns.
- `x_binning`: The number of columns to group into a sum.
- `y`: The top-most row coordinate (zero-based).
- `height`: The number of rows.
- `y_binning`: The number of rows to group into a sum.

5.2.1.2 *PicamRoIs*

Description

PicamRoIs defines a set of non-overlapping Regions of Interest (ROIs.)

Structure Definition

The structure definition for *PicamRoIs* is:

```
typedef struct PicamRoIs
{
    PicamRoi*  roi_array;
    piint  roi_count;
} PicamRoIs;
```

Variable Definitions

The variables required by *PicamRoIs* are:

roi_array: An array of one or more regions.

roi_count: The number of regions.

5.2.1.3 *PicamPulse*

Description

PicamPulse defines a gate pulse.

Structure Definition

The structure definition for *PicamPulse* is:

```
typedef struct PicamPulse
{
    piflt  delay;
    piflt  width;
} PicamPulse;
```

Variable Definitions

The variables required by *PicamPulse* are:

delay: The delay until a gate pulse begins.

width: The width of the gate pulse.

5.2.1.4 *PicamModulation*

Description

[PicamModulation](#) defines a custom intensifier modulation sequence point.

Structure Definition

The structure definition for [PicamModulation](#) is:

```
typedef struct PicamModulation
{
    piflt duration;
    piflt frequency;
    piflt phase;
    piflt output_signal_frequency;
} PicamModulation;
```

Variable Definitions

The variables required by [PicamModulation](#) are:

duration: The time, in mS, the intensifier is modulating.

frequency: The frequency, in MHz, of the intensifier modulation.

phase: The phase, in degrees, of the intensifier with respect to the modulation output signal.

output_signal_frequency: The frequency, in MHz, of the user RF output signal

5.2.1.5 *PicamModulations*

Description

[PicamModulations](#) defines a sequence of intensifier modulation sequence points.

Structure Definition

The structure definition for [PicamModulations](#) is:

```
typedef struct PicamModulations
{
    PicamModulation* modulation_array;
    pint modulation_count;
} PicamModulations;
```

Variable Definitions

The variables required by [PicamModulations](#) are:

modulation_array: An array of one or more sequence points.

modulation_count: The number of sequence points.

5.2.1.6 *PicamStatusPurview*

Description

PicamStatusPurview defines the scope of a status.

Structure Definition

The structure definition for *PicamStatusPurview* is:

```
typedef struct PicamStatusPurview
{
    const puint*  values_array;
    puint  values_count;
} PicamStatusPurview;
```

Variable Definitions

The variables required by *PicamStatusPurview* are:

values_array: The allowable status values.

values_count: The number of allowable status values.

5.2.2 Camera Parameter Constraints

This section provides detailed programming information about the following Camera Parameter Constraint data structures:

- [PicamCollectionConstraint](#)
- [PicamRangeConstraint](#)
- [PicamRoisConstraint](#)
- [PicamPulseConstraint](#)
- [PicamModulationsConstraint](#)

5.2.2.1 *PicamCollectionConstraint*

Description

[PicamCollectionConstraint](#) defines the constraints placed on a variable whose value is selected from a list of predefined values.

Structure Definition

The structure definition for [PicamCollectionConstraint](#) is:

```
typedef struct PicamCollectionConstraint
{
    PicamConstraintScope  scope;
    PicamConstraintSeverity severity;
    const piflt* values_array;
    piint values_count;
} PicamCollectionConstraint;
```

Variable Definitions

The variables required by [PicamCollectionConstraint](#) are:

scope: The scope of the constraint.

severity: The severity of the constraint.

values_array: The allowable values.

values_count: The number of allowable values.

5.2.2.2 *PicamRangeConstraint*

Description

PicamRangeConstraint defines the constraints placed a numeric variable whose value lies within a linear range of numeric values.

Structure Definition

The structure definition for *PicamRangeConstraint* is:

```
typedef struct PicamRangeConstraint
{
    PicamConstraintScope    scope;
    PicamConstraintSeverity severity;
    pibln    empty_set;
    piflt    minimum;
    piflt    maximum;
    piflt    increment;
    const piflt* excluded_values_array;
    piint    excluded_values_count;
    const piflt* outlying_values_array;
    piint    outlying_values_count;
} PicamRangeConstraint;
```

Variable Definitions

The variables required by *PicamRangeConstraint* are:

- `scope`: The scope of the constraint.
- `severity`: The severity of the constraint.
- `empty_set`: Indicates when there are no valid values within the range.
Valid values are:
 - `TRUE`
There are no valid values within the range.
When `TRUE`, only `scope` and `severity` are relevant.
 - `FALSE`
There is at least one valid value within the range.
- `minimum`: The smallest value within the range.
NOTE: `outlying_values_array` may include a smaller value.
- `maximum`: The largest value within the range.
NOTE: `outlying_values_array` may include a larger value.
- `increment`: The numeric gap between consecutive values within the range.
- `excluded_values_array`: The set of values within the range (excluding `minimum` and `maximum`) that is not valid.
This is null when all values within the range are valid.

continued on next page

continued from previous page

- `excluded_values_count`: The number of items within `excluded_values_array`.
This is 0 when there are no excluded values.
- `outlying_values_array`: The set of valid values that lie outside of the range of values.
This is null when no valid values fall outside of the range.
- `outlying_values_count`: The number of items within `outlying_values_array`.
This is 0 when there are no outlying values.

5.2.2.3 *PicamRoisConstraint*

Description

`PicamRoisConstraint` defines the constraints placed on a set of Regions of Interest (ROIs).



NOTE:

Regions of Interest may not overlap.

Structure Definition

The structure definition for `PicamRoisConstraint` is:

```
typedef struct PicamRoisConstraint
{
    PicamConstraintScope    scope;
    PicamConstraintSeverity severity;
    pibln                   empty_set;
    PicamRoisConstraintRulesMask rules;
    piint                    maximum_roi_count;
    PicamRangeConstraint    x_constraint;
    PicamRangeConstraint    width_constraint;
    const piint*             x_binning_limits_array;
    piint                    x_binning_limits_count;
    PicamRangeConstraint    y_constraint;
    PicamRangeConstraint    height_constraint;
    const piint*             y_binning_limits_array;
    piint                    y_binning_limits_count;
} PicamRoisConstraint;
```

Variable Definitions

The variables required by `PicamRoisConstraint` are:

<code>scope</code> :	The scope of the constraint.
<code>severity</code> :	The severity of the constraint.
<code>empty_set</code> :	Indicates when there are no valid Regions of Interest defined. Valid values are: <ul style="list-style-type: none"> • <code>TRUE</code> There are no valid ROIs defined. When <code>TRUE</code>, only <code>scope</code> and <code>severity</code> are relevant. • <code>FALSE</code> There is at least one valid ROI defined.
<code>rules</code> :	Complex set of rules to which a parameter of this type must adhere.
<code>maximum_roi_count</code> :	The maximum number of ROIs permitted.
<code>x_constraint</code> :	The constraint governing the value of <code>PicamRoi.x</code> .
<code>width_constraint</code> :	The constraint governing the value of <code>PicamRoi.width</code> .
<code>x_binning_limits_array</code> :	The list of valid values for <code>PicamRoi.x_binning</code> . NOTE: An additional requirement is that <code>PicamRoi.x_binning</code> must always divide evenly into <code>PicamRoi.width</code> . This is null when no additional limits are required.
<code>x_binning_limits_count</code> :	The number of items in <code>x_binning_limits_array</code> . This is 0 when no additional limits are required.
<code>y_constraint</code> :	The constraint governing the value of <code>PicamRoi.y</code> .
<code>height_constraint</code> :	The constraint governing the value of <code>PicamRoi.height</code> .
<code>y_binning_limits_array</code> :	The list of valid values for <code>PicamRoi.y_binning</code> . NOTE: An additional requirement is that <code>PicamRoi.y_binning</code> must always divide evenly into <code>PicamRoi.height</code> . This is null when no additional limits are required.
<code>y_binning_limits_count</code> :	The number of items in <code>y_binning_limits_array</code> . This is 0 when no additional limits are required.

5.2.2.4 *PicamPulseConstraint*

Description

`PicamPulseConstraint` defines the constraints placed on a valid gate pulse.

Structure Definition

The structure definition for `PicamPulseConstraint` is:

```
typedef struct PicamPulseConstraint
{
    PicamConstraintScope    scope;
    PicamConstraintSeverity severity;
    pibln                   empty_set;
    PicamRangeConstraint    delay_constraint;
    PicamRangeConstraint    width_constraint;
    piflt                   minimum_duration;
    piflt                   maximum_duration;
} PicamPulseConstraint;
```

Variable Definitions

The variables required by `PicamPulseConstraint` are:

`scope`: The scope of the constraint.

`severity`: The severity of the constraint.

`empty_set`: Indicates when there are no valid Pulses defined.

Valid values are:

- TRUE
There are no valid Pulses defined.
When TRUE, only `scope` and `severity` are relevant.
- FALSE
There is at least one valid Pulse defined.

`delay_constraint`: The constraint governing the value of `PicamPulse.delay`.

`width_constraint`: The constraint governing the value of `PicamPulse.width`.

`minimum_duration`: The minimum numeric value for:
[`PicamPulse.delay` + `PicamPulse.width`]

`maximum_duration`: The maximum numeric value for:
[`PicamPulse.delay` + `PicamPulse.width`]

5.2.2.5 *PicamModulationsConstraint*

Description

[PicamModulationsConstraint](#) defines the constraints placed on custom intensifier modulation sequence points.

Structure Definition

The structure definition for [PicamModulationsConstraint](#) is:

```
typedef struct PicamModulationsConstraint
{
    PicamConstraintScope    scope;
    PicamConstraintSeverity severity;
    pibln                   empty_set;
    piint                   maximum_modulation_count;
    PicamRangeConstraint    duration_constraint;
    PicamRangeConstraint    frequency_constraint;
    PicamRangeConstraint    phase_constraint;
    PicamRangeConstraint    output_signal_frequency_constraint;
} PicamModulationsConstraint;
```

Variable Definitions

The variables required by [PicamModulationsConstraint](#) are:

scope:	The scope of the constraint.
severity:	The severity of the constraint.
empty_set:	Indicates when there are no valid modulation points defined. Valid values are: <ul style="list-style-type: none"> • TRUE There are no valid modulation points defined. When TRUE, only scope and severity are relevant. • FALSE There is at least one valid modulation point defined.
maximum_modulation_count:	The maximum number of modulation sequence points.
duration_constraint:	The constraint governing the value of <code>PicamModulation.duration</code> .
frequency_constraint:	The constraint governing the value of <code>PicamModulation.frequency</code> .
phase_constraint:	The constraint governing the value of <code>PicamModulation.phase</code> .
output_signal_frequency_constraint:	The constraint governing the value of <code>PicamModulation.output_signal_frequency</code> .

5.3 Programmers' Reference for Camera Configuration APIs

This section provides a detailed programmers' reference guide for the following APIs:

- Camera Parameter Value APIs
 - `Picam_GetParameterIntegerValue()`
 - `Picam_CanSetParameterIntegerValue()`
 - `Picam_SetParameterIntegerValue()`
 - `Picam_GetParameterLargeIntegerValue()`
 - `Picam_CanSetParameterLargeIntegerValue()`
 - `Picam_SetParameterLargeIntegerValue()`
 - `Picam_GetParameterFloatingPointValue()`
 - `Picam_CanSetParameterFloatingPointValue()`
 - `Picam_SetParameterFloatingPointValue()`
 - `Picam_DestroyRois()`
 - `Picam_GetParameterRoisValue()`
 - `Picam_CanSetParameterRoisValue()`
 - `Picam_SetParameterRoisValue()`
 - `Picam_DestroyPulses()`
 - `Picam_GetParameterPulseValue()`
 - `Picam_CanSetParameterPulseValue()`
 - `Picam_SetParameterPulseValue()`
 - `Picam_DestroyModulations()`
 - `Picam_GetParameterModulationsValue()`
 - `Picam_CanSetParameterModulationsValue()`
 - `Picam_SetParameterModulationsValue()`
 - `Picam_GetParameterIntegerDefaultValue()`
 - `Picam_GetParameterLargeIntegerDefaultValue()`
 - `Picam_GetParameterFloatingPointDefaultValue()`
 - `Picam_GetParameterRoisPointDefaultValue()`
 - `Picam_GetParameterPulseDefaultValue()`
 - `Picam_GetParameterModulationsDefaultValue()`
 - `Picam_RestoreParametersToDefaultValues()`
 - `Picam_CanSetParameterOnline()`
 - `Picam_SetParameterIntegerValueOnline()`
 - `Picam_SetParameterFloatingPointValueOnline()`
 - `Picam_SetParameterPulseValueOnline()`
 - `Picam_CanReadParameter()`
 - `Picam_ReadParameterIntegerValue()`
 - `Picam_ReadParameterFloatingPointValue()`
 - `Picam_CanWaitForStatusParameter()`
 - `Picam_DestroyStatusPurviews()`
 - `Picam_GetStatusParameterPurview()`
 - `Picam_EstimateTimeToStatusParameterValue()`
 - `Picam_WaitForStatusParameterValue()`
- Camera Parameter Information APIs
 - `Picam_DestroyParameters()`
 - `Picam_GetParameters()`
 - `Picam_DoesParameterExist()`
 - `Picam_IsParameterRelevant()`
 - `Picam_GetParameterValueType()`
 - `Picam_GetParameterEnumeratedType()`
 - `Picam_GetParameterValueAccess()`
 - `Picam_GetParameterConstraintType()`

-
- **Camera Parameter Constraints APIs**
 - `Picam_DestroyCollectionConstraints()`
 - `Picam_GetParameterCollectionConstraint()`
 - `Picam_DestroyRangeConstraints()`
 - `Picam_GetParameterRangeConstraint()`
 - `Picam_DestroyRoisConstraints()`
 - `Picam_GetParameterRoisConstraint()`
 - `Picam_DestroyPulseConstraints()`
 - `Picam_GetParameterPulseConstraint()`
 - `Picam_DestroyModulationsConstraints()`
 - `Picam_GetParameterModulationsConstraint()`
 - **Camera Parameter Commitment APIs**
 - `Picam_AreParametersCommitted()`
 - `Picam_CommitParameters()`

5.3.1 Camera Parameter Value APIs

This section provides programming information for APIs used when working with parameter values.

5.3.1.1 *Picam_GetParameterIntegerValue()*

Description

`Picam_GetParameterIntegerValue()` returns the integer value for a specified parameter.

Syntax

The syntax for `Picam_GetParameterIntegerValue()` is:

```
PICAM_API Picam_GetParameterIntegerValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    piint* value);
```

Input Parameters

Input parameters for `Picam_GetParameterIntegerValue()` are:

`camera`: Handle for the camera for which the integer value is being requested.

`parameter`: Specifies the parameter that is to be queried.

Valid parameters are those of type:

- `PicamValueType_Integer`;
- `PicamValueType_Boolean`;
- `PicamValueType_Enumeration`.

Output Parameters

Output parameters for `Picam_GetParameterIntegerValue()` are:

`value`: Pointer to the integer value of the specified parameter.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `model`, or
- `device`.

Stored values for any specific parameter are not necessarily the same for the `device` and `model` instances.

5.3.1.2 *Picam_CanSetParameterIntegerValue()*

Description

`Picam_CanSetParameterIntegerValue()` determines if an integer value is valid for a specified parameter.

Syntax

The syntax for `Picam_CanSetParameterIntegerValue()` is:

```
PICAM_API Picam_CanSetParameterIntegerValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    piint value,  
    pibln* settable);
```

Input Parameters

Input parameters for `Picam_CanSetParameterIntegerValue()` are:

- camera: Handle for the camera for which the value/parameter combination is being validated.
- parameter: Specifies the parameter which is to be tested.
- value: The integer value that is to be tested.

Output Parameters

Output parameters for `Picam_CanSetParameterIntegerValue()` are:

- settable: Pointer to the test results. Indicates if the integer value is a valid value for the specified parameter.
Valid values are:
 - TRUE
Indicates that the integer value is a valid value for the specified parameter.
 - FALSE
Indicates that the integer value is an invalid value for the specified parameter.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `model`, or
- `device`.

5.3.1.3 *Picam_SetParameterIntegerValue()*

Description

`Picam_SetParameterIntegerValue()` sets a parameter to a specified integer value during camera setup.

Syntax

The syntax for `Picam_SetParameterIntegerValue()` is:

```
PICAM_API Picam_SetParameterIntegerValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    piint value);
```

Input Parameters

Input parameters for `Picam_SetParameterIntegerValue()` are:

`camera`: Handle for the camera being configured.

`parameter`: Specifies the parameter that is to be set with an integer value.

Valid parameters are those of type:

- `PicamValueType_Integer`;
- `PicamValueType_Boolean`;
- `PicamValueType_Enumeration`.

`value`: The integer value to which the parameter is to be set.

Output Parameters

There are no output parameters associated with `Picam_SetParameterIntegerValue()`.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `model`;
The `model` parameter may be set independently from the corresponding `device` parameter. However, doing so requires that all parameters be committed to the device prior to starting any data acquisition by calling `Picam_CommitParameters()`.
- `device`.
Setting a `device` parameter automatically sets the corresponding `model` parameter to the same value.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_CommitParameters()`.

5.3.1.4 *Picam_GetParameterLargeIntegerValue()*

Description

`Picam_GetParameterLargeIntegerValue()` returns the current large integer value for a specified parameter.

Syntax

The syntax for `Picam_GetParameterLargeIntegerValue()` is:

```
PICAM_API Picam_GetParameterLargeIntegerValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    pi64s* value);
```

Input Parameters

Input parameters for `Picam_GetParameterLargeIntegerValue()` are:

- `camera`: Handle for the camera for which the large integer value is being requested.
- `parameter`: Specifies the parameter that is to be queried.
Valid parameters are those of type `PicamValueType_LargeInteger`.

Output Parameters

Output parameters for `Picam_GetParameterLargeIntegerValue()` are:

- `value`: Pointer to the large integer value of the specified parameter.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `model`, or
- `device`.

Stored values for any specific parameter are not necessarily the same for the `device` and `model` instances.

5.3.1.5 *Picam_CanSetParameterLargeIntegerValue()*

Description

`Picam_CanSetParameterLargeIntegerValue()` determines if a large integer value is valid for a specified parameter.

Syntax

The syntax for `Picam_CanSetParameterLargeIntegerValue()` is:

```
PICAM_API Picam_CanSetParameterLargeIntegerValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    pi64s value,  
    pibln* settable);
```

Input Parameters

Input parameters for `Picam_CanSetParameterLargeIntegerValue()` are:

- `camera`: Handle for the camera for which the value/parameter combination is being tested.
- `parameter`: Specifies the parameter which is to be tested.
- `value`: The large integer value that is to be tested.

Output Parameters

Output parameters for `Picam_CanSetParameterLargeIntegerValue()` are:

- `settable`: Pointer to the test results. Indicates if the large integer value is a valid value for the specified parameter.
Valid values are:
 - `TRUE`
Indicates that the large integer value is a valid value for the specified parameter.
 - `FALSE`
Indicates that the large integer value is an invalid value for the specified parameter.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `model`, or
- `device`.

5.3.1.6 *Picam_SetParameterLargeIntegerValue()*

Description

`Picam_SetParameterLargeIntegerValue()` sets a parameter to a specified large integer value during camera setup.

Syntax

The syntax for `Picam_SetParameterLargeIntegerValue()` is:

```
PICAM_API Picam_SetParameterLargeIntegerValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    pi64s value);
```

Input Parameters

Input parameters for `Picam_SetParameterLargeIntegerValue()` are:

- camera: Handle for the camera being configured.
- parameter: Specifies the parameter that is to be set with a large integer value.
Valid parameters are those of type
`PicamValueType_LargeInteger`.
- value: The large integer value to which the parameter is to be set.

Output Parameters

There are no output parameters associated with
`Picam_SetParameterLargeIntegerValue()`.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `model`;
The `model` parameter may be set independently from the corresponding `device` parameter. However, doing so requires that all parameters be committed to the device prior to starting any data acquisition by calling
`Picam_CommitParameters()`.
- `device`.
Setting a `device` parameter automatically sets the corresponding `model` parameter to the same value.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_CommitParameters()`

5.3.1.7 *Picam_GetParameterFloatingPointValue()*

Description

`Picam_GetParameterFloatingPointValue()` returns the current floating point value for a specified parameter.

Syntax

The syntax for `Picam_GetParameterFloatingPointValue()` is:

```
PICAM_API Picam_GetParameterFloatingValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    piflt* value);
```

Input Parameters

Input parameters for `Picam_GetParameterFloatingPointValue()` are:

- `camera`: Handle for the camera for which the floating point value is being requested.
- `parameter`: Specifies the parameter that is to be queried.
Valid parameters are those of type `PicamValueType_FloatingPoint`.

Output Parameters

Output parameters for `Picam_GetParameterFloatingPointValue()` are:

- `value`: Pointer to the floating point value of the specified parameter.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `model`, or
- `device`.

Stored values for any specific parameter are not necessarily the same for the `device` and `model` instances.

5.3.1.8 *Picam_CanSetParameterFloatingPointValue()*

Description

Picam_CanSetParameterFloatingPointValue() determines if a floating point value is valid for a specified parameter.

Syntax

The syntax for *Picam_CanSetParameterFloatingPointValue()* is:

```
PICAM_API Picam_CanSetParameterFloatingValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    piflt value,  
    pibln* settable);
```

Input Parameters

Input parameters for *Picam_CanSetParameterFloatingPointValue()* are:

- camera: Handle for the camera for which the value/parameter combination is being validated.
- parameter: Specifies the parameter which is to be tested.
- value: The floating point value that is to be tested.

Output Parameters

Output parameters for *Picam_CanSetParameterFloatingPointValue()* are:

- settable: Pointer to the test results. Indicates if the floating point value is a valid value for the specified parameter.

Valid values are:

- TRUE
Indicates that the floating point value is a valid value for the specified parameter.
- FALSE
Indicates that the floating point value is an invalid value for the specified parameter.

Advanced API Usage

When used in conjunction with Advanced APIs, *camera* may be a handle to either the:

- *model*, or
- *device*.

5.3.1.9 *Picam_SetParameterFloatingPointValue()*

Description

`Picam_SetParameterFloatingPointValue()` sets a parameter to a specified floating point value during camera setup.

Syntax

The syntax for `Picam_SetParameterFloatingPointValue()` is:

```
PICAM_API Picam_SetParameterFloatingValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    piflt value);
```

Input Parameters

Input parameters for `Picam_SetParameterFloatingPointValue()` are:

camera: Handle for the camera being configured.

parameter: Specifies the parameter that is to be set with a floating point value.
Valid parameters are those of type
`PicamValueType_FloatingPoint`.

value: The floating point value to which the parameter is to be set.

Output Parameters

There are no output parameters associated with
`Picam_SetParameterFloatingPointValue()`.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `model`;
The `model` parameter may be set independently from the corresponding `device` parameter. However, doing so requires that all parameters be committed to the device prior to starting any data acquisition by calling
`Picam_CommitParameters()`.
- `device`.
Setting a `device` parameter automatically sets the corresponding `model` parameter to the same value.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_CommitParameters()`

5.3.1.10 *Picam_DestroyRois()*

Description

`Picam_DestroyRois()` releases memory that has been allocated by PICam for use by the array `rois`.

If `rois` is null, calling `Picam_DestroyRois()` has no effect.

Syntax

The syntax for `Picam_DestroyRois()` is:

```
PICAM_API Picam_DestroyRois(  
    const PicamRois* rois);
```

Input Parameters

Input parameters for `Picam_DestroyRois()` are:

`rois`: Pointer to array memory that is to be released.

Output Parameters

There are no output parameters associated with `Picam_DestroyRois()`.

5.3.1.11 *Picam_GetParameterRoisValue()*

Description

`Picam_GetParameterRoisValue()` returns the current value for a specified Rois parameter.

Syntax

The syntax for `Picam_GetParameterRoisValue()` is:

```
PICAM_API Picam_GetParameterRoisValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamRois** value);
```

Input Parameters

Input parameters for `Picam_GetParameterRoisValue()` are:

- `camera`: Handle for the camera for which the Rois parameter value is being requested.
- `parameter`: Specifies the Rois parameter for which the current value is to be returned.
Valid parameters are those of type `PicamValueType_Rois`.

Output Parameters

Output parameters for `Picam_GetParameterRoisValue()` are:

- `value`: Pointer to the memory location in which the value of the specified Rois parameter has been stored.
NOTE: This memory is allocated by PICam and must be released by calling `Picam_DestroyRois()`.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `model`, or
- `device`.

Stored values for any specific parameter are not necessarily the same for the `device` and `model` instances.

Related Structures

For additional information, refer to the following ROI structures:

- `PicamRoi`;
- `PicamRois`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_DestroyRois()`

5.3.1.12 *Picam_CanSetParameterRoisValue()*

Description

`Picam_CanSetParameterRoisValue()` determines if a value is valid for a specified Rois parameter.

Syntax

The syntax for `Picam_CanSetParameterRoisValue()` is:

```
PICAM_API Picam_CanSetParameterRoisValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamRois* value,  
    pibln* settable);
```

Input Parameters

Input parameters for `Picam_CanSetParameterRoisValue()` are:

- camera: Handle for the camera for which the value/parameter combination is being validated.
- parameter: Specifies the Rois parameter.
- value: The value that is to be tested.

Output Parameters

Output parameters for `Picam_CanSetParameterRoisValue()` are:

- settable: Pointer to the test results. Indicates if the value is valid for the specified Rois parameter.
Valid values are:
 - TRUE
Indicates that the value is valid for the specified Rois parameter.
 - FALSE
Indicates that the value is not valid for the specified Rois parameter.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `model`, or
- `device`.

Related Structures

For additional information, refer to the following ROI constraint structures:

- `PicamRoisConstraint`.

5.3.1.13 *Picam_SetParameterRoisValue()*

Description

`Picam_SetParameterRoisValue()` configures an Rois parameter to a specified value during camera setup.

Syntax

The syntax for `Picam_SetParameterRoisValue()` is:

```
PICAM_API Picam_SetParameterRoisValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamRois* value);
```

Input Parameters

Input parameters for `Picam_SetParameterRoisValue()` are:

- camera: Handle for the camera being configured.
- parameter: Specifies the Rois parameter that is to be configured.
Valid parameters are those of type `PicamValueType_Rois`.
- value: The value to which the Rois parameter is to be set.

Output Parameters

There are no output parameters associated with `Picam_SetParameterRoisValue()`.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `model`;
The `model` parameter may be set independently from the corresponding `device` parameter. However, doing so requires that all parameters be committed to the device prior to starting any data acquisition by calling `Picam_CommitParameters()`.
- `device`.
Setting a `device` parameter automatically sets the corresponding `model` parameter to the same value.

Related Structures

For additional information, refer to the following ROI structures:

- `PicamRoi`;
- `PicamRois`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_CommitParameters()`

5.3.1.14 *Picam_DestroyPulses()*

Description

`Picam_DestroyPulses()` releases memory that has been allocated by PICam for use by pulses.

If `pulses` is null, calling `Picam_DestroyPulses()` has no effect.

Syntax

The syntax for `Picam_DestroyPulses()` is:

```
PICAM_API Picam_DestroyPulses(  
    const PicamPulse* pulses);
```

Input Parameters

Input parameters for `Picam_DestroyPulses()` are:

`pulses`: Pointer to array memory that is to be released.

Output Parameters

There are no output parameters associated with `Picam_DestroyPulses()`.

5.3.1.15 *Picam_GetParameterPulseValue()*

Description

Picam_GetParameterPulseValue() returns the current value for a specified Pulse parameter.

Syntax

The syntax for *Picam_GetParameterPulseValue()* is:

```
PICAM_API Picam_GetParameterPulseValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamPulse** value);
```

Input Parameters

Input parameters for *Picam_GetParameterPulseValue()* are:

- camera: Handle for the camera for which the specified pulse parameter value is being requested.
- parameter: Specifies the Pulse parameter for which the current value is to be returned.
Valid parameters are those of type *PicamValueType_Pulse*.

Output Parameters

Output parameters for *Picam_GetParameterPulseValue()* are:

- value: Pointer to the memory location where the value of the specified Pulse parameter has been stored.
NOTE: This memory is allocated by PICam and must be released by calling *Picam_DestroyPulses()*

Advanced API Usage

When used in conjunction with Advanced APIs, *camera* may be a handle to either the:

- *device*, or
- *model*.

Stored values for any specific parameter are not necessarily the same for the *device* and *model* instances.

Related Structures

For additional information, refer to the following Pulse structure:

- *PicamPulse*.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_DestroyPulses()*

5.3.1.16 *Picam_CanSetParameterPulseValue()*

Description

`Picam_CanSetParameterPulseValue()` determines if a value is valid for a specified Pulse parameter.

Syntax

The syntax for `Picam_CanSetParameterPulseValue()` is:

```
PICAM_API Picam_CanSetParameterPulseValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamPulse* value,  
    pibln* settable);
```

Input Parameters

Input parameters for `Picam_CanSetParameterPulseValue()` are:

- camera: Handle for the camera for which the value/parameter combination is being validated.
- parameter: Specifies the Pulse parameter.
- value: The value that is to be tested.

Output Parameters

Output parameters for `Picam_CanSetParameterPulseValue()` are:

- settable: Pointer to the test results. Indicates if the value is valid for the specified Pulse parameter.
Valid values are:
 - TRUE
Indicates that the value is valid for the specified Pulse parameter.
 - FALSE
Indicates that the value is not valid for the specified Pulse parameter.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

Related Structures

For additional information, refer to the following Pulse constraint structure:

- `PicamPulseConstraint`.

5.3.1.17 *Picam_SetParameterPulseValue()*

Description

`Picam_SetParameterPulseValue()` configures a Pulse parameter to a specified value during camera setup.

Syntax

The syntax for `Picam_SetParameterPulseValue()` is:

```
PICAM_API Picam_SetParameterPulseValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamPulse* value);
```

Input Parameters

Input parameters for `Picam_SetParameterPulseValue()` are:

- camera: Handle for the camera being configured.
- parameter: Specifies the Pulse parameter that is to be configured.
Valid parameters are those of type `PicamValueType_Pulse`.
- value: The value to which the Pulse parameter is to be set.

Output Parameters

There are no output parameters associated with `Picam_SetParameterPulseValue()`

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `model`;
The `model` parameter may be set independently from the corresponding `device` parameter. However, doing so requires that all parameters be committed to the device prior to starting any data acquisition by calling `Picam_CommitParameters()`.
- `device`.
Setting a `device` parameter automatically sets the corresponding `model` parameter to the same value.

Related Structures

For additional information, refer to the following Pulse structure:

- `PicamPulse`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_CommitParameters()`

5.3.1.18 *Picam_DestroyModulations()*

Description

`Picam_DestroyModulations()` releases memory that has been allocated by PICam for use by modulations.

If `modulations` is null, calling `Picam_DestroyModulations()` has no effect.

Syntax

The syntax for `Picam_DestroyModulations()` is:

```
PICAM_API Picam_DestroyModulations(  
    const PicamModulations* modulations);
```

Input Parameters

Input parameters for `Picam_DestroyModulations()` are:

`modulations`: Pointer to array memory that is to be released.

Output Parameters

There are no output parameters associated with `Picam_DestroyModulations()`.

5.3.1.19 *Picam_GetParameterModulationsValue()*

Description

Picam_GetParameterModulationsValue() returns the current value for a specified intensifier modulation sequence parameter.

Syntax

The syntax for *Picam_GetParameterModulationsValue()* is:

```
PICAM_API Picam_GetParameterModulationsValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamModulations** value);
```

Input Parameters

Input parameters for *Picam_GetParameterModulationsValue()* are:

- `camera`: Handle for the camera for which the intensifier modulation sequence parameter value is being requested..
- `parameter`: Specifies the intensifier modulation sequence parameter for which the current value is to be returned.
Valid parameters are those of type *PicamValueType_Modulations*.

Output Parameters

Output parameters for *Picam_GetParameterModulationsValue()* are:

- `value`: Pointer to the memory location in which the value of the specified intensifier modulation sequence parameter is stored.

NOTE: This memory is allocated by PICam and must be released by calling *Picam_DestroyModulations()*

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

Stored values for any specific parameter are not necessarily the same for the `device` and `model` instances.

Related Structures

For additional information, refer to the following intensifier modulation sequence structures:

- *PicamModulation*;
- *PicamModulations*.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_DestroyModulations()*.

5.3.1.20 *Picam_CanSetParameterModulationsValue()*

Description

`Picam_CanSetParameterModulationsValue()` determines if a value is valid for a specified intensifier modulation sequence parameter.

Syntax

The syntax for `Picam_CanSetParameterModulationsValue()` is:

```
PICAM_API Picam_CanSetParameterModulationsValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamModulations* value,  
    pibln* settable);
```

Input Parameters

Input parameters for `Picam_CanSetParameterModulationsValue()` are:

- camera: Handle for the camera for which the value/parameter combination is being validated.
- parameter: Specifies the intensifier modulation sequence parameter.
- value: The value that is to be tested.

Output Parameters

Output parameters for `Picam_CanSetParameterModulationsValue()` are:

- settable: Pointer to the test results. Indicates if the value is valid for the specified intensifier modulation sequence parameter.
Valid values are:
 - TRUE
Indicates that the value is valid for the specified intensifier modulation sequence parameter.
 - FALSE
Indicates that the value is not valid for the specified intensifier modulation sequence parameter.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

Related Structures

For additional information, refer to the following intensifier modulation sequence structures:

- `PicamModulation`;
- `PicamModulations`.

5.3.1.21 *Picam_SetParameterModulationsValue()*

Description

`Picam_SetParameterModulationsValue()` configures an intensifier modulation sequence parameter to a specified value during camera setup.

Syntax

The syntax for `Picam_SetParameterModulationsValue()` is:

```
PICAM_API Picam_SetParameterModulationsValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamModulations* value);
```

Input Parameters

Input parameters for `Picam_SetParameterModulationsValue()` are:

- camera: Handle for the camera being configured.
- parameter: Specifies the Pulse parameter that is to be configured.
Valid parameters are those of type `PicamValueType_Pulse`.
- value: The value to which the intensifier modulation sequence parameter is to be set.

Output Parameters

There are no output parameters associated with `Picam_SetParameterModulationsValue()`

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `model`;
The `model` parameter may be set independently from the corresponding `device` parameter. However, doing so requires that all parameters be committed to the device prior to starting any data acquisition by calling `Picam_CommitParameters()`.
- `device`.
Setting a `device` parameter automatically sets the corresponding `model` parameter to the same value.

Related Structures

For additional information, refer to the following intensifier modulation sequence structures:

- `PicamModulation`;
- `PicamModulations`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_CommitParameters()`

5.3.1.22 *Picam_GetParameterIntegerDefaultValue()*

Description

Picam_GetParameterIntegerDefaultValue() returns the integer default value for a specified parameter.

Syntax

The syntax for *Picam_GetParameterIntegerDefaultValue()* is:

```
PICAM_API Picam_GetParameterIntegerDefaultValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    piint* value);
```

Input Parameters

Input parameters for *Picam_GetParameterIntegerDefaultValue()* are:

- camera: Handle for the camera for which the default parameter value is being requested.
- parameter: Specifies the parameter for which the integer default value is to be returned.
Valid parameters are those of type:
 - *PicamValueType_Integer*;
 - *PicamValueType_Boolean*;
 - *PicamValueType_Enumeration*.

Output Parameters

Output parameters for *Picam_GetParameterIntegerDefaultValue()* are:

- value: Pointer to the memory location in which the integer default value for the specified parameter has been stored.

Advanced API Usage

When used in conjunction with Advanced APIs, *camera* may be a handle to either the:

- *device*, or
- *model*.

Both the *device* and *model* share the same default value.

5.3.1.23 *Picam_GetParameterLargeIntegerDefaultValue()*

Description

`Picam_GetParameterLargeIntegerDefaultValue()` returns the large integer default value for a specified parameter.

Syntax

The syntax for `Picam_GetParameterLargeIntegerDefaultValue()` is:

```
PICAM_API Picam_GetParameterLargeIntegerDefaultValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    pi64s* value);
```

Input Parameters

Input parameters for `Picam_GetParameterLargeIntegerDefaultValue()` are:

- `camera`: Handle for the camera for which the default parameter value is being requested.
- `parameter`: Specifies the parameter for which the default value is to be returned. Valid parameter are those of type `PicamValueType_LargeInteger`.

Output Parameters

Output parameters for `Picam_GetParameterLargeIntegerDefaultValue()` are:

- `value`: Pointer to the memory location in which the large integer default value for the specified parameter has been stored.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

Both the `device` and `model` share the same default value.

5.3.1.24 *Picam_GetParameterFloatingPointDefaultValue()*

Description

`Picam_GetParameterFloatingPointDefaultValue()` returns the floating point default value for a specified parameter.

Syntax

The syntax for `Picam_GetParameterFloatingPointDefaultValue()` is:

```
PICAM_API Picam_GetParameterFloatingPointDefaultValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    piflt* value);
```

Input Parameters

Input parameters for `Picam_GetParameterFloatingPointDefaultValue()` are:

- `camera`: Handle for the camera for which the default parameter value is being requested.
- `parameter`: Specifies the parameter for which the default value is to be returned.
Valid parameters are those of type `PicamValueType_FloatingPoint`.

Output Parameters

Output parameters for `Picam_GetParameterFloatingPointDefaultValue()` are:

- `value`: Pointer to the memory location in which the floating point default value for the specified parameter has been stored.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

Both the `device` and `model` share the same default value.

5.3.1.25 *Picam_GetParameterRoisPointDefaultValue()*

Description

Picam_GetParameterRoisPointDefaultValue() returns the default value for a specified Rois parameter.

Syntax

The syntax for *Picam_GetParameterRoisPointDefaultValue()* is:

```
PICAM_API Picam_GetParameterRoisDefaultValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamRois** value);
```

Input Parameters

Input parameters for *Picam_GetParameterRoisPointDefaultValue()* are:

- camera: Handle for the camera for which the default parameter value is being requested.
- parameter: Specifies the Rois parameter for which the default value is to be returned.
Valid parameters are those of type *PicamValueType_Rois*.

Output Parameters

Output parameters for *Picam_GetParameterRoisPointDefaultValue()* are:

- value: Pointer to the memory location in which the default value for the specified Rois parameter has been stored.
NOTE: This memory is allocated by PICam and must be released by calling *Picam_DestroyRois()*.

Advanced API Usage

When used in conjunction with Advanced APIs, *camera* may be a handle to either the:

- *device*, or
- *model*.

Both the *device* and *model* share the same default value.

Related Structures

For additional information, refer to the following ROI structures:

- *PicamRoi*;
- *PicamRois*.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_DestroyRois()*.

5.3.1.26 *Picam_GetParameterPulseDefaultValue()*

Description

`Picam_GetParameterPulseDefaultValue()` returns the default value for a specified Pulse parameter.

Syntax

The syntax for `Picam_GetParameterPulseDefaultValue()` is:

```
PICAM_API Picam_GetParameterPulseDefaultValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamPulse** value);
```

Input Parameters

Input parameters for `Picam_GetParameterPulseDefaultValue()` are:

- `camera`: Handle for the camera for which the default parameter value is being requested.
- `parameter`: Specifies the Pulse parameter for which the default value is to be returned.
Valid parameters are those of type `PicamValueType_Pulse`.

Output Parameters

Output parameters for `Picam_GetParameterPulseDefaultValue()` are:

- `value`: Pointer to the memory location in which the default value for the specified Pulse parameter has been stored.

NOTE: This memory is allocated by PICam and must be released by calling `Picam_DestroyPulses()`.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

Both the `device` and `model` share the same default value.

Related Structures

For additional information, refer to the following Pulse structure:

- `PicamPulse`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_DestroyPulses()`.

5.3.1.27 *Picam_GetParameterModulationsDefaultValue()*

Description

Picam_GetParameterModulationsDefaultValue() returns the default value for a specified intensifier modulation sequence parameter.

Syntax

The syntax for *Picam_GetParameterModulationsDefaultValue()* is:

```
PICAM_API Picam_GetParameterModulationsDefaultValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamModulations** value);
```

Input Parameters

Input parameters for *Picam_GetParameterModulationsDefaultValue()* are:

- camera: Handle for the camera for which the default parameter value is being requested.
- parameter: Specifies the intensifier modulation sequence parameter for which the default value is to be returned.
Valid parameters are those of type *PicamValueType_Modulations*.

Output Parameters

Output parameters for *Picam_GetParameterModulationsDefaultValue()* are:

- value: Pointer to the memory location in which the default value for the specified intensifier modulation sequence parameter has been stored.
NOTE: This memory is allocated by PICam and must be released by calling *Picam_DestroyModulations()*.

Advanced API Usage

When used in conjunction with Advanced APIs, *camera* may be a handle to either the:

- *device*, or
- *model*.

Both the *device* and *model* share the same default value.

Related Structures

For additional information, refer to the following intensifier modulation sequence structures:

- *PicamModulation*;
- *PicamModulations*.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_DestroyModulations()*.

5.3.1.28 *Picam_RestoreParametersToDefaultValues()*

Description

`Picam_RestoreParametersToDefaultValues()` will set all read/write parameters to default values.

Syntax

The syntax for `Picam_RestoreParametersToDefaultValues()` is:

```
PICAM_API Picam_RestoreParametersToDefaultValues(  
                                                PicamHandle camera);
```

Input Parameters

Input parameters for `Picam_RestoreParametersToDefaultValues()` are:

`camera`: Handle for the camera for which parameters are to be restored.

Output Parameters

There are no output parameters associated with `Picam_RestoreParametersToDefaultValues()`

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

5.3.1.29 *Picam_CanSetParameterOnline()*

Description

`Picam_CanSetParameterOnline()` determines if the specified parameter can be configured during data acquisition.

Syntax

The syntax for `Picam_CanSetParameterOnline()` is:

```
PICAM_API Picam_CanSetParameterOnline(  
    PicamHandle camera,  
    PicamParameter parameter,  
    pibln* onlineable);
```

Input Parameters

Input parameters for `Picam_CanSetParameterOnline()` are:

- `camera`: Handle for the camera under test.
- `parameter`: Specifies the parameter for which the ability to be configured during data acquisition is to be determined.

Output Parameters

Output parameters for `Picam_CanSetParameterOnline()` are:

- `onlineable`: Pointer to the test results. Indicates if the specified parameter value can be set during data acquisition.
Valid values are:
 - `TRUE`
Indicates that the specified parameter can be configured during data acquisition.
 - `FALSE`
Indicates that the specified parameter cannot be configured during data acquisition.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

5.3.1.30 *Picam_SetParameterIntegerValueOnline()*

Description

`Picam_SetParameterIntegerValueOnline()` configures the specified parameter with an integer value during data acquisition.



NOTE:

The specified parameter must be capable of being configured during data acquisition.

Refer to `Picam_CanSetParameterOnline()` for additional information.

Syntax

The syntax for `Picam_SetParameterIntegerValueOnline()` is:

```
PICAM_API Picam_SetParameterIntegerValueOnline(  
    PicamHandle camera,  
    PicamParameter parameter,  
    piint value);
```

Input Parameters

Input parameters for `Picam_SetParameterIntegerValueOnline()` are:

`camera`: Handle for the camera being configured.

`parameter`: Specifies the parameter that is to be configured.

Valid parameters are those of type `PicamValueType_Integer`.

`value`: The integer value with which the specified parameter is to be configured.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

`Picam_SetParameterIntegerValueOnline()` effectively sets `parameter` on the camera device.

Output Parameters

There are no output parameters associated with

`Picam_SetParameterIntegerValueOnline()`

Related APIs

For additional information, refer to the following related APIs:

- `Picam_CommitParameters()`.

5.3.1.31 *Picam_SetParameterFloatingPointValueOnline()*

Description

`Picam_SetParameterFloatingPointValueOnline()` configures the specified parameter with a floating point value during data acquisition.



NOTE:

The specified parameter must be capable of being configured during data acquisition.

Refer to `Picam_CanSetParameterOnline()` for additional information.

Syntax

The syntax for `Picam_SetParameterFloatingPointValueOnline()` is:

```
PICAM_API Picam_SetParameterFloatingPointValueOnline(  
    PicamHandle camera,  
    PicamParameter parameter,  
    piflt value);
```

Input Parameters

Input parameters for `Picam_SetParameterFloatingPointValueOnline()` are:

- `camera`: Handle for the camera being configured.
- `parameter`: Specifies the parameter that is to be configured during data acquisition.
Valid parameters are those of type `PicamValueType_FloatingPoint`.
- `value`: The floating point value with which the specified parameter is to be configured.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

`Picam_SetParameterFloatingPointValueOnline()` effectively sets `parameter` on the camera device.

Output Parameters

There are no output parameters associated with

`Picam_SetParameterFloatingPointValueOnline()`

Related APIs

For additional information, refer to the following related APIs:

- `Picam_CommitParameters()`.

5.3.1.32 *Picam_SetParameterPulseValueOnline()*

Description

`Picam_SetParameterPulseValueOnline()` configures the specified Pulse parameter during data acquisition.



NOTE:

The specified parameter must be capable of being configured during data acquisition.

Refer to `Picam_CanSetParameterOnline()` for additional information.

Syntax

The syntax for `Picam_SetParameterPulseValueOnline()` is:

```
PICAM_API Picam_SetParameterPulseValueOnline(  
                                PicamHandle camera,  
                                PicamParameter parameter,  
                                const PicamPulse* value);
```

Input Parameters

Input parameters for `Picam_SetParameterPulseValueOnline()` are:

`camera`: Handle for the camera being configured.

`parameter`: Specifies the Pulse parameter that is to be configured during data acquisition.

Valid parameters are those of type `PicamValueType_Pulse`.

`value`: Pointer to the memory location in which the desired configuration value is stored.

Output Parameters

There are no output parameters associated with

`Picam_SetParameterPulseValueOnline()`.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

`Picam_SetParameterPulseValueOnline()` effectively sets `parameter` on the camera device.

5.3.1.33 *Picam_CanReadParameter()*

Description

`Picam_CanReadParameter()` determines if a parameter value can be read directly from hardware connected to the system.

Syntax

The syntax for `Picam_CanReadParameter()` is:

```
PICAM_API Picam_CanReadParameter(  
    PicamHandle camera,  
    PicamParameter parameter,  
    pibln* readable);
```

Input Parameters

Input parameters for `Picam_CanReadParameter()` are:

- camera: Handle for the camera under test.
- parameter: Specifies the parameter for which the ability to read its value directly from the hardware is to be determined.

Output Parameters

Output parameters for `Picam_CanReadParameter()` are:

- readable: Pointer to the test results. Indicates if the specified parameter value can be read directly from the hardware.
Valid values are:
 - TRUE
Indicates that the value for the specified parameter can be read from the hardware.
 - FALSE
Indicates that the value for the specified parameter cannot be read from the hardware.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- device, or
- model.

5.3.1.34 *Picam_ReadParameterIntegerValue()*

Description

Picam_ReadParameterIntegerValue() returns the integer value for a specified parameter as read directly from hardware connected to the system.



NOTE:

The specified parameter must be capable of being read directly from the hardware.

Refer to *Picam_CanReadParameter()* for additional information.

Syntax

The syntax for *Picam_ReadParameterIntegerValue()* is:

```
PICAM_API Picam_ReadParameterIntegerValue(
                                PicamHandle camera,
                                PicamParameter parameter,
                                piint* value);
```

Input Parameters

Input parameters for *Picam_ReadParameterIntegerValue()* are:

camera: Handle for the camera under test.

parameter: Specifies the parameter that is to have its value read from hardware.

NOTE: The specified parameter must be capable of being read directly from hardware.

Refer to *Picam_CanReadParameter()* for additional information.

Valid parameters are those of type:

- *PicamValueType_Integer;*
- *PicamValueType_Boolean;*
- *PicamValueType_Enumeration.*

Output Parameters

Output parameters for *Picam_ReadParameterIntegerValue()* are:

value: Pointer to the memory location in which the parameter value is stored.

Advanced API Usage

When used in conjunction with Advanced APIs, *camera* may be a handle to either the:

- *device*, or
- *model*.

Picam_ReadParameterIntegerValue() effectively gets *parameter* from the camera device.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_CanReadParameter()*.

5.3.1.35 *Picam_ReadParameterFloatingPointValue()*

Description

`Picam_ReadParameterFloatingPointValue()` returns the floating point value for a specified parameter as read directly from hardware connected to the system.



NOTE:

The specified parameter must be capable of being read directly from the hardware.

Refer to `Picam_CanReadParameter()` for additional information.

Syntax

The syntax for `Picam_ReadParameterFloatingPointValue()` is:

```
PICAM_API Picam_ReadParameterFloatingPointValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    piflt* value);
```

Input Parameters

Input parameters for `Picam_ReadParameterFloatingPointValue()` are:

`camera`: Handle for the camera under test.

`parameter`: Specifies the parameter that is to have its value read from hardware.

NOTE: The specified parameter must be capable of being read directly from hardware.

Refer to `Picam_CanReadParameter()` for additional information.

Valid parameter are those of type:

- `PicamValueType_FloatingPoint`.

Output Parameters

Output parameters for `Picam_ReadParameterFloatingPointValue()` are:

`value`: Pointer to the memory location in which the parameter value is stored.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

`Picam_ReadParameterFloatingPointValue()` effectively gets `parameter` from the camera device.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_CanReadParameter()`.

5.3.1.36 *Picam_CanWaitForStatusParameter()*

Description

`Picam_CanWaitForStatusParameter()` determines if a parameter is a waitable status.

Syntax

The syntax for `Picam_CanWaitForStatusParameter()` is:

```
PICAM_API Picam_CanWaitForStatusParameter(  
    PicamHandle camera,  
    PicamParameter parameter,  
    pibln* waitable);
```

Input Parameters

Input parameters for `Picam_CanWaitForStatusParameter()` are:

camera: Handle for the camera under test.

parameter: Specifies the parameter to check as a waitable status.

Output Parameters

Output parameters for `Picam_CanWaitForStatusParameter()` are:

waitable: Pointer to the test results. Indicates if the specified parameter is a waitable status.

Valid values are:

- TRUE
Indicates that the parameter is a waitable status.
- FALSE
Indicates that the parameter is not a waitable status.

Advanced API Usage

When used in conjunction with Advanced APIs, camera may be a handle to either the:

- device, or
- model.

5.3.1.37 *Picam_DestroyStatusPurviews()*

Description

`Picam_DestroyStatusPurviews()` releases memory that has been allocated by PICam for use by the `purviews_array`.

If the `purviews_array` is null, calling `Picam_DestroyStatusPurviews()` has no effect.

Syntax

The syntax for `Picam_DestroyStatusPurviews()` is:

```
PICAM_API Picam_DestroyStatusPurviews(  
    const PicamStatusPurview* purviews_array);
```

Input Parameters

Input parameters for `Picam_DestroyStatusPurviews()` are:

`purviews_array`: Pointer to array memory that is to be released.

Output Parameters

There are no output parameters associated with `Picam_DestroyStatusPurviews()`.

5.3.1.38 *Picam_GetStatusParameterPurview()*

Description

Picam_GetStatusParameterPurview() returns the scope of a waitable status.

Syntax

The syntax for *Picam_GetStatusParameterPurview()* is:

```
PICAM_API Picam_GetStatusParameterPurview(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamStatusPurview** purview);
```

Input Parameters

Input parameters for *Picam_GetStatusParameterPurview()* are:

camera: Handle for the camera for which the status purview is being requested.

parameter: Specifies the parameter whose status purview is being requested.

NOTE: The specified parameter must be a waitable status.

Refer to *Picam_CanWaitForStatusParameter()* for additional information.

Output Parameters

Output parameters for *Picam_GetStatusParameterPurview()* are:

purview: Pointer to the allocated status purview.

NOTE: This memory is allocated by PICam and must be released by calling *Picam_DestroyStatusPurviews()*.

Advanced API Usage

When used in conjunction with Advanced APIs, *camera* may be a handle to either the:

- *device*, or
- *model*.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_CanWaitForStatusParameter()*
- *Picam_DestroyStatusPurviews()*

5.3.1.39 *Picam_EstimateTimeToStatusParameterValue()*

Description

[Picam_EstimateTimeToStatusParameterValue\(\)](#) returns the estimated time, in milliseconds, for a particular status to be reached.

Syntax

The syntax for [Picam_EstimateTimeToStatusParameterValue\(\)](#) is:

```
PICAM_API Picam_EstimateTimeToStatusParameterValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    piint value,  
    piint* estimated_time);
```

Input Parameters

Input parameters for [Picam_EstimateTimeToStatusParameterValue\(\)](#) are:

camera: Handle for the camera whose time to status will be estimated.

parameter: Handle for the camera whose time to status will be estimated.

NOTE: The specified parameter must be a waitable status.

Refer to [Picam_CanWaitForStatusParameter\(\)](#) for additional information.

value: Specifies the status for which the time is to be estimated.

NOTE: The specified value must be in the status purview.

Refer to [Picam_GetStatusParameterPurview\(\)](#) for additional information.

Output Parameters

Output parameters for [Picam_EstimateTimeToStatusParameterValue\(\)](#) are:

estimated_time: Pointer to the estimated time in milliseconds.

NOTE: If the time cannot be estimated, -1 is returned.

Advanced API Usage

When used in conjunction with Advanced APIs, camera may be a handle to either the:

- device, or
- model.

Related APIs

For additional information, refer to the following related APIs:

- [Picam_CanWaitForStatusParameter\(\)](#)
- [Picam_GetStatusParameterPurview\(\)](#)

5.3.1.40 *Picam_WaitForStatusParameterValue()*

Description

Picam_WaitForStatusParameterValue() waits for a particular status to be reached or until *time_out* milliseconds has elapsed.

PicamError_TimeOutOccurred is returned if *time_out* has elapsed.

Syntax

The syntax for *Picam_WaitForStatusParameterValue()* is:

```
PICAM_API Picam_WaitForStatusParameterValue(  
    PicamHandle camera,  
    PicamParameter parameter,  
    piint value,  
    piint time_out);
```

Input Parameters

Input parameters for *Picam_WaitForStatusParameterValue()* are:

camera: Handle for the camera whose status will be awaited.

parameter: Specifies the parameter whose status will be awaited.

NOTE: The specified parameter must be a waitable status.

Refer to *Picam_CanWaitForStatusParameter()* for additional information.

value: Specifies the status to await.

NOTE: The specified value must be in the status purview.

Refer to *Picam_GetStatusParameterPurview()* for additional information.

time_out: Specifies the time to wait, in milliseconds.

NOTE: Use -1 to wait indefinitely.

Output Parameters

There are no output parameters associated with *Picam_WaitForStatusParameterValue()*.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_CanWaitForStatusParameter()*
- *Picam_GetStatusParameterPurview()*

5.3.2 Camera Parameter Information APIs

This section provides programming information for APIs used to configure and retrieve Camera Parameter information.

5.3.2.1 *Picam_DestroyParameters()*

Description

`Picam_DestroyParameters()` releases memory that has been allocated by PICam for use by `parameter_array`.

If `parameter_array` is null, calling `Picam_DestroyParameters()` has no effect.



NOTE:

`parameter_array` may be a single `PicamParameter` allocated by PICam.

Syntax

The syntax for `Picam_DestroyParameters()` is:

```
PICAM_API Picam_DestroyParameters(  
    const PicamParameter* parameter_array);
```

Input Parameters

Input parameters for `Picam_DestroyParameters()` are:

`parameter_array`: Pointer to array memory that is to be released.

Output Parameters

There are no output parameters associated with `Picam_DestroyParameters()`.

Related Structures

For additional information, refer to the following parameter structure:

- `PicamParameter`.

5.3.2.2 *Picam_GetParameters()*

Description

Picam_GetParameters() returns a list of parameters that are available for a specified camera. The number of parameters is also returned.

Syntax

The syntax for *Picam_GetParameters()* is:

```
PICAM_API Picam_GetParameters(  
    PicamHandle camera,  
    const PicamParameter** parameter_array,  
    piint* parameter_count);
```

Input Parameters

Input parameters for *Picam_GetParameters()* are:

camera: Handle for the camera under test.

Output Parameters

Output parameters for *Picam_GetParameters()* are:

parameter_array: Pointer to the allocated array in which the list of parameters associated with the specified camera is stored.

NOTE: This memory is allocated by PICam and must be released by calling *Picam_DestroyParameters()*.

parameter_count: Pointer to the memory location in which the number of available parameters associated with the specified camera is stored.

Advanced API Usage

When used in conjunction with Advanced APIs, camera may be a handle to either the:

- device, or
- model.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_DestroyParameters()*.

Related Structures

For additional information, refer to the following parameter structure:

- *PicamParameter*.

5.3.2.3 *Picam_DoesParameterExist()*

Description

`Picam_DoesParameterExist()` determines if a specified parameter is available for a specified camera.

Syntax

The syntax for `Picam_DoesParameterExist()` is:

```
PICAM_API Picam_DoesParameterExist(  
    PicamHandle camera,  
    PicamParameter parameter,  
    pibln* exists);
```

Input Parameters

Input parameters for `Picam_DoesParameterExist()` are:

camera: Handle for the camera under test.

parameter: Specifies the parameter for which availability is being determined.

Output Parameters

Output parameters for `Picam_DoesParameterExist()` are:

exists: Pointer to the test results. Indicates if the specified parameter is available on the specified camera.

Valid values are:

- TRUE
Indicates that the specified parameter is available on the specified camera.
- FALSE
Indicates that the specified parameter is not available on the specified camera.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- device, or
- model.

Related Structures

For additional information, refer to the following parameter structure:

- `PicamParameter`.

5.3.2.4 *Picam_IsParameterRelevant()*

Description

`Picam_IsParameterRelevant()` determines if the value of a specified parameter is currently applicable for a specified camera.

Syntax

The syntax for `Picam_IsParameterRelevant()` is:

```
PICAM_API Picam_IsParameterRelevant(  
    PicamHandle camera,  
    PicamParameter parameter,  
    pibln* relevant);
```

Input Parameters

Input parameters for `Picam_IsParameterRelevant()` are:

`camera`: Handle for the camera under test.

`parameter`: Specifies the parameter for which value applicability is being determined.

Output Parameters

Output parameters for `Picam_IsParameterRelevant()` are:

`relevant`: Pointer to the test results. Indicates if the specified parameter value is currently applicable for the specified camera.

Valid values are:

- `TRUE`
Indicates that the specified parameter value is currently applicable for the specified camera.
- `FALSE`
Indicates that the specified parameter value is not currently applicable for the specified camera.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

Related Structures

For additional information, refer to the following parameter structure:

- `PicamParameter`.

5.3.2.5 *Picam_GetParameterValueType()*

Description

`Picam_GetParameterValueType()` returns the data type for a value stored within a specified parameter.

Syntax

The syntax for `Picam_GetParameterValueType()` is:

```
PICAM_API Picam_GetParameterValueType(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamValueType* type);
```

Input Parameters

Input parameters for `Picam_GetParameterValueType()` are:

`camera`: Handle for the camera under test.

`parameter`: Specifies the parameter for which the data type of the stored value is being requested.

Output Parameters

Output parameters for `Picam_GetParameterValueType()` are:

`type`: Pointer to the memory location in which the data type of the specified parameter's value is stored.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

Related Structures

For additional information, refer to the following parameter structure:

- `PicamParameter`;
- `PicamValueType`.

5.3.2.6 *Picam_GetParameterEnumeratedType()*

Description

Picam_GetParameterEnumeratedType() returns the enumeration type for a specified parameter.

Syntax

The syntax for *Picam_GetParameterEnumeratedType()* is:

```
PICAM_API Picam_GetParameterEnumeratedType(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamEnumeratedType* type);
```

Input Parameters

Input parameters for *Picam_GetParameterEnumeratedType()* are:

camera: Handle for the camera under test.

parameter: Specifies the parameter for which the enumeration type is being requested.
Valid parameters are those of type *PicamValueType_Enumeration*.

Output Parameters

Output parameters for *Picam_GetParameterEnumeratedType()* are:

type: Pointer to the memory location in which the enumeration type of the specified parameter is stored.

Advanced API Usage

When used in conjunction with Advanced APIs, *camera* may be a handle to either the:

- *device*, or
- *model*.

Related Structures

For additional information, refer to the following parameter structure:

- *PicamParameter*.

5.3.2.7 *Picam_GetParameterValueAccess()*

Description

`Picam_GetParameterValueAccess()` returns the read/write permissions for the specified parameter.

Syntax

The syntax for `Picam_GetParameterValueAccess()` is:

```
PICAM_API Picam_GetParameterValueAccess(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamValueAccess* access);
```

Input Parameters

Input parameters for `Picam_GetParameterValueAccess()` are:

- `camera`: Handle for the camera under test.
- `parameter`: Specifies the parameter for which read/write permission is being requested.

Output Parameters

Output parameters for `Picam_GetParameterValueAccess()` are:

- `access`: Pointer to the memory location in which the read/write permission for the specified parameter is stored.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

Related Structures

For additional information, refer to the following parameter structure:

- `PicamParameter`.

5.3.2.8 *Picam_GetParameterConstraintType()*

Description

`Picam_GetParameterConstraintType()` returns the type of constraint placed on a specified parameter.

Syntax

The syntax for `Picam_GetParameterConstraintType()` is:

```
PICAM_API Picam_GetParameterConstraintType(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamConstraintType* type);
```

Input Parameters

Input parameters for `Picam_GetParameterConstraintType()` are:

`camera`: Handle for the camera under test.

`parameter`: Specifies the parameter for which constraint information is being requested.

Output Parameters

Output parameters for `Picam_GetParameterConstraintType()` are:

`type`: Pointer to the memory location in which constraint information for the specified parameter is stored.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

Related Structures

For additional information, refer to the following parameter structure:

- `PicamParameter`.

5.3.3 Camera Parameter Constraints APIs

This section provides programming information for APIs used to configure camera parameter constraints.

5.3.3.1 *Picam_DestroyCollectionConstraints()*

Description

`Picam_DestroyCollectionConstraints()` releases memory that has been allocated by PICam for use by `constraint_array`.

If `constraint_array` is null, calling `Picam_DestroyCollectionConstraints()` has no effect.



NOTE:

`constraint_array` may be a single `PicamCollectionConstraint` allocated by PICam.

Syntax

The syntax for `Picam_DestroyCollectionConstraints()` is:

```
PICAM_API Picam_DestroyCollectionConstraints(  
    const PicamCollectionConstraint* constraint_array);
```

Input Parameters

Input parameters for `Picam_DestroyCollectionConstraints()` are:

`constraint_array`: Pointer to array memory that is to be released.

Output Parameters

There are no output parameters associated with `Picam_DestroyCollectionConstraints()`.

Related Structures

For additional information, refer to the following parameter structure:

- `PicamCollectionConstraint`.

5.3.3.2 *Picam_GetParameterCollectionConstraint()*

Description

Picam_GetParameterCollectionConstraint() returns constraint information for a specified constraint category and parameter combination.

Syntax

The syntax for *Picam_GetParameterCollectionConstraint()* is:

```
PICAM_API Picam_GetParameterCollectionConstraint(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamConstraintCategory category,  
    const PicamCollectionConstraint** constraint);
```

Input Parameters

Input parameters for *Picam_GetParameterCollectionConstraint()* are:

- camera: Handle for the camera for which constraint information is being returned.
- parameter: Specifies the parameter for which constraint information is being requested.
- category: Specifies the constraint category for which the list of constraints is being requested.

Output Parameters

Output parameters for *Picam_GetParameterCollectionConstraint()* are:

- constraint: Pointer to the allocated array in which the list of constraints available for the specified constraint category and parameter combination are stored.

NOTE: This memory is allocated by PICam and must be released by calling *Picam_DestroyCollectionConstraints()*.

Advanced API Usage

When used in conjunction with Advanced APIs, camera may be a handle to either the:

- device, or
- model.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_DestroyCollectionConstraints()*.

Related Structures

For additional information, refer to the following parameter structures:

- *PicamParameter*;
- *PicamCollectionConstraint*.

5.3.3.3 *Picam_DestroyRangeConstraints()*

Description

`Picam_DestroyRangeConstraints()` releases memory that has been allocated by PICam for use by `constraint_array`.

If `constraint_array` is null, calling `Picam_DestroyRangeConstraints()` has no effect.



NOTE:

`constraint_array` may be a single `PicamRangeConstraint` allocated by PICam.

Syntax

The syntax for `Picam_DestroyRangeConstraints()` is:

```
PICAM_API Picam_DestroyRangeConstraints(  
    const PicamRangeConstraint* constraint_array);
```

Input Parameters

Input parameters for `Picam_DestroyRangeConstraints()` are:

`constraint_array`: Pointer to array memory that is to be released.

Output Parameters

There are no output parameters associated with `Picam_DestroyRangeConstraints()`.

Related Structures

For additional information, refer to the following parameter structure:

- `PicamRangeConstraint`.

5.3.3.4 *Picam_GetParameterRangeConstraint()*

Description

`Picam_GetParameterRangeConstraint()` returns range constraints for a specified constraint category and parameter combination.

Syntax

The syntax for `Picam_GetParameterRangeConstraint()` is:

```
PICAM_API Picam_GetParameterRangeConstraint(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamConstraintCategory category,  
    const PicamRangeConstraint** constraint);
```

Input Parameters

Input parameters for `Picam_GetParameterRangeConstraint()` are:

- `camera`: Handle for the camera for which range constraints are being returned.
- `parameter`: Specifies the parameter for which range constraint information is being requested.
- `category`: Specifies the constraint category for which range constraint information is being requested.

Output Parameters

Output parameters for `Picam_GetParameterRangeConstraint()` are:

- `constraint`: Pointer to the allocated array in which the range constraints for the specified constraint category and parameter combination are stored.
NOTE: This memory is allocated by PICam and must be released by calling `Picam_DestroyRangeConstraints()`.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_DestroyRangeConstraints()`.

Related Structures

For additional information, refer to the following parameter structures:

- `PicamParameter`;
- `PicamRangeConstraint`.

5.3.3.5 *Picam_DestroyRoisConstraints()*

Description

`Picam_DestroyRoisConstraints()` releases memory that has been allocated by PICam for use by `constraint_array`.

If `constraint_array` is null, calling `Picam_DestroyRoisConstraints()` has no effect.



NOTE:

`constraint_array` may be a single `PicamRoisConstraint` allocated by PICam.

Syntax

The syntax for `Picam_DestroyRoisConstraints()` is:

```
PICAM_API Picam_DestroyRoisConstraints(  
    const PicamRoisConstraint* constraint_array);
```

Input Parameters

Input parameters for `Picam_DestroyRoisConstraints()` are:

`constraint_array`: Pointer to array memory that is to be released.

Output Parameters

There are no output parameters associated with `Picam_DestroyRoisConstraints()`.

Related Structures

For additional information, refer to the following parameter structure:

- `PicamRoisConstraint`.

5.3.3.6 *Picam_GetParameterRoisConstraint()*

Description

Picam_GetParameterRoisConstraint() returns Roi constraints for a specified constraint category and parameter combination.

Syntax

The syntax for *Picam_GetParameterRoisConstraint()* is:

```
PICAM_API Picam_GetParameterRoisConstraint(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamConstraintCategory category,  
    const PicamRoisConstraint** constraint);
```

Input Parameters

Input parameters for *Picam_GetParameterRoisConstraint()* are:

- camera: Handle for the camera for which constraint information is being returned.
- parameter: Specifies the parameter for which Rois constraint information is being requested.
- category: Specifies the constraint category for which Roi constraint information is being requested.

Output Parameters

Output parameters for *Picam_GetParameterRoisConstraint()* are:

- constraint: Pointer to the allocated array in which the Rois constraints for the specified constraint category and parameter combination are stored.
NOTE: This memory is allocated by PICam and must be released by calling *Picam_DestroyRoisConstraints()*.

Advanced API Usage

When used in conjunction with Advanced APIs, camera may be a handle to either the:

- device, or
- model.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_DestroyRoisConstraints()*.

Related Structures

For additional information, refer to the following parameter structures:

- *PicamParameter*;
- *PicamRoisConstraint*.

5.3.3.7 *Picam_DestroyPulseConstraints()*

Description

`Picam_DestroyPulseConstraints()` releases memory that has been allocated by PICam for use by `constraint_array`.

If `constraint_array` is null, calling `Picam_DestroyPulseConstraints()` has no effect.



NOTE:

`constraint_array` may be a single `PicamPulseConstraint` allocated by PICam.

Syntax

The syntax for `Picam_DestroyPulseConstraints()` is:

```
PICAM_API Picam_DestroyPulseConstraints(  
    const PicamPulseConstraint* constraint_array);
```

Input Parameters

Input parameters for `Picam_DestroyPulseConstraints()` are:

`constraint_array`: Pointer to array memory that is to be released.

Output Parameters

There are no output parameters associated with `Picam_DestroyPulseConstraints()`.

Related Structures

For additional information, refer to the following parameter structure:

- `PicamPulseConstraint`.

5.3.3.8 *Picam_GetParameterPulseConstraint()*

Description

Picam_GetParameterPulseConstraint() returns Pulse constraints for a specified constraint category and parameter combination.

Syntax

The syntax for *Picam_GetParameterPulseConstraint()* is:

```
PICAM_API Picam_GetParameterPulseConstraint(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamConstraintCategory category,  
    const PicamPulseConstraint** constraint);
```

Input Parameters

Input parameters for *Picam_GetParameterPulseConstraint()* are:

- camera: Handle for the camera for which constraint information is being returned.
- parameter: Specifies the parameter for which Pulse constraint information is being requested.
- category: Specifies the constraint category for which Pulse constraint information is being requested.

Output Parameters

Output parameters for *Picam_GetParameterPulseConstraint()* are:

- constraint: Pointer to the allocated array in which the Pulse constraints for the specified constraint category and parameter combination are stored.
NOTE: This memory is allocated by PICam and must be released by calling *Picam_DestroyPulseConstraints()*.

Advanced API Usage

When used in conjunction with Advanced APIs, camera may be a handle to either the:

- device, or
- model.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_DestroyPulseConstraints()*.

Related Structures

For additional information, refer to the following parameter structures:

- *PicamParameter*;
- *PicamPulseConstraint*.

5.3.3.9 *Picam_DestroyModulationsConstraints()*

Description

`Picam_DestroyModulationsConstraints()` releases memory that has been allocated by PICam for use by `constraint_array`.

If `constraint_array` is null, calling `Picam_DestroyModulationsConstraints()` has no effect.



NOTE:

`constraint_array` may be a single `PicamModulationsConstraint` allocated by PICam.

Syntax

The syntax for `Picam_DestroyModulationsConstraints()` is:

```
PICAM_API Picam_DestroyModulationConstraints(  
    const PicamModulationsConstraint* constraint_array);
```

Input Parameters

Input parameters for `Picam_DestroyModulationsConstraints()` are:

`constraint_array`: Pointer to array memory that is to be released.

Output Parameters

There are no output parameters associated with `Picam_DestroyModulationsConstraints()`.

Related Structures

For additional information, refer to the following parameter structure:

- `PicamModulationsConstraint`.

5.3.3.10 *Picam_GetParameterModulationsConstraint()*

Description

Picam_GetParameterModulationsConstraint() returns intensifier modulation sequence constraints for a specified constraint category and parameter combination.

Syntax

The syntax for *Picam_GetParameterModulationsConstraint()* is:

```
PICAM_API Picam_GetParameterModulationsConstraint(
                                PicamHandle camera,
                                PicamParameter parameter,
                                PicamConstraintCategory category,
                                const PicamModulationsConstraint** constraint);
```

Input Parameters

Input parameters for *Picam_GetParameterModulationsConstraint()* are:

- camera: Handle for the camera for which constraint information is being returned.
- parameter: Specifies the parameter for which intensifier modulation sequence constraint information is being requested.
- category: Specifies the constraint category for which intensifier modulation sequence constraint information is being requested.

Output Parameters

Output parameters for *Picam_GetParameterModulationsConstraint()* are:

- constraint: Pointer to the allocated array in which the intensifier modulation sequence constraints for the specified constraint category and parameter combination are stored.

NOTE: This memory is allocated by PICam and must be released by calling *Picam_DestroyModulationsConstraints()*.

Advanced API Usage

When used in conjunction with Advanced APIs, camera may be a handle to either the:

- device, or
- model.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_DestroyModulationsConstraints()*.

Related Structures

For additional information, refer to the following parameter structures:

- *PicamParameter*;
- *PicamModulationsConstraint*.

5.3.4 Camera Parameter Commitment APIs

This section provides programming information about APIs used to commit parameter values.

5.3.4.1 *Picam_AreParametersCommitted()*

Description

`Picam_AreParametersCommitted()` determines if the parameter configuration changes have been applied to the specified camera.

Syntax

The syntax for `Picam_AreParametersCommitted()` is:

```
PICAM_API Picam_AreParametersCommitted(  
                                PicamHandle camera,  
                                pibln* committed);
```

Input Parameters

Input parameters for `Picam_AreParametersCommitted()` are:

`camera`: Handle for the camera for which parameter configuration status information is being determined.

Output Parameters

Output parameters for `Picam_AreParametersCommitted()` are:

`committed`: Pointer to the test results. Indicates if parameter configuration changes have been committed to the specified camera.

Valid values are:

- TRUE
Indicates that parameter configuration changes have been committed to the specified camera.
- FALSE
Indicates that parameter configuration changes have not been committed to the specified camera.

Advanced API Usage

When used in conjunction with Advanced APIs, `camera` may be a handle to either the:

- `device`, or
- `model`.

5.3.4.2 *Picam_CommitParameters()*

Description

Picam_CommitParameters() validates parameter values and applies these valid values to the specified camera during system setup and configuration.

- Any parameter that fails to satisfy its required constraint(s) is flagged as invalid and is stored within *failed_parameter_array*.
- The number of invalid parameters is stored in *failed_parameter_count*. If no invalid parameters are detected, this value is 0.

Syntax

The syntax for *Picam_CommitParameters()* is:

```
PICAM_API Picam_CommitParameters(
    PicamHandle camera,
    const PicamParameter** failed_parameter_array,
    piint* failed_parameter_count);
```

Input Parameters

Input parameters for *Picam_CommitParameters()* are:

camera: Handle for the camera for which parameter values are being configured.

Output Parameters

Output parameters for *Picam_CommitParameters()* are:

failed_parameter_array: Pointer to the allocated array in which the list of failed/invalid parameters is stored.

If no invalid parameters are detected, this is a null object.

NOTE: This memory is allocated by PICam and must be released by calling *Picam_DestroyParameters()*.

failed_parameter_count: Pointer to the memory location in which the number of failed/invalid parameters is stored.

Advanced API Usage

When used in conjunction with Advanced APIs, *camera* may be a handle to either the:

- *device*, or
- *model*.

Picam_CommitParameters() systematically configures device with (valid) parameter values that have been stored in *model*.

Chapter 6: Camera Data Acquisition APIs

Once a camera has been configured and the parameters are committed, the camera is ready to acquire data. Data can be acquired either synchronously or asynchronously.

By default, memory is allocated automatically to accommodate the data. This automatic memory is valid until the next acquisition or until the camera is closed.

By default, the data are returned as follows:

- One frame of sensor data containing each region of interest (in the order defined);
- Followed by any metadata for that frame (timestamps followed by frame tracking, gate tracking delay, gate tracking width, and modulation tracking);
- Repeated for each frame in one readout;
- Possibly followed by any padding between readouts.

Configuring the camera such that the total number of readouts is indeterminate will disable automatic data memory management:

- Basic:
 - Instruct the camera to acquire more data than it can exactly acquire;
This is achieved by setting `PicamParameter_ReadoutCount` to a value greater than the value of `PicamParameter_ExactReadoutCountMaximum`.
 - Instruct the camera to readout data non-destructively (for cameras that have this feature)
- Advanced:
 - Instruct the camera to acquire data indefinitely
This is achieved by setting `PicamParameter_ReadoutCount` to 0.

Also, setting a user-allocated buffer with `PicamAdvanced_SetAcquisitionBuffer()` will disable automatic data memory management.

6.1 Data Format


Acquired data is structured as shown in [Figure 6-1](#).



NOTE:

All partitions are specified in bytes.

Figure 6-1: Data Format Diagram



The camera acquires N readouts each a readout stride apart. One readout stride is comprised of M frames each a frame stride apart followed by padding. A frame stride is divided into frame pixel data and frame metadata by a frame size. Frame pixel data contains data for X regions of interest; whose regions are in the order each region was defined. Frame metadata contains any time stamps followed by the frame tracking number, gate tracking (delay and width), and modulation tracking (duration, frequency, phase, and output signal frequency). All formatting information is available as read-only data acquisition parameters.

6.2 Data Type Definitions

This section provides programming information about PICam data definitions.

6.2.1 Data Acquisition Enumerations

This section provides detailed information about the following data acquisition enumerations:

- `PicamAcquisitionErrorsMask`.

6.2.1.1 *PicamAcquisitionErrorsMask*

Data Type

`PicamAcquisitionErrorsMask` is defined as enum.

Description

`PicamAcquisitionErrorsMask` is the set of acquisition error messages.

Enumerator Definitions

Refer to [Table 6-1](#) for enumerator definitions.

Table 6-1: `PicamAcquisitionErrorsMask` Enumerator Definitions

Enumerator	Description
<code>PicamAcquisitionErrorsMask_CameraFaulted</code>	The camera has critically malfunctioned and is in need of service. Further acquisitions are not possible until the camera has been serviced.
<code>PicamAcquisitionErrorsMask_ConnectionLost</code>	The camera was disconnected.
<code>PicamAcquisitionErrorsMask_DataLost</code>	Data has been lost.
<code>PicamAcquisitionErrorsMask_DataNotArriving</code>	Data is no longer arriving from the camera.
<code>PicamAcquisitionErrorsMask_None</code>	No errors have occurred.
<code>PicamAcquisitionErrorsMask_ShutterOverheated</code>	A connected shutter has overheated and is temporarily disabled. Further acquisitions are not possible until the shutter is no longer overheated.

6.3 Data Acquisition Data Structures

This section provides programming information about the following PICam data acquisition data structures:

- `PicamAvailableData`;
- `PicamAcquisitionStatus`.

6.3.1 `PicamAvailableData`

Description

`PicamAvailableData` represents newly acquired data.

Structure Definition

The structure definition for `PicamAvailableData` is:

```
typedef struct PicamAvailableData
{
    void*   initial_readout;
    pi64s   readout_count;
} PicamAvailableData;
```

Variable Definitions

The variables required by `PicamAvailableData` are:

<code>initial_readout</code> :	Pointer to the start of the first available readout.
<code>readout_count</code> :	Indicates how many contiguous readouts are currently available.

6.3.2 PicamAcquisitionStatus

Description

`PicamAcquisitionStatus` reports various status information during data acquisition by the camera.

Structure Definition

The structure definition for `PicamAcquisitionStatus` is:

```
typedef struct PicamAcquisitionStatus
{
    pibln    running;
    PicamAcquisitionErrorsMask errors;
    piflt    readout_rate;
} PicamAcquisitionStatus;
```

Variable Definitions

The variables required by `PicamAcquisitionStatus` are:

`running`: Indicates the data acquisition status/

Valid values are:

- `TRUE`
Indicates an acquisition is in progress.
- `FALSE`
Indicates there is no current data acquisition in progress.

`errors`: Contains any errors that have occurred.

`readout_rate`: The rate of capture in readouts-per-second when acquiring more than one readout.

6.4 Programmers' Reference for Acquisition Control APIs

This section provides programming information for the following acquisition control APIs:

- `Picam_Acquire()`;
- `Picam_StartAcquisition()`;
- `Picam_StopAcquisition()`;
- `Picam_IsAcquisitionRunning()`;
- `Picam_WaitForAcquisitionUpdate()`.

6.4.1 `Picam_Acquire()`

Description

`Picam_Acquire()` performs a specified number of data readouts (specified by `readout_count`) and returns once the acquisition has been completed.



NOTE:

This function cannot be called when cameras are configured for non-destructively readout. This is because the number of readouts acquired is no longer guaranteed to be fixed. As an example, changing the exposure time online will change the number of non-destructive readouts and therefore the total number of readouts acquired.



NOTE:

Parameters must be committed prior to initiating data acquisition. Refer to [Section 5.3.4.2](#), `Picam_CommitParameters()`, on page 156 for additional information.

Data acquisition is successful when:

- The delay between successive readouts does not exceed `readout_time_out`, and
- No errors have occurred.

Data acquisition is immediately halted when:

- The delay between successive readouts exceeds that specified by `readout_time_out`.

The error message `PicamError_TimeOutOccurred` is returned.

- Any other error conditions are detected.
Associated error messages are stored in the `errors` parameter.

Syntax

The syntax for `Picam_Acquire()` is:

```
PICAM_API Picam_Acquire(
    PicamHandle camera,
    pi64s readout_count,
    piint readout_time_out,
    PicamAvailableData* available,
    PicamAcquisitionErrorsMask* errors);
```

Input Parameters

Input parameters for `Picam_Acquire()` are:

- `camera`: Handle for the camera from which data are to be acquired.
- `readout_count`: The number of readouts desired.
Valid values are in the range:
[1...`PicamParameter_ExactReadoutCountMaximum`]
If this value becomes excessively large, this function may fail due to a lack of sufficient memory.
- `readout_time_out`: The time, in mS, to wait between each successive readout.
When specifying an infinite length of time, configure this parameter to **-1**.

Output Parameters

Output parameters for `Picam_Acquire()` are:

- `available`: The output buffer used to store data that has been successfully read out from the specified camera.
In the event of a data acquisition failure, this buffer may contain little to no data.
Data stored in this buffer is valid until:
- The next acquisition cycle is initiated; or
 - The camera is closed.
- `errors`: The parameter used to store any error messages that were raised during data acquisition.

Advanced API Usage

When used in conjunction with Advanced APIs, data in the output buffer `available` is also invalidated when `PicamAdvanced_SetAcquisitionBuffer()` is called.

`Picam_Acquire()` is mutually exclusive with the use of an acquisition-updated callback.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_CommitParameters()`;
- `PicamAdvanced_SetAcquisitionBuffer()`.

6.4.2 Picam_StartAcquisition()

Description

`Picam_StartAcquisition()` asynchronously initiates data acquisition and returns immediately.



NOTE:

Parameters must be committed prior to initiating data acquisition. Refer to [Section 5.3.4.2, `Picam_CommitParameters\(\)`](#), on page 156 for information.

Data acquisition continues until:

- The number of readouts specified by `PicamParameter_ReadoutCount` have been acquired;
- An error occurs which immediately halts data acquisition (refer to [Section 6.4.1, `Picam_Acquire\(\)`](#), on page 162 for additional information); or
- `Picam_StopAcquisition()` is called.



NOTE:

To determine the current data acquisition status, call `Picam_WaitForAcquisitionUpdate()`.

Syntax

The syntax for `Picam_StartAcquisition()` is:

```
PICAM_API Picam_StartAcquisition(  
    PicamHandle camera);
```

Input Parameters

Input parameters for `Picam_StartAcquisition()` are:

`camera`: Handle for the camera for which data acquisition is to be initiated.

Output Parameters

There are no output parameters associated with `Picam_StartAcquisition()`.

Advanced API Usage

When used in conjunction with Advanced APIs, if `PicamParameter_ReadoutCount` = 0, the camera will run continuously until `Picam_StopAcquisition()` is called.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_CommitParameters()`;
- `Picam_Acquire()`;
- `Picam_StopAcquisition()`;
- `Picam_WaitForAcquisitionUpdate()`.

6.4.3 Picam_StopAcquisition()

Description

`Picam_StopAcquisition()` halts an in-progress data acquisition.



NOTE:

[Advanced API Usage ONLY]

If `PicamParameter_ReadoutCount = 0`, the camera will run continuously until `Picam_StopAcquisition()` has been called.

Syntax

The syntax for `Picam_StopAcquisition()` is:

```
PICAM_API Picam_StopAcquisition(  
    PicamHandle camera);
```

Input Parameters

Input parameters for `Picam_StopAcquisition()` are:

`camera`: Handle for the camera from which data acquisition is to be halted.

Output Parameters

There are no output parameters associated with `Picam_StopAcquisition()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_StartAcquisition()`.

6.4.4 Picam_IsAcquisitionRunning()

Description

`Picam_IsAcquisitionRunning()` determines if there is an active data acquisition in process.

Syntax

The syntax for `Picam_IsAcquisitionRunning()` is:

```
PICAM_API Picam_IsAcquisitionRunning(  
                                PicamHandle camera,  
                                pibln* running);
```

Input Parameters

Input parameters for `Picam_IsAcquisitionRunning()` are:

`camera`: Handle for the camera for which the data acquisition status is being determined.

Output Parameters

Output parameters for `Picam_IsAcquisitionRunning()` are:

`running`: Indicates if there is a an active data acquisition in progress.

Valid values are:

- TRUE
Indicates that there is an active data acquisition in process.
- FALSE
Indicates that there is no active data acquisition in process.

6.4.5 Picam_WaitForAcquisitionUpdate()

Description

`Picam_WaitForAcquisitionUpdate()` is used in combination with `Picam_StartAcquisition()` and indicates when:

- New data are available; or
- The camera's status has changed.

Usage

`Picam_WaitForAcquisitionUpdate()` must be continuously called until `PicamAcquisitionStatus.running` returns `FALSE`. This is true regardless of any acquisition errors that may be returned or if `Picam_StopAcquisition()` has been called.

Any errors returned during data acquisition are stored in `PicamAcquisitionStatus.errors` and acquisition is immediately halted.

However, if new data is not available within the time specified by `readout_time_out`:

- The `PicamError_TimeOutOccurred` error is returned;
- Data acquisition will continue; and
- The contents of both the data buffer `available` as well as the `status` data structure are invalid.

Syntax

The syntax for `Picam_WaitForAcquisitionUpdate()` is:

```
PICAM_API Picam_WaitForAcquisitionUpdate(  
                                PicamHandle camera,  
                                piint readout_time_out,  
                                PicamAvailableData* available,  
                                PicamAcquisitionStatus* status);
```

Input Parameters

Input parameters for `Picam_WaitForAcquisitionUpdate()` are:

- `camera`: Handle for the camera from which data is being acquired.
- `readout_time_out`: The time, in mS, to wait between each successive readout. When specifying an infinite length of time, configure this parameter to **-1**.

Output Parameters

Output parameters for `Picam_WaitForAcquisitionUpdate()` are:

- `available`: The output buffer used to store newly acquired data from the specified camera.
Data stored in this buffer is valid until the next `Picam_WaitForAcquisitionUpdate()` call.
- `status`: Pointer to the `PicamAcquisitionStatus` data structure in which acquisition status information is stored.

Advanced API Usage

When used in conjunction with Advanced APIs, data in the output buffer `available` is also invalidated when `PicamAdvanced_SetAcquisitionBuffer()` is called (in the case of the last `Picam_WaitForAcquisitionUpdate()` call.)

`Picam_WaitForAcquisitionUpdate()` is mutually exclusive with the usage of an acquisition-updated callback.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_StartAcquisition();`
- `PicamAdvanced_SetAcquisitionBuffer()`

Related Structures

For additional information, refer to the following related structure definition:

- `PicamAvailableData;`
- `PicamAcquisitionStatus.`

Chapter 7: Advanced Function APIs

This chapter provides programming information about PICam advanced function APIs, including related data definitions and structures which are included in the `picam_advanced.h` file.


Figure 7-1 illustrates a block diagram of the PICam Advanced Function structure.



REFERENCES:

Refer to [Section 2.6.1, Camera Handles](#), for information about the handles used within PICam.

Figure 7-1: PICam Structure - Advanced



4411-0155_0003

7.1 Data Type Definitions

This section provides programming information about the following PICam advanced data definitions:

- Camera Plug and Play Discovery Data Enumerations
 - [PicamDiscoveryAction](#)
- Camera Access Enumerations
 - [PicamHandleType](#)
- Camera Parameter Information Enumerations
 - [PicamDynamicsMask](#)
- Camera Data Acquisition Enumerations
 - [PicamAcquisitionState](#)
 - [PicamAcquisitionStateErrorsMask](#)

7.1.1 Camera Plug and Play Discovery Data Enumerations

This section provides programming information about camera plug and play discovery data enumerations.

7.1.1.1 *PicamDiscoveryAction*

Data Type

[PicamDiscoveryAction](#) is defined as enum.

Description

[PicamDiscoveryAction](#) is the set of discovery states available for a camera.

Enumerator Definitions

Refer to [Table 7-1](#) for enumerator definitions.

Table 7-1: [PicamDiscoveryAction](#) Enumerator Definitions

Enumerator	Description
PicamDiscoveryAction_Faulted	A camera has critically malfunctioned and is in need of service. Any acquisition in progress will be stopped and further acquisition are not possible until the camera has been serviced.
PicamDiscoveryAction_Found	A camera is now available for use.
PicamDiscoveryAction_Lost	A camera is no longer available for use.

7.1.2 Camera Access Enumerations

This section provides programming information about camera access data enumerations.

7.1.2.1 *PicamHandleType*

Data Type

PicamHandleType is defined as enum.

Description

PicamHandleType is the set of camera handle types.

Enumerator Definitions

Refer to [Table 7-2](#) for enumerator definitions.

Table 7-2: *PicamHandleType* Enumerator Definitions

Enumerator	Description
<i>PicamHandleType_CameraDevice</i>	The handle refers to a camera device.
<i>PicamHandleType_CameraModel</i>	The handle refers to a camera model.
<i>PicamHandleType_EMCalibration</i>	The handle refers to a camera opened for EM calibration.

7.1.3 Camera Parameter Information Enumerations

This section provides programming information about camera parameter information data enumerations.

7.1.3.1 *PicamDynamicsMask*

Data Type

PicamDynamicsMask is defined as enum.

Description

PicamDynamicsMask is the set of descriptors for how parameters and their various attributes may or may not change.

Enumerator Definitions

Refer to [Table 7-3](#) for enumerator definitions.

Table 7-3: *PicamDynamicsMask* Enumerator Definitions (Sheet 1 of 2)

Enumerator	Description
<i>PicamDynamicsMask_None</i>	No parameter attributes may change.
<i>PicamDynamicsMask_Value</i>	The parameter value may change.
<i>PicamDynamicsMask_ValueAccess</i>	The parameter value access may change.
<i>PicamDynamicsMask_IsRelevant</i>	The parameter relevance may change.

Table 7-3: PicamDynamicsMask Enumerator Definitions (Sheet 2 of 2)

Enumerator	Description
PicamDynamicsMask_Constraint	The parameter dependent constraints may change.

7.1.4 Camera Data Acquisition Enumerations

This section provides programming information about camera data acquisition enumerations.

7.1.4.1 PicamAcquisitionState

Data Type

[PicamAcquisitionState](#) is defined as enum.

Description

[PicamAcquisitionState](#) is the set of camera states that can be detected during an acquisition.

Enumerator Definitions

Refer to [Table 7-4](#) for enumerator definitions.

Table 7-4: PicamAcquisitionState Enumerator Definitions

Enumerator	Description
PicamAcquisitionState_ReadoutStarted	The camera has begun to readout data.
PicamAcquisitionState_ReadoutEnded	The camera has finished reading out data.

7.1.4.2 PicamAcquisitionStateErrorsMask

Data Type

[PicamAcquisitionStateErrorsMask](#) is defined as enum.

Description

[PicamAcquisitionStateErrorsMask](#) is the set of errors that can occur while detecting acquisition states.

Enumerator Definitions

Refer to [Table 7-5](#) for enumerator definitions.

Table 7-5: PicamAcquisitionStateErrorsMask Enumerator Definitions

Enumerator	Description
PicamAcquisitionStateErrorsMask_None	No error has occurred.
PicamAcquisitionStateErrorsMask_LostCount	One or more state transitions have been missed.

7.2 Data Structures

This section provides programming information about the following PICam data structures:

- Camera Information Data Structures;
 - [PicamPixelLocation](#)
 - [PicamColumnDefect](#)
 - [PicamRowDefect](#)
 - [PicamPixelDefectMap](#)
- Camera Parameter Validation Data Structures;
 - [PicamValidationResult](#)
 - [PicamValidationResults](#)
 - [PicamFailedDependentParameter](#)
 - [PicamDependentValidationResult](#)
- Camera Data Acquisition Data Structures.
 - [PicamAcquisitionBuffer](#)
 - [PicamAcquisitionStateCounters](#)

7.2.1 Camera Information Data Structures

This section provides programming information about structures used to define and describe the camera.

7.2.1.1 *PicamPixelLocation*

Description

[PicamPixelLocation](#) specifies the location of a pixel within the sensor array. A standard zero-based X-Y coordinate system is used where:

- X represents the column number;
- Y represents the row number.

Structure Definition

The structure definition for [PicamPixelLocation](#) is:

```
typedef struct PicamPixelLocation
{
    pi16s x;
    pi16s y;
} PicamPixelLocation;
```

Variable Definitions

The variables required by [PicamPixelLocation](#) are:

- x: The column coordinate.
- y: The row coordinate.

7.2.1.2 *PicamColumnDefect*

Description

`PicamColumnDefect` specifies the location and size of a single defective column on the sensor.

Structure Definition

The structure definition for `PicamPixelLocation` is:

```
typedef struct PicamColumnDefect
{
    PicamPixelLocation start;
    piint height;
} PicamColumnDefect;
```

Variable Definitions

The variables required by `PicamPixelLocation` are:

start: The top-most defective pixel.

height: The number of rows this column defect spans.

7.2.1.3 *PicamRowDefect*

Description

`PicamRowDefect` specifies the location and size of a single defective row on the sensor.

Structure Definition

The structure definition for `PicamRowDefect` is:

```
typedef struct PicamRowDefect
{
    PicamPixelLocation start;
    piint width;
} PicamRowDefect;
```

Variable Definitions

The variables required by `PicamRowDefect` are:

start: The left-most defective pixel.

height: The number of columns this row defect spans.

7.2.1.4 *PicamPixelDefectMap*

Description

[PicamPixelDefectMap](#) is an array in which all defects for a specified sensor are stored.

Structure Definition

The structure definition for [PicamPixelDefectMap](#) is:

```
typedef struct PicamPixelDefectMap
{
    const PicamColumnDefect*  column_defect_array;
                          piint column_defect_count;
    const PicamRowDefect*    row_defect_array;
                          piint row_defect_count;
    const PicamPixelLocation* point_defect_array;
                          piint point_defect_count;
} PicamPixelDefectMap;
```

Variable Definitions

The variables required by [PicamPixelDefectMap](#) are:

<code>column_defect_array</code> :	A set of all column defects. This is null where there are no column defects.
<code>column_defect_count</code> :	The number of items in column_defect_array . This is 0 when there are no defective columns.
<code>row_defect_array</code> :	A set of all row defects. This is null when there are no row defects.
<code>row_defect_count</code> :	The number of items in row_defect_array . This is 0 when there are no defective rows.
<code>point_defect_array</code> :	A set of all single-point defects. This is null when there are no single-point defects.
<code>point_defect_count</code> :	The number of items in point_defect_array . This is 0 when there are no single-point defects.

7.2.2 Camera Parameter Validation Data Structures

This section provides programming information about parameter validation structures.

7.2.2.1 *PicamValidationResult*

Description

PicamValidationResult provides information about the validation status for a single parameter.

Structure Definition

The structure definition for *PicamValidationResult* is:

```
typedef struct PicamValidationResult
{
    pibln    is_valid;
    const PicamParameter*  failed_parameter;
    const PicamConstraintScope* failed_error_constraint_scope;
    const PicamConstraintScope* failed_warning_constraint_scope;
    const PicamParameter*  error_constraining_parameter_array;
    piint    error_constraining_parameter_count;
    const PicamParameter*  warning_constraining_parameter_array;
    piint    warning_constraining_parameter_count;
} PicamValidationResult;
```

Variable Definitions

The variables required by *PicamValidationResult* are:

is_valid: Indicates the validation status for a single parameter.

Valid values are:

- TRUE
Indicates the parameter validation has succeeded.
- FALSE
Indicates the parameter validation has failed.

failed_parameter: The parameter that has failed validation.

This is null when validation has succeeded.

failed_error_constraint: The scope of the error constraint that has failed.

_scope: This is null when:

- Validation has succeeded, or
- Only a warning constraint has failed validation.

failed_warning_constraint_scope: The scope of the warning constraint that has failed.

This is null when:

- Validation has succeeded, or
- Only an error constraint has failed validation.

error_constraining_parameter_array: An array of parameters involved in constraining the failed parameter when a dependent error failed.

This is null otherwise.

continued on next page

continued from previous page

<code>error_constraining_parameter_count:</code>	The number of items in the array of parameters involved in constraining the failed parameter if a dependent error failed. This is 0 otherwise.
<code>warning_constraining_parameter_array:</code>	An array of parameters involved in constraining the failed parameter if a dependent warning failed. This is null otherwise.
<code>warning_constraining_parameter_count:</code>	The number of items in the array of parameters involved in constraining the failed parameter if a dependent warning failed. This is 0 otherwise.

7.2.2.2 *PicamValidationResults*

Description

`PicamValidationResults` provides information about the validation status for multiple parameters.

Structure Definition

The structure definition for `PicamValidationResults` is:

```
typedef struct PicamValidationResults
{
    pibln is_valid;
    const PicamValidationResult* validation_result_array;
    piint validation_result_count;
} PicamValidationResults;
```

Variable Definitions

The variables required by `PicamValidationResults` are:

<code>is_valid:</code>	Indicates the validation status for multiple tested parameters. Valid values are: <ul style="list-style-type: none"> • <code>TRUE</code> Indicates all parameter validations have succeeded. • <code>FALSE</code> Indicates one or more parameter validations has failed.
<code>validation_result_array:</code>	An array containing a result for each parameter that failed validation; null if validation succeeded.
<code>validation_result_count:</code>	The number of failed parameter results; 0 if validation succeeded.

7.2.2.3 *PicamFailedDependentParameter*

Description

[PicamFailedDependentParameter](#) provides information about a parameter that has failed validation and is itself constrained by a second parameter.

Structure Definition

The structure definition for [PicamFailedDependentParameter](#) is:

```
typedef struct PicamFailedDependentParameter
{
    PicamParameter    failed_parameter;
    const PicamConstraintScope* failed_error_constraint_scope;
    const PicamConstraintScope* failed_warning_constraint_scope;
} PicamFailedDependentParameter;
```

Variable Definitions

The variables required by [PicamFailedDependentParameter](#) are:

`failed_parameter`: The parameter whose validation failed and is constrained by another.

`failed_error_constraint_scope`: The scope of the error constraint that failed.
`constraint_scope`: This is null when only a warning constraint failed.

`failed_warning_constraint_scope`: The scope of the warning constraint that failed.
`constraint_scope`: This is null when only an error constraint failed.

7.2.2.4 *PicamDependentValidationResult*

Description

[PicamDependentValidationResult](#) provides information about the failed validation of a parameter that is constrained by a second parameter.

Structure Definition

The structure definition for [PicamDependentValidationResult](#) is:

```
typedef struct PicamDependentValidationResult
{
    pibln    is_valid;
    PicamParameter    constraining_parameter;
    const PicamFailedDependentParameter*    failed_dependent_parameter_array;
    piint    failed_dependent_parameter_count;
} PicamDependentValidationResult;
```

Variable Definitions

The variables required by [PicamDependentValidationResult](#) are:

is_valid: Indicates the validation status for a parameter that is constrained by a second parameter.

Valid values are:

- **TRUE**
Indicates the parameter validation has succeeded.
- **FALSE**
Indicates the parameter validation has failed.

constraining_parameter: The parameter whose value impacts the constraints of another.

failed_dependent_parameter_array: An array containing all parameters whose constraints are dependent on [constraining_parameter](#) and that have failed validation.

This is null when the validation has succeeded.

failed_dependent_parameter_count: The number of items in an array containing all parameters whose constraints are dependent on [constraining_parameter](#) and that have failed validation.
This is 0 when the validation has succeeded.

7.2.3 Camera Data Acquisition Data Structures

This section provides programming information for camera data acquisition structures.

7.2.3.1 *PicamAcquisitionBuffer*

Description

PicamAcquisitionBuffer is a user-allocated buffer into which acquired data is stored.

Structure Definition

The structure definition for *PicamAcquisitionBuffer* is:

```
typedef struct PicamAcquisitionBuffer
{
    void*   memory;
    pi64s   memory_size;
} PicamAcquisitionBuffer;
```

Variable Definitions

The variables required by *PicamAcquisitionBuffer* are:

memory: Pointer to the top of the user-allocated memory location.

memory_size: Number of bytes allocated for use by the user-allocated memory.

7.2.3.2 *PicamAcquisitionStateCounters*

Description

PicamAcquisitionStateCounters counts all acquisition state transitions registered for detection while acquiring.

Structure Definition

The structure definition for *PicamAcquisitionStateCounters* is:

```
typedef struct PicamAcquisitionStateCounters
{
    pi64s   readout_started_count;
    pi64s   readout_ended_count;
} PicamAcquisitionStateCounters;
```

Variable Definitions

The variables required by *PicamAcquisitionStateCounters* are:

readout_started_count: The number of occurrences where the camera has begun to readout.

readout_ended_count: The number of occurrences where the camera has finished readout.

7.3 Callback Functions

This section provides programming information about the following callbacks used within PICam:

- Camera Discovery Callbacks
 - `PicamDiscoveryCallback()`
- Camera Parameter Value Callbacks
 - `PicamIntegerValueChangedCallback()`
 - `PicamLargeIntegerValueChangedCallback()`
 - `PicamFloatingPointValueChangedCallback()`
 - `PicamRoisValueChangedCallback()`
 - `PicamPulseValueChangedCallback()`
 - `PicamModulationsValueChangedCallback()`
 - `PicamWhenStatusParameterValueCallback()`
 - `PicamIsRelevantChangedCallback()`
 - `PicamValueAccessChangedCallback()`
- Camera Parameter Constraints Callbacks
 - `PicamDependentCollectionConstraintChangedCallback()`
 - `PicamDependentRangeConstraintChangedCallback()`
 - `PicamDependentRoisConstraintChangedCallback()`
 - `PicamDependentPulseConstraintChangedCallback()`
 - `PicamDependentModulationsConstraintChangedCallback()`
- Camera Data Acquisition Callbacks
 - `PicamAcquisitionUpdatedCallback()`
 - `PicamAcquisitionStateUpdatedCallback()`

7.3.1 Camera Discovery Callbacks

This section provides programming information about camera discovery callbacks.

7.3.1.1 *PicamDiscoveryCallback()*

Description

`PicamDiscoveryCallback()` is the callback function for camera discovery.

Syntax

The syntax for `PicamDiscoveryCallback()` is:

```
typedef PicamError (PIL_CALL* PicamDiscoveryCallback)
(
    const PicamCameraID* id,
    PicamHandle device,
    PicamDiscoveryAction action );
```

Input Parameters

The input parameters for `PicamDiscoveryCallback()` are:

- `id`: Pointer to the camera that has been discovered.
- `device`: The handle for an open camera device if `id` is open within this process.
This is null otherwise.
- `action`: The type of discovery.

7.3.2 Camera Parameter Value Callbacks

This section provides programming information about camera parameter value callbacks.

7.3.2.1 *PicamIntegerValueChangedCallback()*

Description

[PicamIntegerValueChangedCallback\(\)](#) is the change notification callback function called when a parameter's integer value has been changed.

Syntax

The syntax for [PicamIntegerValueChangedCallback\(\)](#) is:

```
typedef PicamError (PIL_CALL* PicamIntegerValueChangedCallback)
(
    PicamHandle camera,
    PicamParameter parameter,
    piint value );
```

Input Parameters

Input parameters for [PicamIntegerValueChangedCallback\(\)](#) are:

camera: Handle for the camera for which a parameter's integer value has been changed.

parameter: The parameter which has had its integer value changed.

value: The new integer value.

7.3.2.2 *PicamLargeIntegerValueChangedCallback()*

Description

[PicamLargeIntegerValueChangedCallback\(\)](#) is the change notification callback function called when a parameter's large integer value has been changed.

Syntax

The syntax for [PicamLargeIntegerValueChangedCallback\(\)](#) is:

```
typedef PicamError (PIL_CALL* PicamLargeIntegerValueChangedCallback)
(
    PicamHandle camera,
    PicamParameter parameter,
    pi64s value );
```

Input Parameters

Input parameters for [PicamLargeIntegerValueChangedCallback\(\)](#) are:

camera: Handle for the camera for which a parameter's large integer value has been changed.

parameter: The parameter which has had its large integer value changed.

value: The new large integer value.

7.3.2.3 *PicamFloatingPointValueChangedCallback()*

Description

`PicamFloatingPointValueChangedCallback()` is the change notification callback function called when a parameter's floating point value has been changed.

Syntax

The syntax for `PicamFloatingPointValueChangedCallback()` is:

```
typedef PicamError (PIL_CALL* PicamFloatingPointValueChangedCallback)
(
    PicamHandle camera,
    PicamParameter parameter,
    piflt value );
```

Input Parameters

Input parameters for `PicamFloatingPointValueChangedCallback()` are:

camera: Handle for the camera for which a parameter's floating point value has been changed.

parameter: The parameter which has had its floating point value changed.

value: The new floating point value.

7.3.2.4 *PicamRoIsValueValueChangedCallback()*

Description

`PicamRoIsValueValueChangedCallback()` is the change notification callback function called when a parameter's RoIs value has been changed.

Syntax

The syntax for `PicamRoIsValueValueChangedCallback()` is:

```
typedef PicamError (PIL_CALL* PicamRoIsValueValueChangedCallback)
(
    PicamHandle camera,
    PicamParameter parameter,
    const PicamRoIs* value );
```

Input Parameters

Input parameters for `PicamRoIsValueValueChangedCallback()` are:

camera: Handle for the camera for which a parameter's RoIs value has been changed.

parameter: The parameter which has had its RoIs value changed.

value: Pointer the array location in which the new RoIs value is stored.

7.3.2.5 *PicamPulseValueChangedCallback()*

Description

`PicamPulseValueChangedCallback()` is the change notification callback function called when a parameter's gate pulse value has been changed.

Syntax

The syntax for `PicamPulseValueChangedCallback()` is:

```
typedef PicamError (PIL_CALL* PicamPulseValueChangedCallback)
(
    PicamHandle camera,
    PicamParameter parameter,
    const PicamPulse* value );
```

Input Parameters

Input parameters for `PicamPulseValueChangedCallback()` are:

camera: Handle for the camera for which a parameter's gate pulse value has been changed.

parameter: The parameter which has had its gate pulse value changed.

value: Pointer the array in which the new gate pulse value is stored.

7.3.2.6 *PicamModulationsValueChangedCallback()*

Description

`PicamModulationsValueChangedCallback()` is the change notification callback function called when a parameter's intensifier modulation sequence value has been changed.

Syntax

The syntax for `PicamModulationsValueChangedCallback()` is:

```
typedef PicamError (PIL_CALL* PicamModulationsValueChangedCallback)
(
    PicamHandle camera,
    PicamParameter parameter,
    const PicamModulations* value );
```

Input Parameters

Input parameters for `PicamModulationsValueChangedCallback()` are:

camera: Handle for the camera for which a parameter's intensifier modulation sequence value has been changed.

parameter: The parameter which has had its intensifier modulation sequence value changed.

value: Pointer the array in which the new intensifier modulation sequence value is stored.

7.3.2.7 *PicamWhenStatusParameterValueCallback()*

Description

`PicamWhenStatusParameterValueCallback()` is the notification callback function called when a waitable status value has been met or an error has occurred.

Syntax

The syntax for `PicamWhenStatusParameterValueCallback()` is:

```
typedef PicamError (PIL_CALL* PicamWhenStatusParameterValueCallback)
(
    PicamHandle  device,
    PicamParameter  parameter,
    piint  value,
    PicamError  error);
```

Input Parameters

Input parameters for `PicamWhenStatusParameterValueCallback()` are:

device: Handle for the camera device for which a parameter's status value has been met.

parameter: The parameter whose status value has been met.

value: The status value that has been met.

error: Any error that occurred to prevent the status value from being met.

7.3.2.8 *PicamIsRelevantChangedCallback()*

Description

`PicamIsRelevantChangedCallback()` is the change notification callback function called when a parameter's relevance has been changed.

Syntax

The syntax for `PicamIsRelevantChangedCallback()` is:

```
typedef PicamError (PIL_CALL* PicamIsRelevantChangedCallback)
(
    PicamHandle  camera,
    PicamParameter  parameter,
    piibln  relevant );
```

Input Parameters

Input parameters for `PicamIsRelevantChangedCallback()` are:

camera: Handle for the camera for which a parameter's relevance has been changed.

parameter: The parameter which has had its relevance changed.

relevant: The new relevance.

7.3.2.9 *PicamValueAccessChangedCallback()*

Description

`PicamValueAccessChangedCallback()` is the change notification callback function called when a parameter's value access has been changed.

Syntax

The syntax for `PicamValueAccessChangedCallback()` is:

```
typedef PicamError (PIL_CALL* PicamValueAccessChangedCallback)
(
    PicamHandle camera,
    PicamParameter parameter,
    PicamValueAccess access );
```

Input Parameters

Input parameters for `PicamValueAccessChangedCallback()` are:

- camera: Handle for the camera for which a parameter's value access has been changed.
- parameter: The parameter which has had its value access changed.
- access: The new value access.

7.3.3 Camera Parameter Constraints Callbacks

This section provides programming information about camera parameter constraints callbacks.

7.3.3.1 *PicamDependentCollectionConstraintChangedCallback()*

Description

`PicamDependentCollectionConstraintChangedCallback()` is the change notification callback function called when a parameter's dependent collection constraints have been changed.

Syntax

The syntax for `PicamDependentCollectionConstraintChangedCallback()` is:

```
typedef PicamError (PIL_CALL*PicamDependentCollectionConstraint  
    ChangedCallback)  
(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamCollectionConstraint* constraint );
```

Input Parameters

Input parameters for `PicamDependentCollectionConstraintChangedCallback()` are:

- `camera`: Handle for the camera for which a parameter's dependent collection constraints have been changed.
- `parameter`: The parameter which has had its dependent collection constraints changed.
- `constraint`: Pointer to the array in which the new dependent collection constraints are stored.

7.3.3.2 *PicamDependentRangeConstraintChangedCallback()*

Description

`PicamDependentRangeConstraintChangedCallback()` is the change notification callback function called when a parameter's dependent range constraints have been changed.

Syntax

The syntax for `PicamDependentRangeConstraintChangedCallback()` is:

```
typedef PicamError (PIL_CALL*PicamDependentRangeConstraint
    ChangedCallback)
(
    PicamHandle camera,
    PicamParameter parameter,
    const PicamRangeConstraint* constraint );
```

Input Parameters

Input parameters for `PicamDependentRangeConstraintChangedCallback()` are:

- camera: Handle for the camera for which a parameter's dependent range constraints have been changed.
- parameter: The parameter which has had its dependent range constraints changed.
- constraint: Pointer to the array in which the new dependent range constraints are stored.

7.3.3.3 *PicamDependentRoisConstraintChangedCallback()*

Description

`PicamDependentRoisConstraintChangedCallback()` is the change notification callback function called when a parameter's dependent Rois constraints have been changed.

Syntax

The syntax for `PicamDependentRoisConstraintChangedCallback()` is:

```
typedef PicamError (PIL_CALL* PicamDependentRoisConstraintChanged
    Callback)
(
    PicamHandle camera,
    PicamParameter parameter,
    const PicamRoisConstraint* constraint );
```

Input Parameters

Input parameters for `PicamDependentRoisConstraintChangedCallback()` are:

- camera: Handle for the camera for which a parameter's dependent Rois constraints have been changed.
- parameter: The parameter which has had its dependent Rois constraints changed.
- constraint: Pointer to the array in which the new dependent Rois constraints are stored.

7.3.3.4 *PicamDependentPulseConstraintChangedCallback()*

Description

`PicamDependentPulseConstraintChangedCallback()` is the change notification callback function called when a parameter's dependent gate pulse constraints have been changed.

Syntax

The syntax for `PicamDependentPulseConstraintChangedCallback()` is:

```
typedef PicamError (PIL_CALL* PicamDependentPulseConstraintChanged
                  Callback)
(
    PicamHandle camera,
    PicamParameter parameter,
    const PicamPulseConstraint* constraint );
```

Input Parameters

Input parameters for `PicamDependentPulseConstraintChangedCallback()` are:

`camera`: Handle for the camera for which a parameter's dependent gate pulse constraints have been changed.

`parameter`: The parameter which has had its dependent gate pulse constraints changed.

`constraint`: Pointer to the array in which the new dependent gate pulse constraints are stored.

7.3.3.5 *PicamDependentModulationsConstraintChangedCallback()*

Description

`PicamDependentModulationsConstraintChangedCallback()` is the change notification callback function called when a parameter's dependent intensifier modulations sequence constraints have been changed.

Syntax

The syntax for `PicamDependentModulationsConstraintChangedCallback()` is:

```
typedef PicamError (PIL_CALL* PicamDependentModulationsConstraint
                  ChangedCallback)
(
    PicamHandle camera,
    PicamParameter parameter,
    const PicamModulationsConstraint* constraint );
```

Input Parameters

Input parameters for `PicamDependentModulationsConstraintChangedCallback()` are:

`camera`: Handle for the camera for which a parameter's dependent intensifier modulations sequence constraints have been changed.

`parameter`: The parameter which has had its dependent intensifier modulations sequence constraints changed.

`constraint`: Pointer to the array in which the new dependent intensifier modulations sequence constraints are stored.

7.3.4 Camera Data Acquisition Callbacks

This section provides programming information about camera data acquisition callbacks.

7.3.4.1 *PicamAcquisitionUpdatedCallback()*

Description

PicamAcquisitionUpdatedCallback() is the change notification callback function called when a camera's data acquisition status has changed.

Syntax

The syntax for *PicamAcquisitionUpdatedCallback()* is:

```
typedef PicamError (PIL_CALL* PicamAcquisitionUpdatedCallback)
(
    PicamHandle device,
    const PicamAvailableData* available,
    const PicamAcquisitionStatus* status );
```

Input Parameters

Input parameters for *PicamAcquisitionUpdatedCallback()* are:

- device: Handle for the camera which is acquiring data.
- available: Pointer to the array in which newly acquired data are stored.
If no data are available, this is null.
- status: Pointer to the data acquisition status.

7.3.4.2 *PicamAcquisitionStateUpdatedCallback()*

Description

PicamAcquisitionStateUpdatedCallback() is the notification callback function called when a camera has transitioned into the acquisition state requested for detection.

Syntax

The syntax for *PicamAcquisitionStateUpdatedCallback()* is:

```
typedef PicamError (PIL_CALL* PicamAcquisitionStateUpdatedCallback)
(
    PicamHandle device,
    PicamAcquisitionState current,
    const PicamAcquisitionStateCounters* counters,
    PicamAcquisitionStateErrorsMask errors );
```

Input Parameters

Input parameters for *PicamAcquisitionStateUpdatedCallback()* are:

- device: Handle for the device which transitioned into the acquisition state.
- current: Acquisition state whose transition was detected.
- counters: Pointer to the counted transitions at the time of detection.
- errors: Indicates if any errors have occurred.

7.4 Programmers' Reference for Advanced APIs

This section provides a detailed programmers' reference guide for the following advanced APIs:

- Camera Discovery APIs
 - `PicamAdvanced_RegisterForDiscovery()`
 - `PicamAdvanced_UnregisterForDiscovery()`
 - `PicamAdvanced_DiscoverCameras()`
 - `PicamAdvanced_StopDiscoveringCameras()`
 - `PicamAdvanced_IsDiscoveringCameras()`
- Camera Access APIs
 - `PicamAdvanced_OpenCameraDevice()`
 - `PicamAdvanced_CloseCameraDevice()`
 - `PicamAdvanced_GetOpenCameraDevices()`
 - `PicamAdvanced_GetCameraModel()`
 - `PicamAdvanced_GetCameraDevice()`
 - `PicamAdvanced_GetHandleType()`
- Camera Information APIs
 - `PicamAdvanced_GetUserState()`
 - `PicamAdvanced_SetUserState()`
 - `PicamAdvanced_DestroyPixelDefectMaps()`
 - `PicamAdvanced_GetPixelDefectMap()`
- Camera Parameter Value APIs
 - `PicamAdvanced_RegisterForIntegerValueChanged()`
 - `PicamAdvanced_UnregisterForIntegerValueChanged()`
 - `PicamAdvanced_RegisterForExtrinsicIntegerValueChanged()`
 - `PicamAdvanced_UnregisterForExtrinsicIntegerValueChanged()`
 - `PicamAdvanced_RegisterForLargeIntegerValueChanged()`
 - `PicamAdvanced_UnregisterForLargeIntegerValueChanged()`
 - `PicamAdvanced_RegisterForFloatingPointValueChanged()`
 - `PicamAdvanced_UnregisterForFloatingPointValueChanged()`
 - `PicamAdvanced_RegisterForExtrinsicFloatingPointValueChanged()`
 - `PicamAdvanced_UnregisterForExtrinsicFloatingPointValueChanged()`
 - `PicamAdvanced_RegisterForRoisValueChanged()`
 - `PicamAdvanced_UnregisterForRoisValueChanged()`
 - `PicamAdvanced_RegisterForPulseValueChanged()`
 - `PicamAdvanced_UnregisterForPulseValueChanged()`
 - `PicamAdvanced_RegisterForModulationsValueChanged()`
 - `PicamAdvanced_UnregisterForModulationsValueChanged()`
 - `PicamAdvanced_NotifyWhenStatusParameterValue()`
 - `PicamAdvanced_CancelNotifyWhenStatusParameterValue()`
- Camera Parameter Information APIs
 - `PicamAdvanced_RegisterForIsRelevantChanged()`
 - `PicamAdvanced_UnregisterForIsRelevantChanged()`
 - `PicamAdvanced_RegisterForValueAccessChanged()`
 - `PicamAdvanced_UnregisterForValueAccessChanged()`
 - `PicamAdvanced_GetParameterDynamics()`
 - `PicamAdvanced_GetParameterExtrinsicDynamics()`
- Camera Parameter Constraints APIs
 - `PicamAdvanced_GetParameterCollectionConstraints()`
 - `PicamAdvanced_RegisterForDependentCollectionConstraintChanged()`
 - `PicamAdvanced_UnregisterForDependentCollectionConstraintChanged()`
 - `PicamAdvanced_GetParameterRangeConstraints()`
 - `PicamAdvanced_RegisterForDependentRangeConstraintChanged()`
 - `PicamAdvanced_UnregisterForDependentRangeConstraintChanged()`
 - `PicamAdvanced_GetParameterRoisConstraints()`
 - `PicamAdvanced_RegisterForDependentRoisConstraintChanged()`

- `PicamAdvanced_UnregisterForDependentRoisConstraintChanged()`
- `PicamAdvanced_GetParameterPulseConstraints()`
- `PicamAdvanced_RegisterForDependentPulseConstraintChanged()`
- `PicamAdvanced_UnregisterForDependentPulseConstraintChanged()`
- `PicamAdvanced_GetParameterModulationsConstraints()`
- `PicamAdvanced_RegisterForDependentModulationsConstraintChanged()`
- `PicamAdvanced_UnregisterForDependentModulationsConstraintChanged()`
- **Camera Commitment APIs**
 - `Picam_DestroyValidationResult()`
 - `Picam_DestroyValidationResults()`
 - `PicamAdvanced_ValidateParameter()`
 - `PicamAdvanced_ValidateParameters()`
 - `Picam_DestroyDependentValidationResult()`
 - `PicamAdvanced_ValidateDependentParameter()`
 - `PicamAdvanced_CommitParametersToCameraDevice()`
 - `PicamAdvanced_RefreshParameterFromCameraDevice()`
 - `PicamAdvanced_RefreshParametersFromCameraDevice()`
- **Acquisition Setup APIs**
 - `PicamAdvanced_GetAcquisitionBuffer()`
 - `PicamAdvanced_SetAcquisitionBuffer()`
- **Acquisition Notification APIs**
 - `PicamAdvanced_RegisterForAcquisitionUpdated()`
 - `PicamAdvanced_UnregisterForAcquisitionUpdated()`
- **Acquisition State Notification APIs**
 - `PicamAdvanced_CanRegisterForAcquisitionStateUpdated()`
 - `PicamAdvanced_RegisterForAcquisitionStateUpdated()`
 - `PicamAdvanced_UnregisterForAcquisitionStateUpdated()`
- **Acquisition Control APIs**
 - `PicamAdvanced_HasAcquisitionBufferOverrun()`
 - `PicamAdvanced_ClearReadoutCountOnline()`

7.4.1 Camera Discovery APIs

This section provides programming information for advanced camera discovery APIs.

7.4.1.1 *PicamAdvanced_RegisterForDiscovery()*

Description

`PicamAdvanced_RegisterForDiscovery()` registers a function to call when camera discovery is made.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Callback functions are called when any camera state that affects availability changes, such as when:

- A camera is powered on and/or connected to the host computer;
- A connected camera is powered off or disconnected from the host computer;
- A camera is opened in another process;
- A camera is closed in another process.

Callback functions are also called when a camera has suffered a critical malfunction.

Callbacks are called asynchronously from another thread, but are serialized on that thread. This means that additional notifications do not occur simultaneously, but occur after each callback returns.

A camera may be unavailable for multiple reasons. Therefore, although callbacks may repeatedly indicate a camera is lost each time one of the above states change, but the camera is still not available.

Call `PicamAdvanced_UnregisterForDiscovery()` to unregister each callback once it is no longer required.

Syntax

The syntax for `PicamAdvanced_RegisterForDiscovery()` is:

```
PICAM_API PicamAdvanced_RegisterForDiscovery(  
    PicamDiscoveryCallback discover )
```

Input Parameters

Input parameters for `PicamAdvanced_RegisterForDiscovery()` are:

`discover`: The name assigned to the discovery callback function being registered.

Output Parameters

There are no output parameters associated with `PicamAdvanced_RegisterForDiscovery()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_UnregisterForDiscovery()`

7.4.1.2 *PicamAdvanced_UnregisterForDiscovery()*

Description

`PicamAdvanced_UnregisterForDiscovery()` removes the function from the discovery process such that it is no longer called when a camera discovery is made.

Syntax

The syntax for `PicamAdvanced_UnregisterForDiscovery()` is:

```
PICAM_API PicamAdvanced_UnregisterForDiscovery(  
    PicamDiscoveryCallback discover );
```

Input Parameters

Input parameters for `PicamAdvanced_UnregisterForDiscovery()` are:

`discover`: The name assigned to the discovery callback function being unregistered.

Output Parameters

There are no output parameters associated with `PicamAdvanced_UnregisterForDiscovery()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForDiscovery()`

7.4.1.3 *PicamAdvanced_DiscoverCameras()*

Description

`PicamAdvanced_DiscoverCameras()` asynchronously initiates the camera discovery process.

To halt the discovery process, call `PicamAdvanced_StopDiscoveringCameras()`.

Syntax

The syntax for `PicamAdvanced_DiscoverCameras()` is:

```
PICAM_API PicamAdvanced_DiscoverCameras ( void );
```

Input Parameters

There are no input parameters associated with `PicamAdvanced_DiscoverCameras()`.

Output Parameters

There are no output parameters associated with `PicamAdvanced_DiscoverCameras()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_StopDiscoveringCameras()`

7.4.1.4 *PicamAdvanced_StopDiscoveringCameras()*

Description

`PicamAdvanced_StopDiscoveringCameras()` stops the camera discovery process.

Syntax

The syntax for `PicamAdvanced_StopDiscoveringCameras()` is:

```
PICAM_API PicamAdvanced_StopDiscoveringCameras ( void );
```

Input Parameters

There are no input parameters associated with `PicamAdvanced_StopDiscoveringCameras()`.

Output Parameters

There are no output parameters associated with `PicamAdvanced_StopDiscoveringCameras()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_DiscoverCameras()`

7.4.1.5 *PicamAdvanced_IsDiscoveringCameras()*

Description

`PicamAdvanced_IsDiscoveringCameras()` determines if camera discovery is enabled.

Syntax

The syntax for `PicamAdvanced_IsDiscoveringCameras()` is:

```
PICAM_API PicamAdvanced_IsDiscoveringCameras(  
    pibln*  discovering );
```

Input Parameters

There are no input parameters associated with `PicamAdvanced_IsDiscoveringCameras()`.

Output Parameters

Output parameters for `PicamAdvanced_IsDiscoveringCameras()` are:

`discovering`: Indicates if camera discovery is currently enabled.

Valid values are:

- TRUE
Camera discovery is enabled.
- FALSE
Camera discovery is disabled.

7.4.2 Camera Access APIs

This section provides programming information for advanced camera access APIs.

7.4.2.1 *PicamAdvanced_OpenCameraDevice()*

Description

`PicamAdvanced_OpenCameraDevice()` opens the specified camera and returns a handle to the device.

When done, all resources that have been assigned for use by the camera/device must be released by calling:

- `Picam_CloseCamera();` or
- `PicamAdvanced_CloseCameraDevice().`

Syntax

The syntax for `PicamAdvanced_OpenCameraDevice()` is:

```
PICAM_API PicamAdvanced_OpenCameraDevice(  
    const PicamCameraID* id,  
    PicamHandle* device );
```

Input Parameters

Input parameters for `PicamAdvanced_OpenCameraDevice()` are:

`id`: Pointer to the camera id for the camera device to be opened.

Output Parameters

Output parameters for `PicamAdvanced_OpenCameraDevice()` are:

`device`: Pointer to the handle assigned to the camera device that has been opened.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_CloseCamera();`
- `PicamAdvanced_CloseCameraDevice().`

7.4.2.2 *PicamAdvanced_CloseCameraDevice()*

Description

`PicamAdvanced_CloseCameraDevice()` releases all resources associated with the specified device.

Syntax

The syntax for `PicamAdvanced_CloseCameraDevice()` is:

```
PICAM_API PicamAdvanced_CloseCameraDevice(  
                                         PicamHandle device )
```

Input Parameters

Input parameters for `PicamAdvanced_CloseCameraDevice()` are:

`device`: Handle for the camera for which all resources are to be released.

Output Parameters

There are no output parameters associated with `PicamAdvanced_CloseCameraDevice()`.

7.4.2.3 *PicamAdvanced_GetOpenCameraDevices()*

Description

PicamAdvanced_GetOpenCameraDevices() returns an allocated array of open camera device handles.

in *device_array* whose number of items is in *device_count*.

Returns null and 0 (respectively) if no cameras are opened in this process

Syntax

The syntax for *PicamAdvanced_GetOpenCameraDevices()* is:

```
PICAM_API PicamAdvanced_GetOpenCameraDevices(  
    const PicamHandle** device_array,  
    piint* device_count );
```

Input Parameters

There are no input parameters associated with *PicamAdvanced_GetOpenCameraDevices()*.

Output Parameters

Output parameters for *PicamAdvanced_GetOpenCameraDevices()* are:

device_array: Pointer to the array in which the list of handles for open camera devices is stored.

This is null when there are no open camera devices.

NOTE: This memory is allocated by PICam and must be released by calling *Picam_DestroyHandles()*.

device_count: Pointer to the memory location in which the number of open camera devices is stored.

This is 0 when there are no open camera devices.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_DestroyHandles()*

7.4.2.4 *PicamAdvanced_GetCameraModel()*

Description

PicamAdvanced_GetCameraModel() returns the handle for a specified camera model.

Syntax

The syntax for *PicamAdvanced_GetCameraModel()* is:

```
PICAM_API PicamAdvanced_GetCameraModel(  
    PicamHandle camera,  
    PicamHandle* model );
```

Input Parameters

Input parameters for *PicamAdvanced_GetCameraModel()* are:

camera: Specifies the camera model or camera device for which the handle is to be returned.

Output Parameters

Output parameters for *PicamAdvanced_GetCameraModel()* are:

model: Pointer to the memory location in which the handle for the camera model is stored.

7.4.2.5 *PicamAdvanced_GetCameraDevice()*

Description

PicamAdvanced_GetCameraDevice() returns the handle for a specified camera device.

Syntax

The syntax for *PicamAdvanced_GetCameraDevice()* is:

```
PICAM_API PicamAdvanced_GetCameraDevice(  
    PicamHandle camera,  
    PicamHandle* device );
```

Input Parameters

Input parameters for *PicamAdvanced_GetCameraDevice()* are:

camera: Specifies the camera device or camera model for which the handle is to be returned.

Output Parameters

Output parameters for *PicamAdvanced_GetCameraDevice()* are:

device: Pointer to the memory location in which the handle for the camera device is stored.

7.4.2.6 *PicamAdvanced_GetHandleType()*

Description

`PicamAdvanced_GetHandleType()` returns the type of handle for a specified handle.

Syntax

The syntax for `PicamAdvanced_GetHandleType()` is:

```
PICAM_API PicamAdvanced_GetHandleType(  
    PicamHandle  handle,  
    PicamHandleType*  type );
```

Input Parameters

Input parameters for `PicamAdvanced_GetHandleType()` are:

`handle`: Handle for which the handle type is to be determined.

Output Parameters

Output parameters for `PicamAdvanced_GetHandleType()` are:

`type`: The handle type for the specified handle.

7.4.3 Camera Information APIs

This section provides programming information about advanced camera information APIs.

7.4.3.1 *PicamAdvanced_GetUserState()*

Description

PicamAdvanced_GetUserState() returns user-state information for a specified camera.



NOTE:

This API is thread safe.

Syntax

The syntax for *PicamAdvanced_GetUserState()* is:

```
PICAM_API PicamAdvanced_GetUserState(  
                                PicamHandle camera,  
                                void** user_state );
```

Input Parameters

Input parameters for *PicamAdvanced_GetUserState()* are:

camera: Handle for the camera for which user state information is to be returned.

Valid values are:

- device handle;
- model handle.

NOTE: device and model share the same user state.

Output Parameters

Output parameters for *PicamAdvanced_GetUserState()* are:

user_state: Pointer to the memory location where user-state information is stored.

7.4.3.2 *PicamAdvanced_SetUserState()*

Description

`PicamAdvanced_SetUserState()` sets user-state information for a specified camera.



NOTE:

This API is thread safe.

Syntax

The syntax for `PicamAdvanced_SetUserState()` is:

```
PICAM_API PicamAdvanced_SetUserState(  
                                PicamHandle camera,  
                                void* user_state );
```

Input Parameters

Input parameters for `PicamAdvanced_SetUserState()` are:

camera: Handle for the camera for which user state information is to be configured.

Valid values are:

- device handle;
- model handle.

NOTE: device and model share the same user state.

Output Parameters

Output parameters for `PicamAdvanced_SetUserState()` are:

user_state: Pointer to the memory location where user-state information is stored.

7.4.3.3 *PicamAdvanced_DestroyPixelDefectMaps()*

Description

`PicamAdvanced_DestroyPixelDefectMaps()` releases memory that has been allocated by PICam for use by `defect_map_array`.

If `defect_map_array` is null, calling `PicamAdvanced_DestroyPixelDefectMaps()` has no effect.



NOTE:

`defect_map_array` may be a single `PicamPixelDefectMap` allocated by PICam.

Syntax

The syntax for `PicamAdvanced_DestroyPixelDefectMaps()` is:

```
PICAM_API PicamAdvanced_DestroyPixelDefectMaps(  
    const PicamPixelDefectMap* pixel_defect_map_array );
```

Input Parameters

Input parameters for `PicamAdvanced_DestroyPixelDefectMaps()` are:

`pixel_defect_map_array`: Pointer to the array that is to be released.

Output Parameters

There are no output parameters associated with `PicamAdvanced_DestroyPixelDefectMaps()`.

Related Structures

For additional information, refer to the following related structures:

- `PicamPixelDefectMap`

7.4.3.4 *PicamAdvanced_GetPixelDefectMap()*

Description

[PicamAdvanced_GetPixelDefectMap\(\)](#) returns an allocated array/map in which defective pixels information for a specified camera is stored.

Syntax

The syntax for [PicamAdvanced_GetPixelDefectMap\(\)](#) is:

```
PICAM_API PicamAdvanced_GetPixelDefectMap(  
                                PicamHandle camera,  
                                const PicamPixelDefectMap** pixel_defect_map );
```

Input Parameters

Input parameters for [PicamAdvanced_GetPixelDefectMap\(\)](#) are:

camera: Handle for the camera for which [PicamPixelDefectMap](#) is to be returned.

Valid values are:

- device handle;
- model handle.

NOTE: device and model share the same [PicamPixelDefectMap](#).

Output Parameters

Output parameters for [PicamAdvanced_GetPixelDefectMap\(\)](#) are:

pixel_defect_map: Pointer to the [PicamPixelDefectMap](#) array in which defective pixel information is stored.

When no information is available for the specified camera, this is an array describing zero defects.

NOTE: This memory is allocated by PICam and must be released by calling [PicamAdvanced_DestroyPixelDefectMaps\(\)](#).

Related APIs

For additional information, refer to the following related APIs:

- [PicamAdvanced_DestroyPixelDefectMaps\(\)](#)

Related Structures

For additional information, refer to the following related structures:

- [PicamPixelDefectMap](#)

7.4.4 Camera Parameter Value APIs

This section provides programming information for advanced camera parameter value APIs.

7.4.4.1 *PicamAdvanced_RegisterForIntegerValueChanged()*

Description

`PicamAdvanced_RegisterForIntegerValueChanged()` registers a function to call when the integer value for specified camera parameter has been set, even if it is changed as a result of a different parameter's value being changed.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Registered callbacks are called synchronously from within the thread in which associated parameter values are being set (i.e., called prior to returning from the set operation.)



NOTE:

Parameters whose values have changed due to external influences (e.g., representing the status of camera hardware,) do not result in a callback's being called.

Call `PicamAdvanced_UnregisterForIntegerValueChanged()` to unregister each callback once it is not longer required.

Syntax

The syntax for `PicamAdvanced_RegisterForIntegerValueChanged()` is:

```
PICAM_API PicamAdvanced_RegisterForIntegerValueChanged(
    PicamHandle camera,
    PicamParameter parameter,
    PicamIntegerValueChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_RegisterForIntegerValueChanged()` are:

- camera: Handle for the camera for which the callback is being registered.
- parameter: The parameter for which the callback is being registered.
- changed: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_RegisterForIntegerValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_SetParameterIntegerValue()`;
- `PicamAdvanced_UnregisterForIntegerValueChanged()`.

7.4.4.2 *PicamAdvanced_UnregisterForIntegerValueChanged()*

Description

`PicamAdvanced_UnregisterForIntegerValueChanged()` removes the callback function so that it is no longer called when the integer value for a specified parameter is changed.

Syntax

The syntax for `PicamAdvanced_UnregisterForIntegerValueChanged()` is:

```
PICAM_API PicamAdvanced_UnregisterForIntegerValueChanged(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamIntegerValueChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_UnregisterForIntegerValueChanged()` are:

- camera: Handle for the camera for which the callback is being unregistered.
- parameter: The parameter for which the callback is being unregistered.
- changed: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with `PicamAdvanced_UnregisterForIntegerValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForIntegerValueChanged()`.

7.4.4.3 *PicamAdvanced_RegisterForExtrinsicIntegerValueChanged()*

Description

`PicamAdvanced_RegisterForExtrinsicIntegerValueChanged()` registers a function to call when the integer value for specified camera parameter has changed due to external influences (e.g., representing the status of camera hardware).



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Registered callbacks are called asynchronously on another thread.



NOTE:

`PicamAdvanced_UnregisterForExtrinsicIntegerValueChanged()` must be called to unregister each callback once it is no longer required.

Syntax

The syntax for `PicamAdvanced_RegisterForExtrinsicIntegerValueChanged()` is:

```
PICAM_API PicamAdvanced_RegisterForExtrinsicIntegerValueChanged(  
    PicamHandle device,  
    PicamParameter parameter,  
    PicamIntegerValueChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_RegisterForExtrinsicIntegerValueChanged()` are:

- device: Handle for the device for which the callback is being registered.
- parameter: The parameter for which the callback is being registered.
- changed: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with `PicamAdvanced_RegisterForExtrinsicIntegerValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_UnregisterForExtrinsicIntegerValueChanged()`

7.4.4.4 *PicamAdvanced_UnregisterForExtrinsicIntegerValueChanged()*

Description

`PicamAdvanced_UnregisterForExtrinsicIntegerValueChanged()` removes the callback function so that it is no longer called when the integer value for a specified parameter is changed due to external influences (e.g., representing the status of camera hardware).

Syntax

The syntax for `PicamAdvanced_UnregisterForExtrinsicIntegerValueChanged()` is:

```
PICAM_API PicamAdvanced_UnregisterForExtrinsicIntegerValueChanged(  
    PicamHandle  device,  
    PicamParameter  parameter,  
    PicamIntegerValueChangedCallback  changed );
```

Input Parameters

Input parameters for `PicamAdvanced_UnregisterForExtrinsicIntegerValueChanged()` are:

device: Handle for the device for which the callback is being unregistered.

parameter: The parameter for which the callback is being unregistered.

changed: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_UnregisterForExtrinsicIntegerValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForExtrinsicIntegerValueChanged()`

7.4.4.5 *PicamAdvanced_RegisterForLargeIntegerValueChanged()*

Description

`PicamAdvanced_RegisterForLargeIntegerValueChanged()` registers a function to call when the large integer value for specified camera parameter has been set, even if it is changed as a result of a different parameter's value being changed.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Registered callbacks are called synchronously from within the thread in which associated parameter values are being set (i.e., called prior to returning from the set operation.)



NOTE:

Parameters whose values have change due to external influences (e.g., representing the status of camera hardware,) do not result in a callback's being called.

Call `PicamAdvanced_UnregisterForLargeIntegerValueChanged()` to unregister each callback once it is not longer required.

Syntax

The syntax for `PicamAdvanced_RegisterForLargeIntegerValueChanged()` is:

```
PICAM_API PicamAdvanced_RegisterForLargeIntegerValueChanged(  
                                PicamHandle    camera,  
                                PicamParameter  parameter,  
                                PicamLargeIntegerValueChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_RegisterForLargeIntegerValueChanged()` are:

- camera: Handle for the camera for which the callback is being registered.
- parameter: The parameter for which the callback is being registered.
- changed: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with `PicamAdvanced_RegisterForLargeIntegerValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_SetParameterLargeIntegerValue()`;
- `PicamAdvanced_UnregisterForLargeIntegerValueChanged()`.

7.4.4.6 *PicamAdvanced_UnregisterForLargeIntegerValueChanged()*

Description

`PicamAdvanced_UnregisterForLargeIntegerValueChanged()` removes the callback function so that it is no longer called when the large integer value for a specified parameter is changed.

Syntax

The syntax for `PicamAdvanced_UnregisterForLargeIntegerValueChanged()` is:

```
PICAM_API PicamAdvanced_UnregisterForLargeIntegerValueChanged(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamLargeIntegerValueChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_UnregisterForLargeIntegerValueChanged()` are:

- camera: Handle for the camera for which the callback is being unregistered.
- parameter: The parameter for which the callback is being unregistered.
- changed: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_UnregisterForLargeIntegerValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForLargeIntegerValueChanged()`.

7.4.4.7 *PicamAdvanced_RegisterForFloatingPointValueChanged()*

Description

`PicamAdvanced_RegisterForFloatingPointValueChanged()` registers a function to call when the floating point value for specified camera parameter has been set, even if it is changed as a result of a different parameter's value being changed.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Registered callbacks are called synchronously from within the thread in which associated parameter values are being set (i.e., called prior to returning from the set operation.)



NOTE:

Parameters whose values have change due to external influences (e.g., representing the status of camera hardware,) do not result in a callback's being called.

Call `PicamAdvanced_UnregisterForFloatingPointValueChanged()` to unregister each callback once it is not longer required.

Syntax

The syntax for `PicamAdvanced_RegisterForFloatingPointValueChanged()` is:

```
PICAM_API PicamAdvanced_RegisterForFloatingPointValueChanged(  
                                PicamHandle    camera,  
                                PicamParameter  parameter,  
                                PicamFloatingPointValueChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_RegisterForFloatingPointValueChanged()` are:

- camera: Handle for the camera for which the callback is being registered.
- parameter: The parameter for which the callback is being registered.
- changed: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with `PicamAdvanced_RegisterForFloatingPointValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_SetParameterFloatingPointValue()`;
- `PicamAdvanced_UnregisterForFloatingPointValueChanged()`.

7.4.4.8 *PicamAdvanced_UnregisterForFloatingPointValueChanged()*

Description

`PicamAdvanced_UnregisterForFloatingPointValueChanged()` removes the callback function so that it is no longer called when the floating point value for a specified parameter is changed.

Syntax

The syntax for `PicamAdvanced_UnregisterForFloatingPointValueChanged()` is:

```
PICAM_API PicamAdvanced_UnregisterForFloatingPointValueChanged(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamFloatingPointValueChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_UnregisterForFloatingPointValueChanged()` are:

- camera: Handle for the camera for which the callback is being unregistered.
- parameter: The parameter for which the callback is being unregistered.
- changed: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_UnregisterForFloatingPointValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForFloatingPointValueChanged()`.

7.4.4.9 *PicamAdvanced_RegisterForExtrinsicFloatingPointValueChanged()*

Description

`PicamAdvanced_RegisterForExtrinsicFloatingPointValueChanged()` registers a function to call when the floating point value for specified camera parameter has changed due to external influences (e.g., representing the status of camera hardware).



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Registered callbacks are called asynchronously on another thread.



NOTE:

`PicamAdvanced_UnregisterForExtrinsicFloatingPointValueChanged()` must be called to unregister each callback once it is no longer required.

Syntax

The syntax for `PicamAdvanced_RegisterForExtrinsicFloatingPointValueChanged()` is:

```
PICAM_API PicamAdvanced_RegisterForExtrinsicFloatingPointValueChanged(
                                PicamHandle device,
                                PicamParameter parameter,
                                PicamFloatingPointValueChangedCallback changed );
```

Input Parameters

Input parameters for

`PicamAdvanced_RegisterForExtrinsicFloatingPointValueChanged()` are:

- `device`: Handle for the device for which the callback is being registered.
- `parameter`: The parameter for which the callback is being registered.
- `changed`: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_RegisterForExtrinsicFloatingPointValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_UnregisterForExtrinsicFloatingPointValueChanged()`

7.4.4.10 *PicamAdvanced_UnregisterForExtrinsicFloatingPointValueChanged()*

Description

`PicamAdvanced_UnregisterForExtrinsicFloatingPointValueChanged()` removes the callback function so that it is no longer called when the floating point value for a specified parameter is changed due to external influences (e.g., representing the status of camera hardware).

Syntax

The syntax for

`PicamAdvanced_UnregisterForExtrinsicFloatingPointValueChanged()` is:

```
PICAM_API PicamAdvanced_UnregisterForExtrinsicFloatingPointValueChanged(  
                                PicamHandle  device,  
                                PicamParameter parameter,  
                                PicamIntegerValueChangedCallback changed );
```

Input Parameters

Input parameters for

`PicamAdvanced_UnregisterForExtrinsicFloatingPointValueChanged()` are:

device: Handle for the device for which the callback is being unregistered.

parameter: The parameter for which the callback is being unregistered.

changed: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_UnregisterForExtrinsicFloatingPointValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForExtrinsicFloatingPointValueChanged()`

7.4.4.11 *PicamAdvanced_RegisterForRoisValueChanged()*

Description

`PicamAdvanced_RegisterForRoisValueChanged()` registers a function to call when the value of a specified Rois parameter has been set, even if it is changed as a result of a different parameter's value being changed.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Registered callbacks are called synchronously from within the thread in which associated parameter values are being set (i.e., called prior to returning from the set operation.)



NOTE:

Parameters whose values have change due to external influences (e.g., representing the status of camera hardware,) do not result in a callback's being called.

Call `PicamAdvanced_UnregisterForRoisValueChanged()` to unregister each callback once it is not longer required.

Syntax

The syntax for `PicamAdvanced_RegisterForRoisValueChanged()` is:

```
PICAM_API PicamAdvanced_RegisterForRoisValueChanged(
                                PicamHandle   camera,
                                PicamParameter parameter,
                                PicamRoisValueChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_RegisterForRoisValueChanged()` are:

- camera: Handle for the camera for which the callback is being registered.
- parameter: The parameter for which the callback is being registered.
- changed: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with `PicamAdvanced_RegisterForRoisValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_SetParameterRoisValue()`;
- `PicamAdvanced_UnregisterForRoisValueChanged()`.

7.4.4.12 *PicamAdvanced_UnregisterForRoisValueChanged()*

Description

`PicamAdvanced_UnregisterForRoisValueChanged()` removes the callback function so that it is no longer called when the value of a specified Rois parameter is changed.

Syntax

The syntax for `PicamAdvanced_UnregisterForRoisValueChanged()` is:

```
PICAM_API PicamAdvanced_UnregisterForRoisValueChanged(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamRoisValueChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_UnregisterForRoisValueChanged()` are:

- camera: Handle for the camera for which the callback is being unregistered.
- parameter: The parameter for which the callback is being unregistered.
- changed: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with `PicamAdvanced_UnregisterForRoisValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForRoisValueChanged()`.

7.4.4.13 *PicamAdvanced_RegisterForPulseValueChanged()*

Description

`PicamAdvanced_RegisterForPulseValueChanged()` registers a function to call when the value of a specified gate pulse parameter has been set, even if it is changed as a result of a different parameter's value being changed.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Registered callbacks are called synchronously from within the thread in which associated parameter values are being set (i.e., called prior to returning from the set operation.)



NOTE:

Parameters whose values have change due to external influences (e.g., representing the status of camera hardware,) do not result in a callback's being called.

Call `PicamAdvanced_UnregisterForPulseValueChanged()` to unregister each callback once it is not longer required.

Syntax

The syntax for `PicamAdvanced_RegisterForPulseValueChanged()` is:

```
PICAM_API PicamAdvanced_RegisterForPulseValueChanged(
                                PicamHandle  camera,
                                PicamParameter parameter,
                                PicamPulseValueChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_RegisterForPulseValueChanged()` are:

- camera: Handle for the camera for which the callback is being registered.
- parameter: The parameter for which the callback is being registered.
- changed: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with `PicamAdvanced_RegisterForPulseValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_SetParameterPulseValue()`;
- `PicamAdvanced_UnregisterForPulseValueChanged()`.

7.4.4.14 *PicamAdvanced_UnregisterForPulseValueChanged()*

Description

`PicamAdvanced_UnregisterForPulseValueChanged()` removes the callback function so that it is no longer called when the value of a specified gate pulse parameter is changed.

Syntax

The syntax for `PicamAdvanced_UnregisterForPulseValueChanged()` is:

```
PICAM_API PicamAdvanced_UnregisterForPulseValueChanged(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamPulseValueChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_UnregisterForPulseValueChanged()` are:

- camera: Handle for the camera for which the callback is being unregistered.
- parameter: The parameter for which the callback is being unregistered.
- changed: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with `PicamAdvanced_UnregisterForPulseValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForPulseValueChanged()`.

7.4.4.15 *PicamAdvanced_RegisterForModulationsValueChanged()*

Description

`PicamAdvanced_RegisterForModulationsValueChanged()` registers a function to call when the value of a specified intensifier modulation sequence parameter has been set, even if it is changed as a result of a different parameter's value being changed.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Registered callbacks are called synchronously from within the thread in which associated parameter values are being set (i.e., called prior to returning from the set operation.)



NOTE:

Parameters whose values have change due to external influences (e.g., representing the status of camera hardware,) do not result in a callback's being called.

Call `PicamAdvanced_UnregisterForModulationsValueChanged()` to unregister each callback once it is not longer required.

Syntax

The syntax for `PicamAdvanced_RegisterForModulationsValueChanged()` is:

```
PICAM_API PicamAdvanced_RegisterForModulationsValueChanged(
                                PicamHandle  camera,
                                PicamParameter parameter,
                                PicamModulationsValueChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_RegisterForModulationsValueChanged()` are:

- camera: Handle for the camera for which the callback is being registered.
- parameter: The parameter for which the callback is being registered.
- changed: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with `PicamAdvanced_RegisterForModulationsValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_SetParameterModulationsValue()`;
- `PicamAdvanced_UnregisterForModulationsValueChanged()`.

7.4.4.16 *PicamAdvanced_UnregisterForModulationsValueChanged()*

Description

`PicamAdvanced_UnregisterForModulationsValueChanged()` removes the callback function so that it is no longer called when the value of an intensifier modulation sequence parameter is changed.

Syntax

The syntax for `PicamAdvanced_UnregisterForModulationsValueChanged()` is:

```
PICAM_API PicamAdvanced_UnregisterForModulationsValueChanged(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamModulationsValueChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_UnregisterForModulationsValueChanged()` are:

- camera: Handle for the camera for which the callback is being unregistered.
- parameter: The parameter for which the callback is being unregistered.
- changed: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with `PicamAdvanced_UnregisterForModulationsValueChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForModulationsValueChanged()`.

7.4.4.17 *PicamAdvanced_NotifyWhenStatusParameterValue()*

Description

`PicamAdvanced_NotifyWhenStatusParameterValue()` sets a function to call once when the value of a specified status has been met or when an error has occurred.



NOTE:

Multiple functions may be set. When this is the case, the functions are called in the order in which they have been set.

Set callbacks are called asynchronously from within the thread.

Syntax

The syntax for `PicamAdvanced_NotifyWhenStatusParameterValue()` is:

```
PICAM_API PicamAdvanced_NotifyWhenStatusParameterValue(
    PicamHandle device,
    PicamParameter parameter,
    piint value,
    PicamWhenStatusParameterValueCallback when);
```

Input Parameters

Input parameters for `PicamAdvanced_NotifyWhenStatusParameterValue()` are:

`device`: Handle for the device for which the callback is being set.

`parameter`: The parameter for which the callback is being set.

NOTE: The specified parameter must be a waitable status.

Refer to `Picam_CanWaitForStatusParameter()` for additional information.

`value`: The status value to notify when met.

`when`: The name assigned to the callback function being set.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_NotifyWhenStatusParameterValue()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_CanWaitForStatusParameter()`;
- `PicamAdvanced_CancelNotifyWhenStatusParameterValue()`.

7.4.4.18 *PicamAdvanced_CancelNotifyWhenStatusParameterValue()*

Description

PicamAdvanced_CancelNotifyWhenStatusParameterValue() cancels a function to call once when the value of a specified status has been met or when an error has occurred.

Syntax

The syntax for *PicamAdvanced_CancelNotifyWhenStatusParameterValue()* is:

```
PICAM_API PicamAdvanced_CancelNotifyWhenStatusParameterValue(  
    PicamHandle  device,  
    PicamParameter  parameter,  
    piint  value,  
    PicamWhenStatusParameterValueCallback  when);
```

Input Parameters

Input parameters for *PicamAdvanced_CancelNotifyWhenStatusParameterValue()* are:

device: Handle for the device for which the callback is being canceled.

parameter: The parameter for which the callback is being canceled.

NOTE: The specified parameter must be a waitable status.

Refer to *Picam_CanWaitForStatusParameter()* for additional information.

value: The status value to no longer notify when met.

when: The name assigned to the callback function being canceled.

Output Parameters

There are no output parameters associated with

PicamAdvanced_CancelNotifyWhenStatusParameterValue().

Related APIs

For additional information, refer to the following related APIs:

- *Picam_CanWaitForStatusParameter()*;
- *PicamAdvanced_NotifyWhenStatusParameterValue()*.

7.4.5 Camera Parameter Information APIs

This section provides programming information for advanced camera parameter information APIs.

7.4.5.1 *PicamAdvanced_RegisterForIsRelevantChanged()*

Description

`PicamAdvanced_RegisterForIsRelevantChanged()` registers a function to call when the relevance for a parameter has been changed, even if it is changed as a result of a different parameter's value being changed.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Registered callbacks are called synchronously from within the thread in which associated parameter values are being set (i.e., called prior to returning from the set operation.)

Call `PicamAdvanced_UnregisterForIsRelevantChanged()` to unregister each callback once it is no longer required.

Syntax

The syntax for `PicamAdvanced_RegisterForIsRelevantChanged()` is:

```
PICAM_API PicamAdvanced_RegisterForIsRelevantChanged(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamIsRelevantChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_RegisterForIsRelevantChanged()` are:

camera: Handle for the camera for which the callback is being registered.

parameter: The parameter for which the callback is being registered.

changed: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_RegisterForIsRelevantChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_UnregisterForIsRelevantChanged()`.

7.4.5.2 *PicamAdvanced_UnregisterForIsRelevantChanged()*

Description

`PicamAdvanced_UnregisterForIsRelevantChanged()` removes the callback function so that it is no longer called when the relevance for a parameter has been changed.

Syntax

The syntax for `PicamAdvanced_UnregisterForIsRelevantChanged()` is:

```
PICAM_API PicamAdvanced_UnregisterForIsRelevantChanged(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamIsRelevantChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_UnregisterForIsRelevantChanged()` are:

- camera: Handle for the camera for which the callback is being unregistered.
- parameter: The parameter for which the callback is being unregistered.
- changed: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with `PicamAdvanced_UnregisterForIsRelevantChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForIsRelevantChanged()`.

7.4.5.3 *PicamAdvanced_RegisterForValueAccessChanged()*

Description

`PicamAdvanced_RegisterForValueAccessChanged()` registers a function to call when the value access for a parameter has been changed, even if it is changed as a result of a different parameter's value being changed.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Registered callbacks are called synchronously from within the thread in which associated parameter values are being set (i.e., called prior to returning from the set operation.)

Call `PicamAdvanced_UnregisterForValueAccessChanged()` to unregister each callback once it is no longer required.

Syntax

The syntax for `PicamAdvanced_RegisterForValueAccessChanged()` is:

```
PICAM_API PicamAdvanced_RegisterForValueAccessChanged(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamValueAccessChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_RegisterForValueAccessChanged()` are:

`camera`: Handle for the camera for which the callback is being registered.

`parameter`: The parameter for which the callback is being registered.

`changed`: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_RegisterForValueAccessChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_UnregisterForValueAccessChanged()`.

7.4.5.4 *PicamAdvanced_UnregisterForValueAccessChanged()*

Description

`PicamAdvanced_UnregisterForValueAccessChanged()` removes the callback function so that it is no longer called when the value access for a parameter has been changed.

Syntax

The syntax for `PicamAdvanced_UnregisterForValueAccessChanged()` is:

```
PICAM_API PicamAdvanced_UnregisterForValueAccessChanged(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamValueAccessChangedCallback changed );
```

Input Parameters

Input parameters for `PicamAdvanced_UnregisterForValueAccessChanged()` are:

- camera: Handle for the camera for which the callback is being unregistered.
- parameter: The parameter for which the callback is being unregistered.
- changed: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with `PicamAdvanced_UnregisterForValueAccessChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForValueAccessChanged()`.

7.4.5.5 *PicamAdvanced_GetParameterDynamics()*

Description

`PicamAdvanced_GetParameterDynamics()` returns the dynamics for a specified parameter.

Syntax

The syntax for `PicamAdvanced_GetParameterDynamics()` is:

```
PICAM_API PicamAdvanced_GetParameterDynamics(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamDynamicsMask* dynamics );
```

Input Parameters

Input parameters for `PicamAdvanced_GetParameterDynamics()` are:

`camera`: Handle for the camera for which the dynamics information is to be returned.

`parameter`: The parameter for which dynamics information is to be returned.

Output Parameters

Output parameters for `PicamAdvanced_GetParameterDynamics()` are:

`dynamics`: Pointer to the memory location in which the dynamics information is stored.

7.4.5.6 *PicamAdvanced_GetParameterExtrinsicDynamics()*

Description

`PicamAdvanced_GetParameterExtrinsicDynamics()` returns the dynamics for a specified parameter that can change due to external influences (e.g., representing the status of camera hardware).

Syntax

The syntax for `PicamAdvanced_GetParameterExtrinsicDynamics()` is:

```
PICAM_API PicamAdvanced_GetParameterExtrinsicDynamics(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamDynamicsMask* extrinsic);
```

Input Parameters

Input parameters for `PicamAdvanced_GetParameterExtrinsicDynamics()` are:

`camera`: Handle for the camera for which the extrinsic dynamics information is to be returned.

`parameter`: The parameter for which extrinsic dynamics information is to be returned.

Output Parameters

Output parameters for `PicamAdvanced_GetParameterExtrinsicDynamics()` are:

`extrinsic`: Pointer to the memory location in which the extrinsic dynamics information is stored.

7.4.6 Camera Parameter Constraints APIs

This section provides programming information for advanced camera parameter constraint APIs.

7.4.6.1 *PicamAdvanced_GetParameterCollectionConstraints()*

Description

`PicamAdvanced_GetParameterCollectionConstraints()` returns an allocated array in which all collection constraints for a specified camera parameter are stored.

Syntax

The syntax for `PicamAdvanced_GetParameterCollectionConstraints()` is:

```
PICAM_API PicamAdvanced_GetParameterCollectionConstraints(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamCollectionConstraint** constraint_array,  
    piint* constraint_count );
```

Input Parameters

Input parameters for `PicamAdvanced_GetParameterCollectionConstraints()` are:

`camera`: Handle for the camera for which the collection constraint information is to be returned.

`parameter`: The parameter for which collection constraint information is to be returned.

Output Parameters

Output parameters for `PicamAdvanced_GetParameterCollectionConstraints()` are:

`constraint_array`: Pointer to the array in which collection constraint information is stored.

NOTE: This memory is allocated by PICam and must be released by calling `Picam_DestroyCollectionConstraints()`

`constraint_count`: Pointer to the memory location in which the number of constraints is stored.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_DestroyCollectionConstraints()`.

7.4.6.2 *PicamAdvanced_RegisterForDependentCollectionConstraintChanged()*

Description

`PicamAdvanced_RegisterForDependentCollectionConstraintChanged()` registers a function to call when any dependent collection constraint has been changed due to the setting of a DIFFERENT parameter's value.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Registered callbacks are called synchronously from within the thread in which associated parameter values are being set (i.e., called prior to returning from the set operation.)

Call `PicamAdvanced_UnregisterForDependentCollectionConstraintChanged()` to unregister each callback once it is not longer required.

Syntax

The syntax for

`PicamAdvanced_RegisterForDependentCollectionConstraintChanged()` is:

```
PICAM_API PicamAdvanced_RegisterForDependentCollectionConstraint
Changed(
                                PicamHandle  camera,
                                PicamParameter  parameter,
                                PicamDependentCollectionConstraintChangedCallback  changed );
```

Input Parameters

Input parameters for

`PicamAdvanced_RegisterForDependentCollectionConstraintChanged()` are:

- camera: Handle for the camera for which the callback is being registered.
- parameter: The parameter for which the callback is being registered.
- changed: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_RegisterForDependentCollectionConstraintChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_UnregisterForDependentCollectionConstraintChanged()`.

7.4.6.3 *PicamAdvanced_UnregisterForDependentCollectionConstraintChanged()*

Description

`PicamAdvanced_UnregisterForDependentCollectionConstraintChanged()` removes the callback function so that it is no longer called when any dependent collection constraint has been changed.

Syntax

The syntax for

`PicamAdvanced_UnregisterForDependentCollectionConstraintChanged()` is:

```
PICAM_API PicamAdvanced_UnregisterForDependentCollectionConstraint
Changed(
                                PicamHandle camera,
                                PicamParameter parameter,
                                PicamDependentCollectionConstraintChangedCallback changed );
```

Input Parameters

Input parameters for

`PicamAdvanced_UnregisterForDependentCollectionConstraintChanged()` are:

- `camera`: Handle for the camera for which the callback is being unregistered.
- `parameter`: The parameter for which the callback is being unregistered.
- `changed`: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_UnregisterForDependentCollectionConstraintChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForDependentCollectionConstraintChanged()`.

7.4.6.4 *PicamAdvanced_GetParameterRangeConstraints()*

Description

PicamAdvanced_GetParameterRangeConstraints() returns an allocated array in which all range constraints for a specified camera parameter are stored.

Syntax

The syntax for *PicamAdvanced_GetParameterRangeConstraints()* is:

```
PICAM_API PicamAdvanced_GetParameterRangeConstraints(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamRangeConstraint** constraint_array,  
    piint* constraint_count );
```

Input Parameters

Input parameters for *PicamAdvanced_GetParameterRangeConstraints()* are:

camera: Handle for the camera for which the range constraint information is to be returned.

parameter: The parameter for which range constraint information is to be returned.

Output Parameters

Output parameters for *PicamAdvanced_GetParameterRangeConstraints()* are:

constraint_array: Pointer to the array in which range constraint information is stored.

NOTE: This memory is allocated by PICam and must be released by calling *Picam_DestroyRangeConstraints()*

constraint_count: Pointer to the memory location in which the number of constraints is stored.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_DestroyRangeConstraints()*.

7.4.6.5 *PicamAdvanced_RegisterForDependentRangeConstraintChanged()*

Description

[PicamAdvanced_RegisterForDependentRangeConstraintChanged\(\)](#) registers a function to call when any dependent range constraint has been changed due to the setting of a DIFFERENT parameter's value.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Registered callbacks are called synchronously from within the thread in which associated parameter values are being set (i.e., called prior to returning from the set operation.)

Call [PicamAdvanced_UnregisterForDependentRangeConstraintChanged\(\)](#) to unregister each callback once it is not longer required.

Syntax

The syntax for [PicamAdvanced_RegisterForDependentRangeConstraintChanged\(\)](#) is:

```
PICAM_API PicamAdvanced_RegisterForDependentRangeConstraintChanged(  
                                                PicamHandle camera,  
                                                PicamParameter parameter,  
                                                PicamDependentRangeConstraintChangedCallback changed );
```

Input Parameters

Input parameters for

[PicamAdvanced_RegisterForDependentRangeConstraintChanged\(\)](#) are:

camera: Handle for the camera for which the callback is being registered.

parameter: The parameter for which the callback is being registered.

changed: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with

[PicamAdvanced_RegisterForDependentRangeConstraintChanged\(\)](#).

Related APIs

For additional information, refer to the following related APIs:

- [PicamAdvanced_UnregisterForDependentRangeConstraintChanged\(\)](#).

7.4.6.6 *PicamAdvanced_UnregisterForDependentRangeConstraintChanged()*

Description

`PicamAdvanced_UnregisterForDependentRangeConstraintChanged()` removes the callback function so that it is no longer called when any dependent range constraint has been changed.

Syntax

The syntax for `PicamAdvanced_UnregisterForDependentRangeConstraintChanged()` is:

```
PICAM_API PicamAdvanced_UnregisterForDependentRangeConstraintChanged(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamDependentRangeConstraintChangedCallback changed );
```

Input Parameters

Input parameters for

`PicamAdvanced_UnregisterForDependentRangeConstraintChanged()` are:

- camera: Handle for the camera for which the callback is being unregistered.
- parameter: The parameter for which the callback is being unregistered.
- changed: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_UnregisterForDependentRangeConstraintChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForDependentRangeConstraintChanged()`.

7.4.6.7 *PicamAdvanced_GetParameterRoisConstraints()*

Description

[PicamAdvanced_GetParameterRoisConstraints\(\)](#) returns an allocated array in which all Rois constraints for a specified camera parameter are stored.

Syntax

The syntax for [PicamAdvanced_GetParameterRoisConstraints\(\)](#) is:

```
PICAM_API PicamAdvanced_GetParameterRoisConstraints(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamRoisConstraint** constraint_array,  
    piint* constraint_count );
```

Input Parameters

Input parameters for [PicamAdvanced_GetParameterRoisConstraints\(\)](#) are:

camera: Handle for the camera for which the Rois constraint information is to be returned.

parameter: The parameter for which Rois constraint information is to be returned.

Output Parameters

Output parameters for [PicamAdvanced_GetParameterRoisConstraints\(\)](#) are:

constraint_array: Pointer to the array in which Rois constraint information is stored.

NOTE: This memory is allocated by PICam and must be released by calling [Picam_DestroyRoisConstraints\(\)](#)

constraint_count: Pointer to the memory location in which the number of constraints is stored.

Related APIs

For additional information, refer to the following related APIs:

- [Picam_DestroyRoisConstraints\(\)](#)

7.4.6.8 *PicamAdvanced_RegisterForDependentRoisConstraintChanged()*

Description

[PicamAdvanced_RegisterForDependentRoisConstraintChanged\(\)](#) registers a function to call when any dependent Rois constraint has been changed due to the setting of a DIFFERENT parameter's value.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Registered callbacks are called synchronously from within the thread in which associated parameter values are being set (i.e., called prior to returning from the set operation.)

Call [PicamAdvanced_UnregisterForDependentRoisConstraintChanged\(\)](#) to unregister each callback once it is no longer required.

Syntax

The syntax for [PicamAdvanced_RegisterForDependentRoisConstraintChanged\(\)](#) is:

```
PICAM_API PicamAdvanced_RegisterForDependentRoisConstraintChanged(  
                                PicamHandle camera,  
                                PicamParameter parameter,  
                                PicamDependentRoisConstraintChangedCallback changed );
```

Input Parameters

Input parameters for [PicamAdvanced_RegisterForDependentRoisConstraintChanged\(\)](#) are:

- camera: Handle for the camera for which the callback is being registered.
- parameter: The parameter for which the callback is being registered.
- changed: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with

[PicamAdvanced_RegisterForDependentRoisConstraintChanged\(\)](#).

Related APIs

For additional information, refer to the following related APIs:

- [PicamAdvanced_UnregisterForDependentRoisConstraintChanged\(\)](#)

7.4.6.9 *PicamAdvanced_UnregisterForDependentRoisConstraintChanged()*

Description

`PicamAdvanced_UnregisterForDependentRoisConstraintChanged()` removes the callback function so that it is no longer called when any dependent Rois constraint has been changed.

Syntax

The syntax for `PicamAdvanced_UnregisterForDependentRoisConstraintChanged()` is:

```
PICAM_API PicamAdvanced_UnregisterForDependentRoisConstraintChanged(  
                                PicamHandle camera,  
                                PicamParameter parameter,  
                                PicamDependentRoisConstraintChangedCallback changed );
```

Input Parameters

Input parameters for

`PicamAdvanced_UnregisterForDependentRoisConstraintChanged()` are:

- camera: Handle for the camera for which the callback is being unregistered.
- parameter: The parameter for which the callback is being unregistered.
- changed: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_UnregisterForDependentRoisConstraintChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForDependentRoisConstraintChanged()`.

7.4.6.10 *PicamAdvanced_GetParameterPulseConstraints()*

Description

PicamAdvanced_GetParameterPulseConstraints() returns an allocated array in which all Pulse constraints for a specified camera parameter are stored.

Syntax

The syntax for *PicamAdvanced_GetParameterPulseConstraints()* is:

```
PICAM_API PicamAdvanced_GetParameterPulseConstraints(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamPulseConstraint** constraint_array,  
    piint* constraint_count );
```

Input Parameters

Input parameters for *PicamAdvanced_GetParameterPulseConstraints()* are:

- camera: Handle for the camera for which the Pulse constraint information is to be returned.
- parameter: The parameter for which Pulse constraint information is to be returned.

Output Parameters

Output parameters for *PicamAdvanced_GetParameterPulseConstraints()* are:

- constraint_array: Pointer to the array in which Pulse constraint information is stored.
NOTE: This memory is allocated by PICam and must be released by calling *Picam_DestroyPulseConstraints()*
- constraint_count: Pointer to the memory location in which the number of constraints is stored.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_DestroyPulseConstraints()*

7.4.6.11 *PicamAdvanced_RegisterForDependentPulseConstraintChanged()*

Description

[PicamAdvanced_RegisterForDependentPulseConstraintChanged\(\)](#) registers a function to call when any dependent Pulse constraint has been changed due to the setting of a DIFFERENT parameter's value.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Registered callbacks are called synchronously from within the thread in which associated parameter values are being set (i.e., called prior to returning from the set operation.)

Call [PicamAdvanced_UnregisterForDependentPulseConstraintChanged\(\)](#) to unregister each callback once it is not longer required.

Syntax

The syntax for [PicamAdvanced_RegisterForDependentPulseConstraintChanged\(\)](#) is:

```
PICAM_API PicamAdvanced_RegisterForDependentPulseConstraintChanged(  
                                PicamHandle camera,  
                                PicamParameter parameter,  
                                PicamDependentPulseConstraintChangedCallback changed );
```

Input Parameters

Input parameters for

[PicamAdvanced_RegisterForDependentPulseConstraintChanged\(\)](#) are:

camera: Handle for the camera for which the callback is being registered.

parameter: The parameter for which the callback is being registered.

changed: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with

[PicamAdvanced_RegisterForDependentPulseConstraintChanged\(\)](#).

Related APIs

For additional information, refer to the following related APIs:

- [PicamAdvanced_UnregisterForDependentPulseConstraintChanged\(\)](#).

7.4.6.12 *PicamAdvanced_UnregisterForDependentPulseConstraintChanged()*

Description

`PicamAdvanced_UnregisterForDependentPulseConstraintChanged()` removes the callback function so that it is no longer called when any dependent Pulse constraint has been changed.

Syntax

The syntax for `PicamAdvanced_UnregisterForDependentPulseConstraintChanged()` is:

```
PICAM_API PicamAdvanced_UnregisterForDependentPulseConstraintChanged(  
    PicamHandle camera,  
    PicamParameter parameter,  
    PicamDependentPulseConstraintChangedCallback changed );
```

Input Parameters

Input parameters for

`PicamAdvanced_UnregisterForDependentPulseConstraintChanged()` are:

- camera: Handle for the camera for which the callback is being unregistered.
- parameter: The parameter for which the callback is being unregistered.
- changed: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_UnregisterForDependentPulseConstraintChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForDependentPulseConstraintChanged()`.

7.4.6.13 *PicamAdvanced_GetParameterModulationsConstraints()*

Description

PicamAdvanced_GetParameterModulationsConstraints() returns an allocated array in which all Modulation constraints for a specified camera parameter are stored.

Syntax

The syntax for *PicamAdvanced_GetParameterModulationsConstraints()* is:

```
PICAM_API PicamAdvanced_GetParameterModulationsConstraints(  
    PicamHandle camera,  
    PicamParameter parameter,  
    const PicamModulationsConstraint** constraint_array,  
    piint* constraint_count );
```

Input Parameters

Input parameters for *PicamAdvanced_GetParameterModulationsConstraints()* are:

- camera: Handle for the camera for which the Modulation constraint information is to be returned.
- parameter: The parameter for which Modulation constraint information is to be returned.

Output Parameters

Output parameters for *PicamAdvanced_GetParameterModulationsConstraints()* are:

- constraint_array: Pointer to the array in which Modulation constraint information is stored.
NOTE: This memory is allocated by PICam and must be released by calling *Picam_DestroyModulationsConstraints()*
- constraint_count: Pointer to the memory location in which the number of constraints is stored.

Related APIs

For additional information, refer to the following related APIs:

- *Picam_DestroyModulationsConstraints()*

7.4.6.14 *PicamAdvanced_RegisterForDependentModulationsConstraintChanged()*

Description

`PicamAdvanced_RegisterForDependentModulationsConstraintChanged()` registers a function to call when any dependent Modulation constraint has been changed due to the setting of a DIFFERENT parameter's value.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Registered callbacks are called synchronously from within the thread in which associated parameter values are being set (i.e., called prior to returning from the set operation.)

Call `PicamAdvanced_UnregisterForDependentModulationsConstraintChanged()` to unregister each callback once it is not longer required.

Syntax

The syntax for

`PicamAdvanced_RegisterForDependentModulationsConstraintChanged()` is:

```
PICAM_API PicamAdvanced_RegisterForDependentModulationsConstraint
Changed(
                                PicamHandle  camera,
                                PicamParameter  parameter,
                                PicamDependentModulationsConstraintChangedCallback  changed );
```

Input Parameters

Input parameters for

`PicamAdvanced_RegisterForDependentModulationsConstraintChanged()` are:

- camera: Handle for the camera for which the callback is being registered.
- parameter: The parameter for which the callback is being registered.
- changed: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_RegisterForDependentModulationsConstraintChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_UnregisterForDependentModulationsConstraintChanged()`

7.4.6.15 *PicamAdvanced_UnregisterForDependentModulationsConstraintChanged()*

Description

`PicamAdvanced_UnregisterForDependentModulationsConstraintChanged()` removes the callback function so that it is no longer called when any dependent Modulation constraint has been changed.

Syntax

The syntax for

`PicamAdvanced_UnregisterForDependentModulationsConstraintChanged()` is:

```
PICAM_API PicamAdvanced_UnregisterForDependentModulations
          ConstraintChanged(
                                PicamHandle camera,
                                PicamParameter parameter,
                                PicamDependentModulationsConstraintChangedCallback changed );
```

Input Parameters

Input parameters for

`PicamAdvanced_UnregisterForDependentModulationsConstraintChanged()` are:

- camera: Handle for the camera for which the callback is being unregistered.
- parameter: The parameter for which the callback is being unregistered.
- changed: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_UnregisterForDependentModulationsConstraintChanged()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForDependentModulationsConstraintChanged()`.

7.4.7 Camera Commitment APIs

This section provides programming information for advanced camera commitment APIs.

7.4.7.1 *Picam_DestroyValidationResult()*

Description

`Picam_DestroyValidationResult()` releases memory that has been allocated by PICam for use by `result`.

If `result` is null, calling `Picam_DestroyValidationResult()` has no effect.

Syntax

The syntax for `Picam_DestroyValidationResult()` is:

```
PICAM_API Picam_DestroyValidationResult(  
    const PicamValidationResult* result );
```

Input Parameters

Input parameters for `Picam_DestroyValidationResult()` are:

`result`: Pointer to the array that is to be released.

Output Parameters

There are no output parameters associated with `Picam_DestroyValidationResult()`.

Related Structures

For additional information, refer to the following related structures:

- `PicamValidationResult`.

7.4.7.2 *Picam_DestroyValidationResults()*

Description

`Picam_DestroyValidationResults()` releases memory that has been allocated by PICam for use by `results`.

If `results` is null, calling `Picam_DestroyValidationResults()` has no effect.

Syntax

The syntax for `Picam_DestroyValidationResults()` is:

```
PICAM_API Picam_DestroyValidationResults(  
    const PicamValidationResults* results );
```

Input Parameters

Input parameters for `Picam_DestroyValidationResults()` are:

`results`: Pointer to the array that is to be released.

Output Parameters

There are not output parameters associated with `Picam_DestroyValidationResults()`.

Related Structures

For additional information, refer to the following related structures:

- `PicamValidationResults`.

7.4.7.3 *PicamAdvanced_ValidateParameter()*

Description

`PicamAdvanced_ValidateParameter()` validates a single, specified parameter against all associated constraints and returns the results.

Syntax

The syntax for `PicamAdvanced_ValidateParameter()` is:

```
PICAM_API PicamAdvanced_ValidateParameter(  
    PicamHandle  model,  
    PicamParameter  parameter,  
    const PicamValidationResult**  result );
```

Input Parameters

Input parameters for `PicamAdvanced_ValidateParameter()` are:

`model`: Handle for the model for which the parameter is being validated.

`parameter`: The parameter being validated.

Output Parameters

Output parameters for `PicamAdvanced_ValidateParameter()` are:

`result`: Pointer to the array in which the validation results for all constraints are stored.

NOTE: This memory is allocated by PICam and must be released by calling `Picam_DestroyValidationResult()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_DestroyValidationResult()`.

Related Structures

For additional information, refer to the following related structures:

- `PicamValidationResult`.

7.4.7.4 *PicamAdvanced_ValidateParameters()*

Description

`PicamAdvanced_ValidateParameters()` validates all parameters against all associated constraints and returns the results.

Syntax

The syntax for `PicamAdvanced_ValidateParameters()` is:

```
PICAM_API PicamAdvanced_ValidateParameters(  
                                PicamHandle model,  
                                const PicamValidationResults** results );
```

Input Parameters

Input parameters for `PicamAdvanced_ValidateParameters()` are:

`model`: Handle for the model for which all parameters are being validated.

Output Parameters

Output parameters for `PicamAdvanced_ValidateParameters()` are:

`result`: Pointer to the array in which the validation results for all constraints are stored.

NOTE: This memory is allocated by PICam and must be released by calling `Picam_DestroyValidationResults()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_DestroyValidationResults()`.

Related Structures

For additional information, refer to the following related structures:

- `PicamValidationResults`.

7.4.7.5 *Picam_DestroyDependentValidationResult()*

Description

`Picam_DestroyDependentValidationResult()` releases memory that has been allocated by PICam for use by `result`.

If `result` is null, calling `Picam_DestroyDependentValidationResult()` has no effect.

Syntax

The syntax for `Picam_DestroyDependentValidationResult()` is:

```
PICAM_API Picam_DestroyDependentValidationResult(  
    const PicamDependentValidationResult* result );
```

Input Parameters

Input parameters for `Picam_DestroyDependentValidationResult()` are:

`result`: Pointer to the array that is to be released.

Output Parameters

There are no output parameters associated with `Picam_DestroyDependentValidationResult()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamDependentValidationResult`.

7.4.7.6 *PicamAdvanced_ValidateDependentParameter()*

Description

`PicamAdvanced_ValidateDependentParameter()` validates all parameters of a specified model whose constraints are dependent on a specified parameter.

Syntax

The syntax for `PicamAdvanced_ValidateDependentParameter()` is:

```
PICAM_API PicamAdvanced_ValidateDependentParameter(  
    PicamHandle  model,  
    PicamParameter  parameter,  
    const PicamDependentValidationResult**  result );
```

Input Parameters

Input parameters for `PicamAdvanced_ValidateDependentParameter()` are:

`model`: Handle for the model for which all dependent parameters are being validated.

`parameter`: The parameter on which all constraints being validated are dependent.

Output Parameters

Output parameters for `PicamAdvanced_ValidateDependentParameter()` are:

`result`: Pointer to the array in which the validation results for all constraints are stored.

NOTE: This memory is allocated by PICam and must be released by calling `Picam_DestroyDependentValidationResult()`.

Related APIs

For additional information, refer to the following related APIs:

- `Picam_DestroyDependentValidationResult()`

7.4.7.7 *PicamAdvanced_CommitParametersToCameraDevice()*

Description

`PicamAdvanced_CommitParametersToCameraDevice()` attempts to configure a camera device with the set of parameter values stored in `model`.



NOTE:

If this action leads to a camera device error, the action fails and the camera device remains untouched.

Syntax

The syntax for `PicamAdvanced_CommitParametersToCameraDevice()` is:

```
PICAM_API PicamAdvanced_CommitParametersToCameraDevice(  
    PicamHandle model );
```

Input Parameters

Input parameters for `PicamAdvanced_CommitParametersToCameraDevice()` are:

`model`: Handle for the model for which all parameters are to be committed.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_CommitParametersToCameraDevice()`.

7.4.7.8 *PicamAdvanced_RefreshParameterFromCameraDevice()*

Description

`PicamAdvanced_RefreshParameterFromCameraDevice()` updates a single parameter's value stored in `model` with the value from the connected camera device.

Syntax

The syntax for `PicamAdvanced_RefreshParameterFromCameraDevice()` is:

```
PICAM_API PicamAdvanced_RefreshParameterFromCameraDevice(  
    PicamHandle  model,  
    PicamParameter  parameter );
```

Input Parameters

Input parameters for `PicamAdvanced_RefreshParameterFromCameraDevice()` are:

`model`: Handle for the model for which the parameter's value is to be overwritten.

`parameter`: The parameter for which the value is to be overwritten.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_RefreshParameterFromCameraDevice()`.

7.4.7.9 *PicamAdvanced_RefreshParametersFromCameraDevice()*

Description

`PicamAdvanced_RefreshParametersFromCameraDevice()` updates all parameter values stored in `model` with values from the connected camera device.

Syntax

The syntax for `PicamAdvanced_RefreshParametersFromCameraDevice()` is:

```
PICAM_API PicamAdvanced_RefreshParametersFromCameraDevice(  
    PicamHandle model );
```

Input Parameters

Input parameters for `PicamAdvanced_RefreshParametersFromCameraDevice()` are:

`model`: Handle for the model for which all parameter values are to be overwritten.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_RefreshParametersFromCameraDevice()`.

7.4.8 Acquisition Setup APIs

This section provides programming information about advanced acquisition setup APIs.

7.4.8.1 *PicamAdvanced_GetAcquisitionBuffer()*

Description

`PicamAdvanced_GetAcquisitionBuffer()` returns the user-allocated buffer to be used during data acquisition.

Syntax

The syntax for `PicamAdvanced_GetAcquisitionBuffer()` is:

```
PICAM_API PicamAdvanced_GetAcquisitionBuffer(  
    PicamHandle device,  
    PicamAcquisitionBuffer* buffer );
```

Input Parameters

Input parameters for `PicamAdvanced_GetAcquisitionBuffer()` are:

`device`: Handle for the device to which the data acquisition buffer is allocated.

Output Parameters

Output parameters for `PicamAdvanced_GetAcquisitionBuffer()` are:

`buffer`: Pointer to the user-allocated data acquisition buffer.
If no buffer has been created/allocated, this points to a null buffer with zero size.

7.4.8.2 *PicamAdvanced_SetAcquisitionBuffer()*

Description

`PicamAdvanced_SetAcquisitionBuffer()` assigns a user-allocated buffer to a specific device.

Syntax

The syntax for `PicamAdvanced_SetAcquisitionBuffer()` is:

```
PICAM_API PicamAdvanced_SetAcquisitionBuffer(  
                                PicamHandle device,  
                                const PicamAcquisitionBuffer* buffer );
```

Input Parameters

Input parameters for `PicamAdvanced_SetAcquisitionBuffer()` are:

`device`: Handle for the device to which the data acquisition buffer is to be allocated.

Output Parameters

Output parameters for `PicamAdvanced_SetAcquisitionBuffer()` are:

`buffer`: Pointer to the user-allocated data acquisition buffer.
To clear this buffer, point to null with zero size.
This buffer can be used to create a circular buffer.

7.4.9 Acquisition Notification APIs

This section provides programming information for advanced acquisition notification APIs.

7.4.9.1 *PicamAdvanced_RegisterForAcquisitionUpdated()*

Description

`PicamAdvanced_RegisterForAcquisitionUpdated()` registers a function to call during data acquisition when:

- New data are available, or
- A change in acquisition status has occurred.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Callbacks are called asynchronously from another thread, but are serialized on that thread. This means that additional notifications do not occur simultaneously, but occur after each callback returns.

Call `PicamAdvanced_UnregisterForAcquisitionUpdated()` to unregister each callback once it is no longer required.

Syntax

The syntax for `PicamAdvanced_RegisterForAcquisitionUpdated()` is:

```
PICAM_API PicamAdvanced_RegisterForAcquisitionUpdated(  
                                                PicamHandle device,  
                                                PicamAcquisitionUpdatedCallback updated );
```

Input Parameters

Input parameters for `PicamAdvanced_RegisterForAcquisitionUpdated()` are:

device: Handle for the device for which the callback is being registered.

changed: The name assigned to the callback function being registered.

Output Parameters

There are no output parameters associated with `PicamAdvanced_RegisterForAcquisitionUpdated()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_UnregisterForAcquisitionUpdated()`.

7.4.9.2 *PicamAdvanced_UnregisterForAcquisitionUpdated()*

Description

`PicamAdvanced_UnregisterForAcquisitionUpdated()` removes the callback function so that it is no longer called during data acquisition.

Syntax

The syntax for `PicamAdvanced_UnregisterForAcquisitionUpdated()` is:

```
PICAM_API PicamAdvanced_UnregisterForAcquisitionUpdated(  
                                PicamHandle device,  
                                PicamAcquisitionUpdatedCallback updated );
```

Input Parameters

Input parameters for `PicamAdvanced_UnregisterForAcquisitionUpdated()` are:

device: Handle for the device for which the callback is being unregistered.

changed: The name assigned to the callback function being unregistered.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_UnregisterForAcquisitionUpdated()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForAcquisitionUpdated()`

7.4.10 Acquisition State Notification APIs

This section provides programming information for advanced acquisition state notification APIs.

7.4.10.1 *PicamAdvanced_CanRegisterForAcquisitionStateUpdated()*

Description

`PicamAdvanced_CanRegisterForAcquisitionStateUpdated()` determines if an acquisition state can be detected.

Syntax

The syntax for `PicamAdvanced_CanRegisterForAcquisitionStateUpdated()` is:

```
PICAM_API PicamAdvanced_CanRegisterForAcquisitionStateUpdated(  
                                PicamHandle device,  
                                PicamAcquisitionState state,  
                                pibln* detectable );
```

Input Parameters

Input parameters for `PicamAdvanced_CanRegisterForAcquisitionStateUpdated()` are:

`device`: Handle for the device under test.

`state`: Specifies the acquisition state to be queried for detectability.

Output Parameters

Output parameters for `PicamAdvanced_CanRegisterForAcquisitionStateUpdated()` are:

`detectable`: Pointer to the test results. Indicates if the specified acquisition state is detectable.

Valid values are:

- `TRUE`
Indicates the specified acquisition state is detectable.
- `FALSE`
Indicates the specified acquisition state is not detectable.

7.4.10.2 *PicamAdvanced_RegisterForAcquisitionStateUpdated()*

Description

`PicamAdvanced_RegisterForAcquisitionStateUpdated()` registers a function to call during data acquisition when the camera transitions to an acquisition state.



NOTE:

Multiple functions may be registered. When this is the case, the functions are called in the order in which they have been registered.

Callbacks are called asynchronously from another thread, but are serialized on that thread. This means that additional notifications do not occur simultaneously, but occur after each callback returns.

Syntax

The syntax for `PicamAdvanced_RegisterForAcquisitionStateUpdated()` is:

```
PICAM_API PicamAdvanced_RegisterForAcquisitionStateUpdated(  
    PicamHandle device,  
    PicamAcquisitionState state,  
    PicamAcquisitionStateUpdatedCallback updated );
```

Input Parameters

Input parameters for `PicamAdvanced_RegisterForAcquisitionStateUpdated()` are:

`device`: Handle for the device for which the callback is being registered.

`state`: Specifies the acquisition state to detect.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_RegisterForAcquisitionStateUpdated()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_UnregisterForAcquisitionStateUpdated()`

7.4.10.3 *PicamAdvanced_UnregisterForAcquisitionStateUpdated()*

Description

`PicamAdvanced_UnregisterForAcquisitionStateUpdated()` removes the callback function so that it is no longer called during data acquisition.

Syntax

The syntax for `PicamAdvanced_UnregisterForAcquisitionStateUpdated()` is:

```
PICAM_API PicamAdvanced_UnregisterForAcquisitionStateUpdated(  
                                PicamHandle  device,  
                                PicamAcquisitionState  state,  
                                PicamAcquisitionStateUpdatedCallback  updated );
```

Input Parameters

Input parameters for `PicamAdvanced_UnregisterForAcquisitionStateUpdated()` are:

device: Handle for the device for which the callback is being unregistered.

state: Specifies the acquisition state to detect no longer.

Output Parameters

There are no output parameters associated with

`PicamAdvanced_UnregisterForAcquisitionStateUpdated()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamAdvanced_RegisterForAcquisitionStateUpdated()`

7.4.11 Acquisition Control APIs

This section provides programming information for advance acquisition control APIs.

7.4.11.1 *PicamAdvanced_HasAcquisitionBufferOverrun()*

Description

`PicamAdvanced_HasAcquisitionBufferOverrun()` determines if a user-allocated circular buffer has overflowed.

Syntax

The syntax for `PicamAdvanced_HasAcquisitionBufferOverrun()` is:

```
PICAM_API PicamAdvanced_HasAcquisitionBufferOverrun(  
    PicamHandle device,  
    pibln*   overrun );
```

Input Parameters

Input parameters for `PicamAdvanced_HasAcquisitionBufferOverrun()` are:

device: Handle for the device for which the status of the associated user-allocated circular buffer is being tested.

Output Parameters

Output parameters for `PicamAdvanced_HasAcquisitionBufferOverrun()` are:

overran: Pointer to the results.

Indicates if the user-allocated circular data buffer has overflowed.

Valid values are:

- **TRUE**
Indicates that the buffer has overflowed.
- **FALSE**
Indicates that the buffer has not overflowed.

7.4.11.2 *PicamAdvanced_ClearReadoutCountOnline()*

Description

`PicamAdvanced_ClearReadoutCountOnline()` tries to set the readout count to 0 [zero] while the camera is running.

Syntax

The syntax for `PicamAdvanced_ClearReadoutCountOnline()` is:

```
PICAM_API PicamAdvanced_ClearReadoutCountOnline(  
                                PicamHandle device,  
                                pibln* cleared );
```

Input Parameters

Input parameters for `PicamAdvanced_ClearReadoutCountOnline()` are:

`device`: Handle for the device for which the status of the online readout count is being tested.

Output Parameters

Output parameters for `PicamAdvanced_ClearReadoutCountOnline()` are:

`cleared`: Pointer to the results.

Indicates if the online readout count has been cleared.

Valid values are:

- TRUE
Indicates that the online readout count has been cleared.
- FALSE
Indicates that the online readout count has not been cleared.

This page is intentionally blank.

Chapter 8: EM Calibration APIs

This chapter provides information about the EM gain calibration APIs. All functions, data definitions, and structures are located in the `picam_em_calibration.h` file.



NOTE:

The information and APIs described within this chapter are **NOT** applicable to emICCD cameras.

8.1 EM Calibration Applications

Each ProEM camera is factory-calibrated for linear EM Gain. Over time, however, aging of the EMCCD array may degrade gain linearity. Because aging appears to be a strong function of the amount of charge that flows through the multiplication register, users who consistently operate the camera at high gain at high light levels may need to recalibrate EM gain more frequently than those who are looking at lower light levels at lower gain.

To compensate for aging, each ProEM includes a built-in shutter (either manual or electro-mechanical) and a light source that allows users to perform on-demand EM Gain Calibration using a calibration application. Once the EM gain calibration has been performed, the gain value entered in the software by the user will be the actual multiplication gain applied to the input signal.



CAUTION!

When calibrating a ProEM camera with a manual shutter (e.g., ProEM:1600,) the shutter **MUST** be closed manually before launching any calibration program.

This is not necessary for a camera with an internal electro-mechanical shutter because the program will automatically close the shutter before beginning the calibration.

PICam users have two options available when creating an EM Calibration application:

- Build the sample code `EMGainCalibration.exe` that is included with PICam. This option requires the least amount of development time and overhead since `EMGainCalibration.exe` is a fully functional application once it has been built.



REFERENCES:

Refer to [Appendix B, EM Gain Calibration Code Sample](#), for additional information about using the sample code.

- Create a custom EM Calibration application using the API routines, structures, and callbacks described in this chapter.

When building a custom application, the `EMGainCalibration.exe` sample code included with PICam is a good resource for the developer when learning about the EM Calibration API library.

8.2 Structure Definitions

This section provides programming information about PICam structure definitions.

8.2.1 EM Calibration Structures

This section provides detailed programming information about the following EM Calibration data structures:

- `PicamEMCalibrationDate`

8.2.1.1 *PicamEMCalibrationDate*

Description

`PicamEMCalibrationDate` specifies the calibration date.

Structure Definition

The structure definition for `PicamEMCalibrationDate` is:

```
typedef struct PicamEMCalibrationDate
{
    pint  year;
    pint  month;
    pint  day;
} PicamEMCalibrationDate;
```

Variable Definitions

The variables required by `PicamEMCalibrationDate` are:

- `year`: The year as an integer (e.g., 2011.)
- `month`: The month as an integer.
Valid values are from [1...12], inclusive.
For example, 3 = March.
- `day`: The day of the month as an integer.
Valid values are from [1...31], inclusive.

8.3 Callback Functions

This section provides programming information about callback functions used by PICam

8.3.1 EM Calibration

This section provides information about the following callback functions:

- `PicamEMCalibrationCallback()`.

8.3.1.1 *PicamEMCalibrationCallback()*

Description

`PicamEMCalibrationCallback()` is the callback function for EM calibration progress and/or cancellation.

Syntax

The syntax for `PicamEMCalibrationCallback()` is:

```
typedef pibln (PIL_CALL* PicamEMCalibrationCallback)(  
    PicamHandle  calibration,  
    piflt  progress,  
    void*  user_state );
```

Input Parameters

The input parameters for `PicamEMCalibrationCallback()` are:

`calibration`: Handle for the camera which is being calibrated.

`progress`: This is the percentage of calibration completion.
Valid values are [0...100], inclusive.

`user_state`: User-supplied data provided when calibration is started.

Return Values

Return values for `PicamEMCalibrationCallback()` are:

`TRUE`: Calibration continues.

`FALSE`: Cancels the calibration.

8.4 Programmers' Reference for EM Calibration APIs

This section provides a detailed programmers' reference guide for the following EM Calibration APIs:

- EM Calibration Access APIs
 - `PicamEMCalibration_OpenCalibration()`
 - `PicamEMCalibration_CloseCalibration()`
 - `PicamEMCalibration_GetOpenCalibrations()`
 - `PicamEMCalibration_GetCameraID()`
- EM Calibration Parameter Value APIs
 - `PicamEMCalibration_GetCalibrationDate()`
 - `PicamEMCalibration_ReadSensorTemperatureReading()`
 - `PicamEMCalibration_ReadSensorTemperatureStatus()`
 - `PicamEMCalibration_GetSensorTemperatureSetPoint()`
 - `PicamEMCalibration_SetSensorTemperatureSetPoint()`
- EM Calibration Parameter Constraints APIs
 - `PicamEMCalibration_GetSensorTemperatureSetPointConstraint()`
- EM Calibration APIs
 - `PicamEMCalibration_SetSensorTemperatureSetPoint()`

8.4.1 EM Calibration Access APIs

This section provides programming information about EM Calibration Access APIs.

8.4.1.1 *PicamEMCalibration_OpenCalibration()*

Description

`PicamEMCalibration_OpenCalibration()` opens a camera for calibration and returns a handle to it.



NOTE:

Opening a camera for calibration is mutually exclusive with opening it for normal usage.

Syntax

The syntax for `PicamEMCalibration_OpenCalibration()` is:

```
PICAM_API PicamEMCalibration_OpenCalibration(  
    const PicamCameraID* id,  
    PicamHandle* calibration );
```

Input Parameters

Input parameters for `PicamEMCalibration_OpenCalibration()` are:

`id`: Pointer to the camera id for the camera being calibrated.

Output Parameters

Output parameters for `PicamEMCalibration_OpenCalibration()` are:

`calibration`: Pointer to the handle assigned to the camera that will be calibrated.

Related APIs

For additional information, refer to the following related APIs:

- `PicamEMCalibration_CloseCalibration()`.

8.4.1.2 *PicamEMCalibration_CloseCalibration()*

Description

`PicamEMCalibration_CloseCalibration()` releases all resources that have been associated with a specified calibration process.

Syntax

The syntax for `PicamEMCalibration_CloseCalibration()` is:

```
PICAM_API PicamEMCalibration_CloseCalibration(  
    PicamHandle calibration );
```

Input Parameters

Input parameters for `PicamEMCalibration_CloseCalibration()` are:

`calibration`: Pointer to the handle for the calibration process for which resources are to be released.

Output Parameters

There are no output parameters associated with `PicamEMCalibration_CloseCalibration()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamEMCalibration_OpenCalibration()`.

8.4.1.3 *PicamEMCalibration_GetOpenCalibrations()*

Description

`PicamEMCalibration_GetOpenCalibrations()` returns an allocated array of open calibration handles.

Syntax

The syntax for `PicamEMCalibration_GetOpenCalibrations()` is:

```
PICAM_API PicamEMCalibration_GetOpenCalibrations(  
    const PicamHandle**  calibrations_array,  
    piint* calibrations_count );
```

Input Parameters

There are no input parameters associated with `PicamEMCalibration_GetOpenCalibrations()`.

Output Parameters

Output parameters for `PicamEMCalibration_GetOpenCalibrations()` are:

- `calibrations_array`: Pointer to the array of handles to open calibration processes.
Returns null when there are no open calibration processes.
NOTE: This memory is allocated by PICam and must be released by calling `Picam_DestroyHandles()`
- `calibrations_count`: Pointer to the memory location in which the number of open calibration processes is stored.
Returns 0 when there are no open calibration processes.

8.4.1.4 *PicamEMCalibration_GetCameraID()*

Description

`PicamEMCalibration_GetCameraID()` returns the camera id associated with a specified calibration process.

Syntax

The syntax for `PicamEMCalibration_GetCameraID()` is:

```
PICAM_API PicamEMCalibration_GetCameraID(  
    PicamHandle calibration,  
    PicamCameraID* id );
```

Input Parameters

Input parameters for `PicamEMCalibration_GetCameraID()` are:

`calibration`: Handle associated with the calibration process for which the associated camera is to be determined.

Output Parameters

Output parameters for `PicamEMCalibration_GetCameraID()` are:

`id`: Pointer to the ID of the camera associated with the specified calibration process.

8.4.2 EM Calibration Parameter Value APIs

This section provides programming information about EM Calibration Parameter Value APIs.

8.4.2.1 *PicamEMCalibration_GetCalibrationDate()*

Description

`PicamEMCalibration_GetCalibrationDate()` returns the date of the most recent successful calibration.

Syntax

The syntax for `PicamEMCalibration_GetCalibrationDate()` is:

```
PICAM_API PicamEMCalibration_GetCalibrationDate(  
    PicamHandle calibration,  
    PicamEMCalibrationDate* value );
```

Input Parameters

Input parameters for `PicamEMCalibration_GetCalibrationDate()` are:

`calibration`: Handle of the camera for which the calibration date is to be determined.

Output Parameters

Output parameters for `PicamEMCalibration_GetCalibrationDate()` are:

`value`: Pointer to the calibration date.

Related Structures

For additional information, refer to the following related APIs:

- `PicamEMCalibrationDate`.

8.4.2.2 *PicamEMCalibration_ReadSensorTemperatureReading()*

Description

`PicamEMCalibration_ReadSensorTemperatureReading()` returns the current sensor temperature, in degrees Celsius, for a specified camera.

Syntax

The syntax for `PicamEMCalibration_ReadSensorTemperatureReading()` is:

```
PICAM_API PicamEMCalibration_ReadSensorTemperatureReading(  
    PicamHandle calibration,  
    piflt* value );
```

Input Parameters

Input parameters for `PicamEMCalibration_ReadSensorTemperatureReading()` are:

calibration: Handle of the camera for which the sensor temperature is to be determined.

Output Parameters

Output parameters for `PicamEMCalibration_ReadSensorTemperatureReading()` are:

value: Pointer to the memory location in which the sensor temperature is stored.

8.4.2.3 *PicamEMCalibration_ReadSensorTemperatureStatus()*

Description

`PicamEMCalibration_ReadSensorTemperatureStatus()` returns the status of the current sensor temperature for a specified camera.



NOTE:

Calibration cannot begin until the status of the current sensor temperature is **locked**.

Syntax

The syntax for `PicamEMCalibration_ReadSensorTemperatureStatus()` is:

```
PICAM_API PicamEMCalibration_ReadSensorTemperatureStatus(  
    PicamHandle calibration,  
    PicamSensorTemperatureStatus* value );
```

Input Parameters

Input parameters for `PicamEMCalibration_ReadSensorTemperatureStatus()` are:

calibration: Handle of the camera for which the status of the sensor temperature is to be determined.

Output Parameters

Output parameters for `PicamEMCalibration_ReadSensorTemperatureStatus()` are:

value: Pointer to the memory location in which the status information is stored.

Related Structures

For additional information, refer to the following related APIs:

- `PicamSensorTemperatureStatus`.

8.4.2.4 *PicamEMCalibration_GetSensorTemperatureSetPoint()*

Description

`PicamEMCalibration_GetSensorTemperatureSetPoint()` returns the temperature setpoint that has been programmed for a specified camera.

Syntax

The syntax for `PicamEMCalibration_GetSensorTemperatureSetPoint()` is:

```
PICAM_API PicamEMCalibration_GetSensorTemperatureSetPoint(  
    PicamHandle calibration,  
    piflt* value );
```

Input Parameters

Input parameters for `PicamEMCalibration_GetSensorTemperatureSetPoint()` are:

calibration: Handle of the camera for which the programmed temperature setpoint is to be determined.

Output Parameters

Output parameters for `PicamEMCalibration_GetSensorTemperatureSetPoint()` are:

value: Pointer to the memory location in which the setpoint information is stored.

Related APIs

For additional information, refer to the following related APIs:

- `PicamEMCalibration_SetSensorTemperatureSetPoint()`.

8.4.2.5 *PicamEMCalibration_SetSensorTemperatureSetPoint()*

Description

`PicamEMCalibration_SetSensorTemperatureSetPoint()` configures the sensor temperature setpoint for a specified camera to a specified value.

Syntax

The syntax for `PicamEMCalibration_SetSensorTemperatureSetPoint()` is:

```
PICAM_API PicamEMCalibration_SetSensorTemperatureSetPoint(  
    PicamHandle calibration,  
    piflt value );
```

Input Parameters

Input parameters for `PicamEMCalibration_SetSensorTemperatureSetPoint()` are:

`calibration`: Handle of the camera for which the temperature setpoint is to be programmed.

`value`: The desired temperature setpoint, in degrees Celsius.

Output Parameters

There are no output parameters associated with

`PicamEMCalibration_SetSensorTemperatureSetPoint()`.

8.4.3 EM Calibration Parameter Constraints APIs

This section provides programming information about EM Calibration Parameter Constraint APIs.

8.4.3.1 *PicamEMCalibration_GetSensorTemperatureSetPointConstraint()*

Description

[PicamEMCalibration_GetSensorTemperatureSetPointConstraint\(\)](#) returns an allocated constraint in which the set of valid temperature setpoints, in degrees Celsius, for a specified camera is stored.

Syntax

The syntax for [PicamEMCalibration_GetSensorTemperatureSetPointConstraint\(\)](#) is:

```
PICAM_API PicamEMCalibration_GetSensorTemperatureSetPointConstraint(  
    PicamHandle calibration,  
    const PicamRangeConstraint** constraint );
```

Input Parameters

Input parameters for

[PicamEMCalibration_GetSensorTemperatureSetPointConstraint\(\)](#) are:

calibration: Handle for the camera for which the valid range of temperature setpoints is to be returned.

Output Parameters

Output parameters for

[PicamEMCalibration_GetSensorTemperatureSetPointConstraint\(\)](#) are:

constraint: Pointer to the allocated constraint in which the set of valid temperature setpoints is stored.

NOTE: This memory is allocated by PICam and must be released by calling [Picam_DestroyCollectionConstraints\(\)](#)

Related APIs

For additional information, refer to the following related APIs:

- [Picam_DestroyHandles\(\)](#).

Related Structures

For additional information, refer to the following related APIs:

- [PicamRangeConstraint](#).

8.4.4 EM Calibration APIs

This section provides programming information about EM Calibration APIs.

8.4.4.1 *PicamEMCalibration_Calibrate()*

Description

`PicamEMCalibration_Calibrate()` calibrates the EM Gain for a specified camera.



NOTE:

Calibration cannot begin until the status of the current sensor temperature is **locked**.



NOTE:

If calibration is cancelled (via the use of the callback function `PicamEMCalibrationCallback()`) this function returns `PicamError_OperationCanceled`.

Syntax

The syntax for `PicamEMCalibration_Calibrate()` is:

```
PICAM_API PicamEMCalibration_Calibrate(
                                PicamHandle calibration,
                                PicamEMCalibrationCallback callback,
                                void* user_state );
```

Input Parameters

Input parameters for `PicamEMCalibration_Calibrate()` are:

`calibration`: Handle for the camera for which EM calibration is to be performed.

`callback`: Optional Callback function.

Specifying a Callback provides additional functionality, such as:

- The ability to cancel a calibration process;
- The ability to obtain calibration progress information.

`user_state`: [optional]

When used, allows the caller to provide user-defined data to the callback function.

Output Parameters

There are no output parameters associated with `PicamEMCalibration_Calibrate()`.

Related APIs

For additional information, refer to the following related APIs:

- `PicamEMCalibration_ReadSensorTemperatureStatus()`;
- `PicamEMCalibrationCallback()`.

This page is intentionally blank.

Appendix A: Available Parameters



NOTES:

- Parameters are listed using a truncated version of their names (i.e., the `PicamParameter_` prefix has been dropped.)

For example, the parameter named `PicamParameter_ExposureTime` is listed as `ExposureTime`.

- An asterisk indicates that the parameter does not apply to all members of a camera family.

Table A-1: Symbol Key for Table A-2

Value Types		Constraint Types	
F = Floating Point	M = Modulations	R = Range	M = Modulation
E = Enumeration	R = Region of Interest	C = Collection	P = Pulse
B = Boolean	P = Pulse	Ri = Region of Interest	
I = Integer	L = Large Integer		

Table A-2: Parameter Information and Camera Support (Sheet 1 of 6)

Parameter Name	Read Only	Value Type	Constraint Type	NIRvana LN	PI-MAX 3/4	PI-MTE	PIoNIR/NIRvana/ST	PIXIS	ProEM/ + / -HS	PyLoN	PyLoN-IR	SOPHIA	Quad-RO
Shutter Timing													
ActiveShutter		E	C									✓	
ExposureTime		F	R	✓		✓	✓	✓	✓	✓	✓	✓	✓
ExternalShutterStatus	✓	E										✓	
ExternalShutterType	✓	E										✓	
InactiveShutterTimingModeResult	✓	E										✓*	
InternalShutterStatus	✓	E										✓	
InternalShutterType	✓	E										✓	
ShutterClosingDelay		F	R	✓		✓	✓	✓	✓	✓	✓	✓	✓
ShutterDelayResolution		F	C	✓		✓	✓	✓	✓	✓	✓	✓	✓

Table A-2: Parameter Information and Camera Support (Sheet 2 of 6)

Parameter Name	Read Only	Value Type	Constraint Type	NIRvana LN	PI-MAX 3/4	PI-MTE	PIoNIR/NIRvana/ST	PIXIS	ProEM/ + / -HS	PyLoN	PyLoN-IR	SOPHIA	Quad-RO
ShutterOpeningDelay		F	R	✓			✓		✓	✓	✓	✓	
ShutterTimingMode		E	C	✓		✓	✓	✓	✓	✓	✓	✓	✓
Gating													
DifEndingGate		P	P		✓*								
DifStartingGate		P	P		✓*								
GatingMode		E	C		✓								
RepetitiveGate		P	P		✓								
SequentialEndingGate		P	P		✓								
SequentialGateStepCount		L	R		✓								
SequentialGateStepIterations		L	R		✓								
SequentialStartingGate		P	P		✓								
Intensifier													
BracketGating		B	C		✓*								
CustomModulationSequence		M	M		✓*								
EMICcdGain		I	R		✓*								
EMICcdGainControlMode		E	C		✓*								
EnableIntensifier		B	C		✓								
EnableModulation		B	C		✓*								
GatingSpeed	✓	E			✓								
IntensifierDiameter	✓	F			✓								
IntensifierGain		I	R		✓								
IntensifierOptions	✓	E			✓								
IntensifierStatus	✓	E			✓								
ModulationDuration		F	R		✓*								
ModulationFrequency		F	R		✓*								
PhosphorDecayDelay		F	R		✓								
PhosphorDecayDelayResolution		F	C		✓								
PhosphorType	✓	E			✓								
PhotocathodeSensitivity	✓	E			✓								
RepetitiveModulationPhase		F	R		✓*								

Table A-2: Parameter Information and Camera Support (Sheet 3 of 6)

Parameter Name	Read Only	Value Type	Constraint Type	NIRvana LN	PI-MAX 3/4	PI-MTE	PIoNIR/NIRvana/ST	PIXIS	ProEM/ + / -HS	PyLoN	PyLoN-IR	SOPHIA	Quad-RO
SequentialEndingModulationPhase		F	R		✓*								
SequentialStartingModulationPhase		F	R		✓*								
Analog to Digital Conversion													
AdcAnalogGain		E	C		✓		✓	✓	✓	✓	✓	✓	✓
AdcBitDepth		I	C	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
AdcEMGain		I	R		✓*				✓				
AdcQuality		E	C		✓*	✓*		✓*	✓*	✓			
AdcSpeed		F	C	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
CorrectPixelBias		B	C		✓*				✓	✓		✓	
Hardware I/O													
AnticipateTrigger		B	C		✓*								
AuxOutput		P	P		✓								
DelayFromPreTrigger		F	R		✓*								
EnableModulationOutputSignal		B	C		✓*								
EnableSyncMaster		B	C		✓								
InvertOutputSignal		B	C	✓	✓		✓		✓	✓	✓	✓	
InvertOutputSignal2		B	C									✓	
ModulationOutputSignalAmplitude		F	R		✓*								
ModulationOutputSignalFrequency		F	R		✓*								
OutputSignal		E	C	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
OutputSignal2		E	C									✓	
SyncMaster2Delay		F	R		✓								
TriggerCoupling		E	C		✓								
TriggerDetermination		E	C	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
TriggerFrequency		F	R		✓								
TriggerResponse		E	C	✓	✓*	✓	✓	✓	✓	✓	✓	✓	✓
TriggerSource		E	C		✓								
TriggerTermination		E	C		✓								
TriggerThreshold		F	R		✓								

Table A-2: Parameter Information and Camera Support (Sheet 5 of 6)

Appendix B: EM Gain Calibration Code Sample



CAUTION!

The information provided within this appendix is **NOT** applicable to emICCD cameras.

The `EMGainCalibration.exe` file is sample code included with PICam which, when built, allows PICam users to perform an EM Gain Calibration that may occasionally be required by ProEM systems.



NOTE:

Users with access to LightField do not need to build the sample code in order to perform an EM Gain Calibration. LightField includes `EMGainCalibration.exe` as part of its normal installation. The fully-functional executable file is located within the standard LightField installation directory (i.e., where `PrincetonInstruments.LightField.exe` is stored).

`EMGainCalibration.exe` is an excellent alternative for developers who do not need to create a custom EM gain calibration application. Once the sample code has been built, it can be included as part of the standard customer installation process and allows users to perform EM calibration on an as needed basis.



CAUTION!

When calibrating a ProEM camera with a manual shutter (e.g., ProEM:1600,) the shutter **MUST** be closed manually before launching the calibration program.

This is not necessary for a camera with an internal electro-mechanical shutter because the program will automatically close the shutter before beginning the calibration.

B.1 EM Gain Calibration Procedure

Perform the following procedure to perform an EM Gain calibration:

1. If a data acquisition program is running (e.g., custom PICam application, LightField, etc.) close it.
2. Verify that the ProEM camera that is to be calibrated is the only ProEM camera connected to the host computer and that it is turned on.
3. If the camera has a manual shutter, verify that it is closed. If necessary, close it.


4. Launch `EMGainCalibration.exe`.

**CAUTION!**

Do not operate the camera while EM gain calibration is in process.

5. When the **EM Gain Calibration** dialog is displayed, the default temperature for the camera is shown in the **Target** field. See [Figure B-1](#).

Figure B-1: Typical EM Gain Calibration Dialog



4411-0155_0004

If the camera typically operates at a different temperature, manually adjust it as necessary.

6. Once the **Current** temperature reaches the **Target** temperature specified:

- The internal shutter closes;

**NOTE:**

When using a manual shutter, it must be closed **prior** to initiating the calibration procedure.

- The internal light illuminates the sensor;
- A series of data frames is acquired;
- The calibration map is then calculated.

**NOTE:**

Wait until the calibration has completed before launching the data acquisition program. It may take up to 10 minutes for the calibration to be completed.

Appendix C: Firmware Upgrade/Restore

This appendix provides the procedures to upgrade and restore a GigE camera's firmware.



NOTE:

It is strongly recommended that cameras be upgraded one at a time to avoid confusion.

C.1 Firmware Upgrade Procedure

Perform the following procedure to upgrade a GigE camera's firmware to be compatible with PICam 4.x:


1. On the Host Computer, navigate to the following directory:

`c:\program files\princeton instruments\picam\firmware`

2. Double-click on the `Firmware_Upgrade.exe` file to launch the upgrade tool.

The **IP Engine Selection** dialog is displayed. See [Figure C-1](#).

Figure C-1: Firmware Upgrade: Typical IP Engine Selection Dialog



4411-0155_0006


3. Within the **Available IP Engines** field, select the desired IP Engine from the list of available devices.

**NOTE:**

Each IP Engine listed represents one camera.

See [Figure C-2](#).

Figure C-2: Firmware Upgrade: Selecting Device to be Upgraded



4411-0155_0007



4. Once selected, click  to begin the automated firmware upgrade process. The **Updating** dialog is displayed, similar to that shown in [Figure C-3](#).

Figure C-3: Firmware Upgrade: Typical Updating Dialog




4411-0155_0008

As the upgrade proceeds, the dialog displays appropriate messages, and the progress bar provides a visual indication.

5. Once the upgrade is complete, click **OK** on the **Updating** dialog. See [Figure C-4](#).

Figure C-4: Firmware Upgrade: Upgrade Complete



4411-0155_0009

6. Finally, cycle power to the camera to complete the Firmware Upgrade.

C.2 Restore Firmware

In the unlikely event that PICam 3.x firmware must be restored onto a camera, this section provides detailed information about using the Princeton Instruments provided Firmware Restore tool.

C.2.1 Precautions

Unlike the firmware upgrade procedure that requires no special preparation or precautions, restoring PICam 3.x firmware requires some planning to avoid unnecessary complications.

It is strongly recommended that PICam 3.x firmware be restored on all affected GigE cameras before uninstalling PICam 4.x from the host computer. The Firmware Restore Tool is not included with PICam 3.x installations. Uninstalling PICam 4.x from the host computer will completely remove it making it unavailable for use afterward.



NOTE:

If it is anticipated that additional GigE cameras will require a firmware restore after PICam 4.x has been uninstalled, move the Firmware Restore Tool into a non-PICam directory on the host computer prior to unistalling PICam 4.x.


Additionally, Pleora EBUS Support is required in order to run the Firmware Restore Tool. Do not uninstall Pleora EBUS Support.

C.2.2 Procedure

Perform the following procedure:

1. On the Host Computer, navigate to the following directory:
`c:\program files\princeton instruments\picam\firmware`
2. Double-click on the `Firmware_Restore.exe` file to launch the firmware restore tool.
 The **IP Engine Selection** dialog is displayed. See [Figure C-5](#).

Figure C-5: Firmware Restore: Typical IP Engine Selection Dialog




3. Within the **Available IP Engines** field, select the desired IP Engine from the list of available devices. See [Figure C-6](#).



NOTE:

Each IP Engine listed represents one camera.

Figure C-6: Firmware Restore: Selecting Device to be Restored





4. Once selected, click  to begin restoring the firmware. The **Updating** dialog is displayed, similar to that shown in [Figure C-7](#).

Figure C-7: Firmware Restore: Typical Updating Dialog




4411-0155_0010

As the firmware restoration proceeds, the dialog displays appropriate messages, and the progress bar provides a visual indication.

5. Once the firmware has been restored, click **OK** on the **Updating** dialog. See [Figure C-8](#).

Figure C-8: Firmware Restore: Complete



4411-0155_0011

6. Finally, cycle power to the camera to complete the Firmware Restore.

This page is intentionally blank.

Appendix D: Debugging GigE Cameras

Beginning with PICam 4.x, Princeton Instruments' GigE cameras incorporate a **Camera Heartbeat** that enables the camera and PICam to coordinate communication with one another. As long as a Heartbeat signal has been received by the camera before the specified Heartbeat Timeout has expired, the Camera will continue to be controlled exclusively by PICam.

Under normal operation, the implementation of a Camera Heartbeat is completely transparent to end-users of GigE cameras. However, developers must be keenly aware of how the Heartbeat Timeout impacts camera availability during successive debugging sessions when a PICam process either crashes or is intentionally killed.

Once PICam tasks have been completed, an orderly cessation of the Camera Heartbeat is initiated, and the communication channel between PICam and the camera is closed. At this point, if desired, the camera can again be controlled by PICam or by another program/device.

If, however, PICam halts unexpectedly (e.g., it crashes, the process is killed,) the camera will continue to wait for the next incoming Heartbeat signal until such time as the Heartbeat Timeout has expired. While waiting, the camera remains unavailable to other processes, devices, and programs.



NOTE:

The primary symptom of an expired Heartbeat is a continuous string of unexpected errors being received.

Only after the Heartbeat Timeout has expired without an incoming Heartbeat signal having been detected will the camera close its communication channel with PICam and become available to other processes or devices. At this point, PICam will need to be reinitialized/restarted.

D.1 Debugging

The introduction of the **Camera Heartbeat** presents additional challenges to developers during the Debugging stage of software development. When a breakpoint is reached, PICam execution halts, and no additional Camera Heartbeats are sent to the camera. If the configured Heartbeat Timeout is too short, it will expire, the camera will close its communication channel with PICam requiring it to be reinitialized, and thus preventing subsequent sections of code from easily being executed, examined, and debugged.

The solution to this dilemma is to extend the timeout period sufficiently by configuring the Heartbeat Timeout for an appropriately large value (e.g., 5 minutes.) Increasing the timeout permits the executed code to be examined/debugged while the camera waits patiently for the next incoming Heartbeat signal from PICam. As long as execution of the next chunk of code has been manually initiated before the Heartbeat Timeout expires, PICam will send another Heartbeat signal to the camera (which, in turn, resets the Heartbeat Timeout timer,) and the next chunk of code executes until it reaches the next breakpoint.

D.1.1 Timeout Period Considerations

When deciding on an appropriate timeout period, achieving a balance between having adequate time to review/debug each section of code while not consistently timing out can be tricky.

If too long of a timeout period has been selected and the PICam process crashes or is subsequently killed (a typical action following any debugging session,) the GigE camera will remain unavailable to a future debugging session until the Heartbeat Timeout has expired.

In order to immediately release the camera following a crashed/killed process, cycling its power will clean up any processes that have been abnormally terminated. However, unless the Heartbeat Timeout is programmed for a shorter time period, if the program experiences a subsequent abnormal termination, the camera will again remain unavailable to future debug sessions, and power will need to be cycled again.

D.1.2 Following Debugging

Once debugging has concluded, be sure to reset the value of the `PICAM_GIGE_TIMEOUT` environment variable to a more appropriate timeout.

D.2 Timeout Configuration

The Heartbeat Timeout, in milliseconds, is defined by the environment variable `PICAM_GIGE_TIMEOUT`.

Valid values are:

- Minimum: 500 ms
- Default: 2000 ms
- Maximum: 4,294,967,295 ms (approximately 49.7 days)

This page is intentionally blank.



www.princetoninstruments.com

info@princetoninstruments.com

USA +1 877 474 2286 | FRANCE +33 (1) 60 86 03 65 | GERMANY +49 (0) 89 660 7793

UK & IRELAND +44 (0) 1628 472 346 | SINGAPORE +65 6408 6240 | CHINA +86 10 659 16460 | JAPAN +81 (3) 5639 2741